

Table of contents

\\CMIC

Physics

T2MES

T2MES_low_res

tfl_loc_multi_iPAT@c			
GRE	FA		120
GRE	FA		60
highres_3D_of_tubes			
BH_FW_4e_GT			
pre_MOLLI_5s(3s)3s_256			
post_MOLLI_4s(1s)3s(1s)2s_256			
T2map	0-25-55	new	SSFP
pre		FB	T1T2
pre	FB	T1T2	high
post		FB	T1T2
pre_MOLLI_5s(3s)3s_256			
post_MOLLI_4s(1s)3s(1s)2s_256			
T2map	0-25-55	new	SSFP
pre		FB	T1T2
pre	FB	T1T2	high
post		FB	T1T2
pre_MOLLI_5s(3s)3s_256			
post_MOLLI_4s(1s)3s(1s)2s_256			
T2map	0-25-55	new	SSFP
pre		FB	T1T2
pre	FB	T1T2	high
post		FB	T1T2
pre_MOLLI_5s(3s)3s_256			
post_MOLLI_4s(1s)3s(1s)2s_256			
T2map	0-25-55	new	SSFP
pre		FB	T1T2
pre	FB	T1T2	high
post		FB	T1T2
BH_FW_4e_GT			
se	te	=	10
se	te	=	20
se	te	=	30
se	te	=	40
se	te	=	50
se	te	=	60
se	te	=	80
se	te	=	100
se	te	=	125
se	te	=	150
IR		GRE	TI=20
IR		GRE	TI=50
IR		GRE	TI=100
IR		GRE	TI=200
IR		GRE	TI=400
IR		GRE	TI=600
IR		GRE	TI=800
IR		GRE	TI=1000
IR		GRE	TI=1300
IR		GRE	TI=1700

IR	GRE
BH_FW_4e2i_GT	
highres_3D_of_tubes_x5	
FAMAP120HIGHBWTP_90	
FAMAP60HIGHBWTP_45	

TI=2100

\\CMIC\Physics\T2MES\T2MES_low_res\tfl_loc_multi_iPAT@c

TA: 4.1 s PM: ISO Voxel size: 2.1×2.1×8.0 mmPAT: 2 Rel. SNR: 1.00 : tfl

Properties

Prio recon	Off
Load images to viewer	On
Inline movie	Off
Auto store images	On
Load images to stamp segments	On
Load images to graphic segments	On
Auto open inline display	Off
Auto close inline display	Off
Start measurement without further preparation	On
Wait for user to start	On
Start measurements	Single measurement

Routine

Slice group	1
Slices	3
Dist. factor	300 %
Position	Isocenter
Orientation	Transversal
Phase enc. dir.	A >> P
Slice group	2
Slices	1
Dist. factor	20 %
Position	Isocenter
Orientation	Sagittal
Phase enc. dir.	A >> P
Slice group	3
Slices	3
Dist. factor	200 %
Position	Isocenter
Orientation	Coronal
Phase enc. dir.	R >> L
AutoAlign	---
Phase oversampling	0 %
FoV read	400 mm
FoV phase	91.7 %
Slice thickness	8.0 mm
TR	500.00 ms
TE	1.46 ms
Averages	1
Concatenations	7
Filter	Distortion Corr.(2D), Prescan Normalize, Elliptical filter
Coil elements	BO1-3;SP1-3

Contrast - Common

TR	500.00 ms
TE	1.46 ms
TD	0 ms
Magn. preparation	None
Flip angle	8 deg
Fat suppr.	None
Wrap-up Magn.	None

Contrast - Dynamic

Averages	1
Averaging mode	Long term
Reconstruction	Magnitude
Measurements	1
Multiple series	Off

Resolution - Common

FoV read	400 mm
FoV phase	91.7 %
Slice thickness	8.0 mm
Base resolution	192
Phase resolution	50 %
Phase partial Fourier	Off
Trajectory	Cartesian
Interpolation	Off

Resolution - iPAT

PAT mode	GRAPPA
Accel. factor PE	2
Ref. lines PE	24
Reference scan mode	Integrated

Resolution - Filter Image

Image Filter	Off
Distortion Corr.	On
Mode	2D
Unfiltered images	Off
Prescan Normalize	On
Unfiltered images	Off
Normalize	Off
B1 filter	Off

Resolution - Filter Rawdata

Raw filter	Off
Elliptical filter	On
POCS	Off

Geometry - Common

Slice group	1
Slices	3
Dist. factor	300 %
Position	Isocenter
Orientation	Transversal
Phase enc. dir.	A >> P
Slice group	2
Slices	1
Dist. factor	20 %
Position	Isocenter
Orientation	Sagittal
Phase enc. dir.	A >> P
Slice group	3
Slices	3
Dist. factor	200 %
Position	Isocenter
Orientation	Coronal
Phase enc. dir.	R >> L
FoV read	400 mm
FoV phase	91.7 %
Slice thickness	8.0 mm
TR	500.00 ms
Multi-slice mode	Sequential
Series	Ascending
Concatenations	7

Geometry - AutoAlign

Slice group	1
Position	Isocenter

Geometry - AutoAlign

Orientation	Transversal
Phase enc. dir.	A >> P
Slice group	2
Position	Isocenter
Orientation	Sagittal
Phase enc. dir.	A >> P
Slice group	3
Position	Isocenter
Orientation	Coronal
Phase enc. dir.	R >> L
AutoAlign	---
Initial Position	Isocenter
L	0.0 mm
P	0.0 mm
H	0.0 mm
Initial Rotation	0.00 deg
Initial Orientation	Coronal

Geometry - Saturation

Fat suppr.	None
Wrap-up Magn.	None
Special sat.	None

Geometry - Navigator**System - Miscellaneous**

Positioning mode	ISO
Table position	H
Table position	0 mm
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Coil Combine Mode	Adaptive Combine
Save uncombined	Off
Matrix Optimization	Off
Coil Focus	Flat
AutoAlign	---
Coil Select Mode	Default

System - Adjustments

B0 Shim mode	Tune up
Adjust with body coil	Off
Confirm freq. adjustment	Off
Assume Dominant Fat	Off
Assume Silicone	Off
Adjustment Tolerance	Auto

System - Adjust Volume

Position	Isocenter
Orientation	Transversal
Rotation	0.00 deg
A >> P	263 mm
R >> L	350 mm
F >> H	350 mm
Reset	Off

System - Tx/Rx

Frequency 1H	63.673772 MHz
Correction factor	1
Gain	High
Img. Scale Cor.	1.000
Reset	Off
? Ref. amplitude 1H	0.000 V

Physio - Signal1

1st Signal/Mode	None
TR	500.00 ms
Concatenations	7
Segments	56

Physio - Cardiac

Tagging	None
Magn. preparation	None
Fat suppr.	None
Dark blood	On
Dark blood thickness	200 %
FoV read	400 mm
FoV phase	91.7 %
Phase resolution	50 %
Cine	Off
Trajectory	Cartesian
Dummy heartbeats	1

Physio - PACE

Resp. control	Off
Concatenations	7

Inline - Common

Subtract	Off
Measurements	1
StdDev	Off
Save original images	On

Inline - Cardiac

Inline Evaluation	Off
Magn. preparation	None
Contrasts	1
TE	1.46 ms
TR	500.00 ms
Save original images	On

Inline - MIP

MIP-Sag	Off
MIP-Cor	Off
MIP-Tra	Off
MIP-Time	Off
Save original images	On

Inline - Composing

Inline Composing	Off
Distortion Corr.	On
Mode	2D
Unfiltered images	Off

Sequence - Part 1

Introduction	Off
Dimension	2D
Reordering	Linear
Asymmetric echo	Weak
Contrasts	1
Flow comp.	No
Optimization	Min. TE
Multi-slice mode	Sequential
Echo spacing	3.7 ms
Sequence type	Gre
Bandwidth	473 Hz/Px

Sequence - Part 2

Define	Shots
Shots per slice	1
Segments	56
RF pulse type	Fast
Gradient mode	Fast
Excitation	Slice-sel.
Flip angle mode	Constant
RF spoiling	On
Phase Enc. Rewinder	On
Cine	Off

Sequence - Assistant

Mode	Off
Allowed delay	0 s

\\CMIC\Physics\T2MES\T2MES_low_res\GRE FA 120

TA: 18:00 PM: FIX Voxel size: 1.9×1.9×8.0 mmRel. SNR: 1.00 : fl

Properties

Prio recon	Off
Load images to viewer	On
Inline movie	Off
Auto store images	On
Load images to stamp segments	Off
Load images to graphic segments	Off
Auto open inline display	Off
Auto close inline display	Off
Start measurement without further preparation	Off
Wait for user to start	Off
Start measurements	Single measurement

Routine

Slice group	1
Slices	1
Position	Isocenter
Orientation	Transversal
Phase enc. dir.	A >> P
AutoAlign	---
FoV read	360 mm
FoV phase	75.0 %
Slice thickness	8.0 mm
TR	10000.00 ms
TE	10.00 ms
Filter	Distortion Corr.(2D)
Coil elements	BO1-3;SP2-4

Contrast - Common

TR	10000.00 ms
TE	10.00 ms
Magn. preparation	None
Flip angle	120 deg

Resolution - Common

FoV read	360 mm
FoV phase	75.0 %
Slice thickness	8.0 mm
Base resolution	192
Phase resolution	75 %

Resolution - Filter Image

Image Filter	Off
Distortion Corr.	On
Mode	2D
Unfiltered images	Off
Prescan Normalize	Off
Normalize	Off
B1 filter	Off

Resolution - Filter Rawdata

Raw filter	Off
Elliptical filter	Off

Geometry - Common

Slice group	1
Slices	1
Position	Isocenter
Orientation	Transversal
Phase enc. dir.	A >> P

Geometry - Common

FoV read	360 mm
FoV phase	75.0 %
Slice thickness	8.0 mm
TR	10000.00 ms

Geometry - AutoAlign

Slice group	1
Position	Isocenter
Orientation	Transversal
Phase enc. dir.	A >> P
AutoAlign	---
Initial Position	Isocenter
L	0.0 mm
P	0.0 mm
H	0.0 mm
Initial Rotation	0.00 deg
Initial Orientation	Transversal

System - Miscellaneous

Positioning mode	FIX
Table position	H
Table position	0 mm
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Coil Combine Mode	Adaptive Combine
Save uncombined	Off
Matrix Optimization	Off
Coil Focus	Flat
AutoAlign	---
Coil Select Mode	Default

System - Adjustments

B0 Shim mode	Cardiac
Adjust with body coil	On
Confirm freq. adjustment	Off
Assume Dominant Fat	Off
Assume Silicone	Off
Adjustment Tolerance	Auto

System - Adjust Volume

! Position	Isocenter
! Orientation	Transversal
! Rotation	0.00 deg
! A >> P	100 mm
! R >> L	100 mm
! F >> H	50 mm
Reset	Off

System - Tx/Rx

Frequency 1H	63.673772 MHz
Correction factor	1
Gain	High
Img. Scale Cor.	1.000
Reset	Off
? Ref. amplitude 1H	0.000 V

Sequence - Part 1

Bandwidth	500 Hz/Px
-----------	-----------

Sequence - Special

Adiab IR Pulse	Tan/Tanh
RF Phase Spoiling	ON

Sequence - Assistant

Mode	Off
------	-----

\\CMIC\Physics\T2MES\T2MES_low_res\GRE FA 60

TA: 18:00 PM: FIX Voxel size: 1.9×1.9×8.0 mmRel. SNR: 1.00 : fl

Properties

Prio recon	Off
Load images to viewer	On
Inline movie	Off
Auto store images	On
Load images to stamp segments	Off
Load images to graphic segments	Off
Auto open inline display	Off
Auto close inline display	Off
Start measurement without further preparation	On
Wait for user to start	Off
Start measurements	Single measurement

Routine

Slice group	1
Slices	1
Position	Isocenter
Orientation	Transversal
Phase enc. dir.	A >> P
AutoAlign	---
FoV read	360 mm
FoV phase	75.0 %
Slice thickness	8.0 mm
TR	10000.00 ms
TE	10.00 ms
Filter	Distortion Corr.(2D)
Coil elements	BO1-3;SP2-4

Contrast - Common

TR	10000.00 ms
TE	10.00 ms
Magn. preparation	None
Flip angle	60 deg

Resolution - Common

FoV read	360 mm
FoV phase	75.0 %
Slice thickness	8.0 mm
Base resolution	192
Phase resolution	75 %

Resolution - Filter Image

Image Filter	Off
Distortion Corr.	On
Mode	2D
Unfiltered images	Off
Prescan Normalize	Off
Normalize	Off
B1 filter	Off

Resolution - Filter Rawdata

Raw filter	Off
Elliptical filter	Off

Geometry - Common

Slice group	1
Slices	1
Position	Isocenter
Orientation	Transversal
Phase enc. dir.	A >> P

Geometry - Common

FoV read	360 mm
FoV phase	75.0 %
Slice thickness	8.0 mm
TR	10000.00 ms

Geometry - AutoAlign

Slice group	1
Position	Isocenter
Orientation	Transversal
Phase enc. dir.	A >> P
AutoAlign	---
Initial Position	Isocenter
L	0.0 mm
P	0.0 mm
H	0.0 mm
Initial Rotation	0.00 deg
Initial Orientation	Transversal

System - Miscellaneous

Positioning mode	FIX
Table position	H
Table position	0 mm
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Coil Combine Mode	Adaptive Combine
Save uncombined	Off
Matrix Optimization	Off
Coil Focus	Flat
AutoAlign	---
Coil Select Mode	Default

System - Adjustments

B0 Shim mode	Cardiac
Adjust with body coil	On
Confirm freq. adjustment	Off
Assume Dominant Fat	Off
Assume Silicone	Off
Adjustment Tolerance	Auto

System - Adjust Volume

! Position	Isocenter
! Orientation	Transversal
! Rotation	0.00 deg
! A >> P	100 mm
! R >> L	100 mm
! F >> H	50 mm
Reset	Off

System - Tx/Rx

Frequency 1H	63.673772 MHz
Correction factor	1
Gain	High
Img. Scale Cor.	1.000
Reset	Off
? Ref. amplitude 1H	0.000 V

Sequence - Part 1

Bandwidth	500 Hz/Px
-----------	-----------

Sequence - Special

Adiab IR Pulse	Tan/Tanh
RF Phase Spoiling	ON

Sequence - Assistant

Mode	Off
------	-----

\\CMIC\Physics\T2MES\T2MES_low_res\highres_3D_of_tubes

TA: 0:58 PM: FIX Voxel size: 1.0×1.0×1.0 mmPAT: 2 Rel. SNR: 1.00 : tfl

Properties

Prio recon	Off
Load images to viewer	On
Inline movie	Off
Auto store images	On
Load images to stamp segments	On
Load images to graphic segments	On
Auto open inline display	Off
Auto close inline display	Off
Start measurement without further preparation	On
Wait for user to start	Off
Start measurements	Single measurement

Routine

Slab group	1
Slabs	1
Dist. factor	20 %
Position	Isocenter
Orientation	Transversal
Phase enc. dir.	A >> P
AutoAlign	---
Phase oversampling	0 %
Slice oversampling	20.0 %
Slices per slab	120
FoV read	256 mm
FoV phase	64.8 %
Slice thickness	1.00 mm
TR	400.90 ms
TE	1.77 ms
Averages	1
Concatenations	1
Filter	None
Coil elements	BO1-3;SP2-4

Contrast - Common

TR	400.90 ms
TE	1.77 ms
Magn. preparation	None
Flip angle	12 deg
Fat suppr.	None
Wrap-up Magn.	None

Contrast - Dynamic

Averages	1
Averaging mode	Short term
Reconstruction	Magnitude
Measurements	1
Multiple series	Off

Resolution - Common

FoV read	256 mm
FoV phase	64.8 %
Slice thickness	1.00 mm
Base resolution	256
Phase resolution	100 %
Slice resolution	100 %
Phase partial Fourier	Off
Slice partial Fourier	Off
Trajectory	Cartesian
Interpolation	Off

Resolution - iPAT

PAT mode	GRAPPA
Accel. factor PE	2
Ref. lines PE	24
Accel. factor 3D	1
Reference scan mode	Integrated

Resolution - Filter Image

Image Filter	Off
Distortion Corr.	Off
Prescan Normalize	Off
Normalize	Off
B1 filter	Off

Resolution - Filter Rawdata

Raw filter	Off
Elliptical filter	Off
POCS	Off

Geometry - Common

Slab group	1
Slabs	1
Dist. factor	20 %
Position	Isocenter
Orientation	Transversal
Phase enc. dir.	A >> P
Slice oversampling	20.0 %
Slices per slab	120
FoV read	256 mm
FoV phase	64.8 %
Slice thickness	1.00 mm
TR	400.90 ms
Multi-slice mode	Sequential
Series	Descending
Concatenations	1

Geometry - AutoAlign

Slab group	1
Position	Isocenter
Orientation	Transversal
Phase enc. dir.	A >> P
AutoAlign	---
Initial Position	Isocenter
L	0.0 mm
P	0.0 mm
H	0.0 mm
Initial Rotation	0.00 deg
Initial Orientation	Transversal

Geometry - Saturation

Fat suppr.	None
Wrap-up Magn.	None
Special sat.	None

Geometry - Navigator**System - Miscellaneous**

Positioning mode	FIX
Table position	H
Table position	0 mm
MSMA	S - C - T

System - Miscellaneous

Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Coil Combine Mode	Sum of Squares
Save uncombined	Off
Matrix Optimization	Off
Coil Focus	Flat
AutoAlign	---
Coil Select Mode	Default

System - Adjustments

B0 Shim mode	Cardiac
Adjust with body coil	On
Confirm freq. adjustment	Off
Assume Dominant Fat	Off
Assume Silicone	Off
Adjustment Tolerance	Auto

System - Adjust Volume

! Position	Isocenter
! Orientation	Transversal
! Rotation	0.00 deg
! A >> P	100 mm
! R >> L	100 mm
! F >> H	40 mm
Reset	Off

System - Tx/Rx

Frequency 1H	63.673772 MHz
Correction factor	1
Gain	Low
Img. Scale Cor.	1.000
Reset	Off
? Ref. amplitude 1H	0.000 V

Physio - Signal1

1st Signal/Mode	None
TR	400.90 ms
Concatenations	1
Segments	95

Physio - Cardiac

Tagging	None
Magn. preparation	None
Fat suppr.	None
Dark blood	Off
FoV read	256 mm
FoV phase	64.8 %
Phase resolution	100 %
Cine	Off
Trajectory	Cartesian
Dummy heartbeats	1

Physio - PACE

Resp. control	Off
Concatenations	1

Inline - Common

Subtract	Off
Measurements	1
StdDev	Off
Save original images	On

Inline - Cardiac

Inline Evaluation	Off
Magn. preparation	None
Contrasts	1
TE	1.77 ms
TR	400.90 ms
Save original images	On

Inline - MIP

MIP-Sag	Off
MIP-Cor	Off
MIP-Tra	Off
MIP-Time	Off
Save original images	On

Inline - Composing

Inline Composing	Off
Distortion Corr.	Off

Sequence - Part 1

Introduction	Off
Dimension	3D
Elliptical scanning	Off
Reordering	Linear
Asymmetric echo	Off
Contrasts	1
Flow comp.	No
Optimization	Min. TE TR
Multi-slice mode	Sequential
Echo spacing	4.2 ms
Sequence type	Gre
Bandwidth	592 Hz/Px

Sequence - Part 2

Define	Shots
Shots per slice	1
Segments	95
RF pulse type	Fast
Gradient mode	Fast
Excitation	Slab-sel.
Flip angle mode	Constant
RF spoiling	On
Phase Enc. Rewinder	On
Cine	Off

Sequence - Assistant

Mode	Off
Allowed delay	0 s

\\CMIC\Physics\T2MES\T2MES_low_res\BH_FW_4e_GT

TA: 8.5 s PM: FIX Voxel size: 1.4×1.4×8.0 mmPAT: Off Rel. SNR: 1.00 : fl

Properties

Prio recon	Off
Load images to viewer	Off
Inline movie	Off
Auto store images	On
Load images to stamp segments	Off
Load images to graphic segments	On
Auto open inline display	Off
Auto close inline display	Off
Start measurement without further preparation	On
Wait for user to start	Off
Start measurements	Single measurement

Routine

Slice group	1
Slices	1
Dist. factor	20 %
Position	Isocenter
Orientation	Transversal
Phase enc. dir.	A >> P
AutoAlign	---
Phase oversampling	0 %
FoV read	360 mm
FoV phase	75.0 %
Slice thickness	8.0 mm
TR	208.00 ms
TE 1	1.56 ms
TE 2	3.88 ms
TE 3	6.2 ms
TE 4	8.52 ms
Averages	1
Concatenations	1
Filter	Distortion Corr.(2D)
Coil elements	BO1-3;SP2-4

Contrast - Common

TR	208.00 ms
TE 1	1.56 ms
TE 2	3.88 ms
TE 3	6.2 ms
TE 4	8.52 ms
Magn. preparation	None
Flip angle	12 deg
Fat suppr.	None
Wrap-up Magn.	None

Contrast - Dynamic

Averages	1
Averaging mode	Short term
Reconstruction	Magnitude
Measurements	1
Multiple series	Off

Resolution - Common

FoV read	360 mm
FoV phase	75.0 %
Slice thickness	8.0 mm
Base resolution	256
Phase resolution	75 %
Phase partial Fourier	Off

Resolution - Common

Trajectory	Cartesian
Interpolation	Off

Resolution - iPAT

PAT mode	None
----------	------

Resolution - Filter Image

Image Filter	Off
Distortion Corr.	On
Mode	2D
Unfiltered images	Off
Prescan Normalize	Off
Normalize	Off
B1 filter	Off

Resolution - Filter Rawdata

Raw filter	Off
Elliptical filter	Off
POCS	Off

Geometry - Common

Slice group	1
Slices	1
Dist. factor	20 %
Position	Isocenter
Orientation	Transversal
Phase enc. dir.	A >> P
FoV read	360 mm
FoV phase	75.0 %
Slice thickness	8.0 mm
TR	208.00 ms
Multi-slice mode	Sequential
Series	Interleaved
Concatenations	1

Geometry - AutoAlign

Slice group	1
Position	Isocenter
Orientation	Transversal
Phase enc. dir.	A >> P
AutoAlign	---
Initial Position	Isocenter
L	0.0 mm
P	0.0 mm
H	0.0 mm
Initial Rotation	0.00 deg
Initial Orientation	Transversal

Geometry - Saturation

Fat suppr.	None
Wrap-up Magn.	None
Special sat.	None

Geometry - Navigator**System - Miscellaneous**

Positioning mode	FIX
Table position	H
Table position	0 mm

System - Miscellaneous

MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Coil Combine Mode	Adaptive Combine
Save uncombined	Off
Matrix Optimization	Off
Coil Focus	Flat
AutoAlign	---
Coil Select Mode	Default

System - Adjustments

B0 Shim mode	Standard
Adjust with body coil	On
Confirm freq. adjustment	Off
Assume Dominant Fat	Off
Assume Silicone	Off
Adjustment Tolerance	Auto

System - Adjust Volume

! Position	Isocenter
! Orientation	Transversal
! Rotation	0.00 deg
! A >> P	100 mm
! R >> L	100 mm
! F >> H	25 mm
Reset	Off

System - Tx/Rx

Frequency 1H	63.673772 MHz
Correction factor	1
Gain	High
Img. Scale Cor.	1.000
Reset	Off
? Ref. amplitude 1H	0.000 V

Physio - Signal1

1st Signal/Mode	ECG/Trigger
Average cycle	947 ± 122 ms
Average cycle	No Signal ms
Captured cycle	947 ± 122 ms
Acquisition window	830 ms
Trigger pulse	1
Trigger delay	622 ms
TR	208.00 ms
Concatenations	1
Segments	20
Phases	1
Adaptive Triggering	Off

Physio - Cardiac

Tagging	None
Magn. preparation	None
Fat suppr.	None
Dark blood	Off
FoV read	360 mm
FoV phase	75.0 %
Phase resolution	75 %
Cine	Off
Trajectory	Cartesian
Dummy heartbeats	1

Physio - PACE

Resp. control	Off
---------------	-----

Physio - PACE

Concatenations	1
----------------	---

Inline - Common

Subtract	Off
Measurements	1
StdDev	Off
Save original images	On

Inline - Cardiac

Inline Evaluation	Off
Magn. preparation	None
Contrasts	4
TE 1	1.56 ms
TE 2	3.88 ms
TE 3	6.2 ms
TE 4	8.52 ms
TR	208.00 ms
Save original images	On

Inline - MIP

MIP-Sag	Off
MIP-Cor	Off
MIP-Tra	Off
MIP-Time	Off
Save original images	On

Inline - Composing

Inline Composing	Off
Distortion Corr.	On
Mode	2D
Unfiltered images	Off

Sequence - Part 1

Introduction	Off
Dimension	2D
Reordering	Linear
Asymmetric echo	Off
Contrasts	4
Flow comp. 1	No
Readout mode	Monopolar
Optimization	Min. TE TR
Multi-slice mode	Sequential
Echo spacing	10.4 ms
Sequence type	Gre
Bandwidth 1	977 Hz/Px
Bandwidth 2	977 Hz/Px
Bandwidth 3	977 Hz/Px
Bandwidth 4	977 Hz/Px

Sequence - Part 2

Define	Segments
Segments	20
RF pulse type	Fast
Gradient mode	Fast
Excitation	Slice-sel.
Flip angle mode	Constant
RF spoiling	On
Phase Enc. Rewinder	On
Cine	Off

Sequence - Special

FatWater Separation	On
Multi-echo Images	On

Sequence - Special

In-Opp Phase Images	Off
Frequency Map	On
T2* Map	Off
Motion Correction	Off
MoCo Averaging Mode	Complex MoCo
MoCo Images Only?	Off
Frequency Offset	0 Hz
No. of Interleaves	0

Sequence - Assistant

Mode	Off
Allowed delay	0 s

\\CMIC\Physics\T2MES\T2MES_low_res\pre_MOLLI_5s(3s)3s_256

TA: 0:11 PM: FIX Voxel size: 1.4×1.4×8.0 mmPAT: 2 Rel. SNR: 1.00 : tfi

Properties

Prio recon	Off
Load images to viewer	On
Inline movie	Off
Auto store images	On
Load images to stamp segments	Off
Load images to graphic segments	On
Auto open inline display	Off
Auto close inline display	Off
Start measurement without further preparation	On
Wait for user to start	Off
Start measurements	Single measurement

Resolution - iPAT

Accel. factor PE	2
Ref. lines PE	36
Reference scan mode	GRE/separate

Resolution - Filter Image

Image Filter	Off
Distortion Corr.	On
Mode	2D
Unfiltered images	Off
Prescan Normalize	Off
Normalize	Off
B1 filter	Off

Routine

Slice group	1
Slices	1
Dist. factor	20 %
Position	Isocenter
Orientation	Transversal
Phase enc. dir.	A >> P
AutoAlign	---
Phase oversampling	0 %
FoV read	360 mm
FoV phase	75.0 %
Slice thickness	8.0 mm
TR	280.56 ms
TE	1.12 ms
Averages	1
Concatenations	1
Filter	Distortion Corr.(2D)
Coil elements	BO1-3;SP2-4

Resolution - Filter Rawdata

Raw filter	Off
Elliptical filter	Off
POCS	Off

Geometry - Common

Slice group	1
Slices	1
Dist. factor	20 %
Position	Isocenter
Orientation	Transversal
Phase enc. dir.	A >> P
FoV read	360 mm
FoV phase	75.0 %
Slice thickness	8.0 mm
TR	280.56 ms
Multi-slice mode	Sequential
Series	Interleaved
Concatenations	1

Contrast - Common

TR	280.56 ms
TE	1.12 ms
Magn. preparation	Non-sel. IR T1map
T1	180 ms
Flip angle	35 deg
Fat suppr.	None
Wrap-up Magn.	None

Geometry - AutoAlign

Slice group	1
Position	Isocenter
Orientation	Transversal
Phase enc. dir.	A >> P
AutoAlign	---
Initial Position	Isocenter
L	0.0 mm
P	0.0 mm
H	0.0 mm
Initial Rotation	0.00 deg
Initial Orientation	Transversal

Contrast - Dynamic

Averages	1
Averaging mode	Short term
Reconstruction	Magnitude
Measurements	1
Multiple series	Off

Resolution - Common

FoV read	360 mm
FoV phase	75.0 %
Slice thickness	8.0 mm
Base resolution	256
Phase resolution	75 %
Phase partial Fourier	7/8
Trajectory	Cartesian
Interpolation	Off

Geometry - Saturation

Fat suppr.	None
Wrap-up Magn.	None
Special sat.	None

Geometry - Navigator**System - Miscellaneous**

Positioning mode	FIX
Table position	H
Table position	0 mm
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P

Resolution - iPAT

PAT mode	GRAPPA
----------	--------

System - Miscellaneous

Transversal	F >> H
Coil Combine Mode	Adaptive Combine
Save uncombined	Off
Matrix Optimization	Off
Coil Focus	Flat
AutoAlign	---
Coil Select Mode	Default

System - Adjustments

B0 Shim mode	Standard
Adjust with body coil	On
Confirm freq. adjustment	Off
Assume Dominant Fat	Off
Assume Silicone	Off
Adjustment Tolerance	Auto

System - Adjust Volume

! Position	Isocenter
! Orientation	Transversal
! Rotation	0.00 deg
! A >> P	100 mm
! R >> L	100 mm
! F >> H	25 mm
Reset	Off

System - Tx/Rx

Frequency 1H	63.673772 MHz
Correction factor	1
Gain	High
Img. Scale Cor.	1.000
Reset	Off
? Ref. amplitude 1H	0.000 V

Physio - Signal1

1st Signal/Mode	ECG/Trigger
Average cycle	947 ± 122 ms
Average cycle	No Signal ms
Captured cycle	947 ± 122 ms
Acquisition window	700 ms
Trigger pulse	1
Trigger delay	418 ms
TR	280.56 ms
Concatenations	1
Segments	72
Phases	1
Adaptive Triggering	Off

Physio - Cardiac

Tagging	None
Magn. preparation	Non-sel. IR T1map
TI	180 ms
Fat suppr.	None
Dark blood	Off
FoV read	360 mm
FoV phase	75.0 %
Phase resolution	75 %
Cine	Off
Trajectory	Cartesian
Dummy heartbeats	0
Motion Correction	Standard

Physio - PACE

Resp. control	Off
Concatenations	1

Sequence - Part 1

Introduction	Off
Dimension	2D
Reordering	Linear
Asymmetric echo	Weak
Contrasts	1
Optimization	Min. TE TR
Multi-slice mode	Sequential
Sequence type	Trufi
Bandwidth	1085 Hz/Px

Sequence - Part 2

Define	Shots
Shots per slice	1
Segments	72
Trufi delta freq.	0 Hz
RF pulse type	Fast
Gradient mode	Fast
Excitation	Slice-sel.
Flip angle mode	Constant
Cine	Off

Sequence - Special

Partition coeff map	Off
ECV map	Off
Synth ECV map	Off
T1 scout mode	Off
Error map	On
Synth PSIR	Off
Periods in seconds	On
16 bit images	Off
T1 sampling scheme	Native

Sequence - Assistant

Mode	Off
Allowed delay	0 s

\\CMIC\Physics\T2MES\T2MES_low_res\post_MOLLI_4s(1s)3s(1s)2s_256

TA: 0:12 PM: FIX Voxel size: 1.4×1.4×8.0 mmPAT: 2 Rel. SNR: 1.00 : tfi

Properties

Prio recon	Off
Load images to viewer	On
Inline movie	Off
Auto store images	On
Load images to stamp segments	Off
Load images to graphic segments	On
Auto open inline display	Off
Auto close inline display	Off
Start measurement without further preparation	On
Wait for user to start	Off
Start measurements	Single measurement

Resolution - iPAT

Accel. factor PE	2
Ref. lines PE	36
Reference scan mode	GRE/separate

Resolution - Filter Image

Image Filter	Off
Distortion Corr.	On
Mode	2D
Unfiltered images	Off
Prescan Normalize	Off
Normalize	Off
B1 filter	Off

Routine

Slice group	1
Slices	1
Dist. factor	20 %
Position	Isocenter
Orientation	Transversal
Phase enc. dir.	A >> P
AutoAlign	---
Phase oversampling	0 %
FoV read	360 mm
FoV phase	75.0 %
Slice thickness	8.0 mm
TR	360.56 ms
TE	1.12 ms
Averages	1
Concatenations	1
Filter	Distortion Corr.(2D)
Coil elements	BO1-3;SP2-4

Resolution - Filter Rawdata

Raw filter	Off
Elliptical filter	Off
POCS	Off

Geometry - Common

Slice group	1
Slices	1
Dist. factor	20 %
Position	Isocenter
Orientation	Transversal
Phase enc. dir.	A >> P
FoV read	360 mm
FoV phase	75.0 %
Slice thickness	8.0 mm
TR	360.56 ms
Multi-slice mode	Sequential
Series	Interleaved
Concatenations	1

Contrast - Common

TR	360.56 ms
TE	1.12 ms
Magn. preparation	Non-sel. IR T1map
T1	260 ms
Flip angle	35 deg
Fat suppr.	None
Wrap-up Magn.	None

Geometry - AutoAlign

Slice group	1
Position	Isocenter
Orientation	Transversal
Phase enc. dir.	A >> P
AutoAlign	---
Initial Position	Isocenter
L	0.0 mm
P	0.0 mm
H	0.0 mm
Initial Rotation	0.00 deg
Initial Orientation	Transversal

Contrast - Dynamic

Averages	1
Averaging mode	Short term
Reconstruction	Magnitude
Measurements	1
Multiple series	Off

Resolution - Common

FoV read	360 mm
FoV phase	75.0 %
Slice thickness	8.0 mm
Base resolution	256
Phase resolution	75 %
Phase partial Fourier	7/8
Trajectory	Cartesian
Interpolation	Off

Geometry - Saturation

Fat suppr.	None
Wrap-up Magn.	None
Special sat.	None

Geometry - Navigator**System - Miscellaneous**

Positioning mode	FIX
Table position	H
Table position	0 mm
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P

Resolution - iPAT

PAT mode	GRAPPA
----------	--------

System - Miscellaneous

Transversal	F >> H
Coil Combine Mode	Adaptive Combine
Save uncombined	Off
Matrix Optimization	Off
Coil Focus	Flat
AutoAlign	---
Coil Select Mode	Default

System - Adjustments

B0 Shim mode	Cardiac
Adjust with body coil	On
Confirm freq. adjustment	Off
Assume Dominant Fat	Off
Assume Silicone	Off
Adjustment Tolerance	Auto

System - Adjust Volume

! Position	Isocenter
! Orientation	Transversal
! Rotation	0.00 deg
! A >> P	100 mm
! R >> L	100 mm
! F >> H	25 mm
Reset	Off

System - Tx/Rx

Frequency 1H	63.673772 MHz
Correction factor	1
Gain	High
Img. Scale Cor.	1.000
Reset	Off
? Ref. amplitude 1H	0.000 V

Physio - Signal1

1st Signal/Mode	ECG/Trigger
Average cycle	947 ± 122 ms
Average cycle	No Signal ms
Captured cycle	947 ± 122 ms
Acquisition window	831 ms
Trigger pulse	1
Trigger delay	470 ms
TR	360.56 ms
Concatenations	1
Segments	72
Phases	1
Adaptive Triggering	Off

Physio - Cardiac

Tagging	None
Magn. preparation	Non-sel. IR T1map
TI	260 ms
Fat suppr.	None
Dark blood	Off
FoV read	360 mm
FoV phase	75.0 %
Phase resolution	75 %
Cine	Off
Trajectory	Cartesian
Dummy heartbeats	0
Motion Correction	Standard

Physio - PACE

Resp. control	Off
Concatenations	1

Sequence - Part 1

Introduction	Off
Dimension	2D
Reordering	Linear
Asymmetric echo	Weak
Contrasts	1
Optimization	Min. TE TR
Multi-slice mode	Sequential
Sequence type	Trufi
Bandwidth	1085 Hz/Px

Sequence - Part 2

Define	Shots
Shots per slice	1
Segments	72
Trufi delta freq.	0 Hz
RF pulse type	Fast
Gradient mode	Fast
Excitation	Slice-sel.
Flip angle mode	Constant
Cine	Off

Sequence - Special

Partition coeff map	On
ECV map	Off
Synth ECV map	On
T1 scout mode	Off
Error map	On
Synth PSIR	On
Periods in seconds	On
16 bit images	Off
T1 sampling scheme	Post Gd

Sequence - Assistant

Mode	Off
Allowed delay	0 s

\\CMIC\Physics\T2MES\T2MES_low_res\T2map 0-25-55 new SSFP

TA: 0:12 PM: FIX Voxel size: 1.9×1.9×8.0 mmPAT: 2 Rel. SNR: 1.00 : tti

Properties

Prio recon	Off
Load images to viewer	On
Inline movie	Off
Auto store images	On
Load images to stamp segments	Off
Load images to graphic segments	On
Auto open inline display	Off
Auto close inline display	Off
Start measurement without further preparation	On
Wait for user to start	Off
Start measurements	Single measurement

Routine

Slice group	1
Slices	1
Dist. factor	20 %
Position	Isocenter
Orientation	Transversal
Phase enc. dir.	A >> P
AutoAlign	---
Phase oversampling	0 %
FoV read	360 mm
FoV phase	75.0 %
Slice thickness	8.0 mm
TR	207.98 ms
TE	1.06 ms
Averages	1
Concatenations	1
Filter	Distortion Corr.(2D)
Coil elements	BO1-3;SP2-4

Contrast - Common

TR	207.98 ms
TE	1.06 ms
Magn. preparation	T2 prep. adiab.
T2 prep. duration 1	0 ms
T2 prep. duration 2	25 ms
T2 prep. duration 3	55 ms
Flip angle	70 deg
Fat suppr.	None
Wrap-up Magn.	None

Contrast - Dynamic

Averages	1
Averaging mode	Short term
Reconstruction	Magnitude
Measurements	1
Multiple series	Off

Resolution - Common

FoV read	360 mm
FoV phase	75.0 %
Slice thickness	8.0 mm
Base resolution	192
Phase resolution	83 %
Phase partial Fourier	7/8
Trajectory	Cartesian
Interpolation	Off

Resolution - iPAT

PAT mode	GRAPPA
Accel. factor PE	2
Ref. lines PE	32
Reference scan mode	GRE/separate

Resolution - Filter Image

Image Filter	Off
Distortion Corr.	On
Mode	2D
Unfiltered images	Off
Prescan Normalize	Off
Normalize	Off
B1 filter	Off

Resolution - Filter Rawdata

Raw filter	Off
Elliptical filter	Off
POCS	Off

Geometry - Common

Slice group	1
Slices	1
Dist. factor	20 %
Position	Isocenter
Orientation	Transversal
Phase enc. dir.	A >> P
FoV read	360 mm
FoV phase	75.0 %
Slice thickness	8.0 mm
TR	207.98 ms
Multi-slice mode	Sequential
Series	Interleaved
Concatenations	1

Geometry - AutoAlign

Slice group	1
Position	Isocenter
Orientation	Transversal
Phase enc. dir.	A >> P
AutoAlign	---
Initial Position	Isocenter
L	0.0 mm
P	0.0 mm
H	0.0 mm
Initial Rotation	0.00 deg
Initial Orientation	Transversal

Geometry - Saturation

Fat suppr.	None
Wrap-up Magn.	None
Special sat.	None

Geometry - Navigator**System - Miscellaneous**

Positioning mode	FIX
Table position	H
Table position	0 mm
MSMA	S - C - T
Sagittal	R >> L

System - Miscellaneous

Coronal	A >> P
Transversal	F >> H
Coil Combine Mode	Adaptive Combine
Save uncombined	Off
Matrix Optimization	Off
Coil Focus	Flat
AutoAlign	---
Coil Select Mode	Default

System - Adjustments

B0 Shim mode	Standard
Adjust with body coil	On
Confirm freq. adjustment	Off
Assume Dominant Fat	Off
Assume Silicone	Off
Adjustment Tolerance	Auto

System - Adjust Volume

! Position	Isocenter
! Orientation	Transversal
! Rotation	0.00 deg
! A >> P	100 mm
! R >> L	100 mm
! F >> H	25 mm
Reset	Off

System - Tx/Rx

Frequency 1H	63.673772 MHz
Correction factor	1
Gain	High
Img. Scale Cor.	1.000
Reset	Off
? Ref. amplitude 1H	0.000 V

Physio - Signal1

1st Signal/Mode	ECG/Trigger
Average cycle	947 ± 122 ms
Average cycle	No Signal ms
Captured cycle	947 ± 122 ms
Acquisition window	831 ms
Trigger pulse	1
Trigger delay	623 ms
TR	207.98 ms
Concatenations	1
Segments	60
Phases	1
Adaptive Triggering	Off

Physio - Cardiac

Tagging	None
Magn. preparation	T2 prep. adiab.
T2 prep. duration 1	0 ms
T2 prep. duration 2	25 ms
T2 prep. duration 3	55 ms
Fat suppr.	None
Dark blood	Off
FoV read	360 mm
FoV phase	75.0 %
Phase resolution	83 %
Cine	Off
Trajectory	Cartesian
Dummy heartbeats	0
Motion Correction	None

Physio - PACE

Resp. control	Off
Concatenations	1

Sequence - Part 1

Introduction	Off
Dimension	2D
Reordering	Linear
Asymmetric echo	Weak
Contrasts	1
Optimization	Min. TE TR
Multi-slice mode	Sequential
Sequence type	Trufi
Bandwidth	1184 Hz/Px

Sequence - Part 2

Define	Shots
Shots per slice	1
Segments	60
Trufi delta freq.	0 Hz
RF pulse type	Fast
Gradient mode	Fast
Excitation	Slice-sel.
Flip angle mode	Constant
Cine	Off

Sequence - Special

Periods in seconds	On
16 bit images	Off

Sequence - Assistant

Mode	Off
Allowed delay	0 s

\\CMIC\Physics\T2MES\T2MES_low_res\pre FB T1T2

TA: 0:48 PM: FIX Voxel size: 1.4×1.4×8.0 mmPAT: 2 Rel. SNR: 1.00 : tfi

Properties

Prio recon	Off
Load images to viewer	On
Inline movie	Off
Auto store images	On
Load images to stamp segments	Off
Load images to graphic segments	On
Auto open inline display	Off
Auto close inline display	Off
Start measurement without further preparation	On
Wait for user to start	Off
Start measurements	Single measurement

Routine

Slice group	1
Slices	1
Dist. factor	20 %
Position	Isocenter
Orientation	Transversal
Phase enc. dir.	A >> P
AutoAlign	---
Phase oversampling	0 %
FoV read	360 mm
FoV phase	75.0 %
Slice thickness	8.0 mm
TR	824.00 ms
TE	1.26 ms
Averages	3
Concatenations	1
Filter	Distortion Corr.(2D)
Coil elements	BO1-3;SP2-4

Contrast - Common

TR	824.00 ms
TE	1.26 ms
Magn. preparation	Non-sel. SR perf
T1	550 ms
T2 prep. duration 1	0 ms
T2 prep. duration 2	0 ms
T2 prep. duration 3	0 ms
T2 prep. duration 4	0 ms
T2 prep. duration 5	0 ms
T2 prep. duration 6	0 ms
T2 prep. duration 7	0 ms
T2 prep. duration 8	55 ms
T2 prep. duration 9	55 ms
Flip angle	100 deg
Fat suppr.	None
Wrap-up Magn.	None

Contrast - Dynamic

Averages	3
Averaging mode	Short term
Reconstruction	Magnitude
Measurements	1
Multiple series	Off

Resolution - Common

FoV read	360 mm
FoV phase	75.0 %

Resolution - Common

Slice thickness	8.0 mm
Base resolution	256
Phase resolution	75 %
Phase partial Fourier	7/8
Trajectory	Cartesian
Interpolation	Off

Resolution - iPAT

PAT mode	GRAPPA
Accel. factor PE	2
Ref. lines PE	36
Reference scan mode	GRE/separate

Resolution - Filter Image

Image Filter	Off
Distortion Corr.	On
Mode	2D
Unfiltered images	Off
Prescan Normalize	Off
Normalize	Off
B1 filter	Off

Resolution - Filter Rawdata

Raw filter	Off
Elliptical filter	Off
POCS	Off

Geometry - Common

Slice group	1
Slices	1
Dist. factor	20 %
Position	Isocenter
Orientation	Transversal
Phase enc. dir.	A >> P
FoV read	360 mm
FoV phase	75.0 %
Slice thickness	8.0 mm
TR	824.00 ms
Multi-slice mode	Sequential
Series	Interleaved
Concatenations	1

Geometry - AutoAlign

Slice group	1
Position	Isocenter
Orientation	Transversal
Phase enc. dir.	A >> P
AutoAlign	---
Initial Position	Isocenter
L	0.0 mm
P	0.0 mm
H	0.0 mm
Initial Rotation	0.00 deg
Initial Orientation	Transversal

Geometry - Saturation

Fat suppr.	None
Wrap-up Magn.	None
Special sat.	None

Geometry - Navigator**System - Miscellaneous**

Positioning mode	FIX
Table position	H
Table position	0 mm
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Coil Combine Mode	Adaptive Combine
Save uncombined	Off
Matrix Optimization	Off
Coil Focus	Flat
AutoAlign	---
Coil Select Mode	Default

System - Adjustments

B0 Shim mode	Cardiac
Adjust with body coil	On
Confirm freq. adjustment	Off
Assume Dominant Fat	Off
Assume Silicone	Off
Adjustment Tolerance	Auto

System - Adjust Volume

! Position	Isocenter
! Orientation	Transversal
! Rotation	0.00 deg
! A >> P	100 mm
! R >> L	100 mm
! F >> H	50 mm
Reset	Off

System - Tx/Rx

Frequency 1H	63.673772 MHz
Correction factor	1
Gain	High
Img. Scale Cor.	1.000
Reset	Off
? Ref. amplitude 1H	0.000 V

Physio - Signal1

1st Signal/Mode	ECG/Trigger
Average cycle	947 ± 122 ms
Average cycle	No Signal ms
Captured cycle	947 ± 122 ms
Acquisition window	830 ms
Trigger pulse	1
Trigger delay	0 ms
TR	824.00 ms
Concatenations	1
Segments	72
Phases	1
Adaptive Triggering	Off

Physio - Cardiac

Tagging	None
Magn. preparation	Non-sel. SR perf
T1	550 ms
T2 prep. duration 1	0 ms
T2 prep. duration 2	0 ms
T2 prep. duration 3	0 ms
T2 prep. duration 4	0 ms
T2 prep. duration 5	0 ms

Physio - Cardiac

T2 prep. duration 6	0 ms
T2 prep. duration 7	0 ms
T2 prep. duration 8	55 ms
T2 prep. duration 9	55 ms
Fat suppr.	None
Dark blood	Off
FoV read	360 mm
FoV phase	75.0 %
Phase resolution	75 %
Cine	Off
Trajectory	Cartesian
Dummy heartbeats	0
Motion Correction	Standard

Physio - PACE

Resp. control	Off
Concatenations	1

Inline - Common

Subtract	Off
Measurements	1
StdDev	Off
Motion Correction	Standard
Save original images	On

Inline - Cardiac

Inline Evaluation	T1 map
Magn. preparation	Non-sel. SR perf
Num. of preps	9
Sampling duration 1	1 beats
Sampling duration 2	1 beats
Sampling duration 3	1 beats
Sampling duration 4	1 beats
Sampling duration 5	1 beats
Sampling duration 6	1 beats
Sampling duration 7	1 beats
Sampling duration 8	1 beats
Sampling duration 9	1 beats
T2 prep. duration 1	0 ms
T2 prep. duration 2	0 ms
T2 prep. duration 3	0 ms
T2 prep. duration 4	0 ms
T2 prep. duration 5	0 ms
T2 prep. duration 6	0 ms
T2 prep. duration 7	0 ms
T2 prep. duration 8	55 ms
T2 prep. duration 9	55 ms
Contrasts	1
TE	1.26 ms
TR	824.00 ms
Recovery duration 1	1 beats
Recovery duration 2	1 beats
Recovery duration 3	0 beats
Recovery duration 4	0 beats
Recovery duration 5	0 beats
Recovery duration 6	0 beats
Recovery duration 7	0 beats
Recovery duration 8	1 beats
Recovery duration 9	1 beats
Motion Correction	Standard
Save original images	On

Inline - MIP

MIP-Sag	Off
---------	-----

Inline - MIP

MIP-Cor	Off
MIP-Tra	Off
MIP-Time	Off
Save original images	On

Inline - Composing

Inline Composing	Off
Distortion Corr.	On
Mode	2D
Unfiltered images	Off

Sequence - Part 1

Introduction	Off
Dimension	2D
Reordering	Linear
Asymmetric echo	Weak
Contrasts	1
Optimization	Min. TE
Multi-slice mode	Sequential
Echo spacing	2.9 ms
Sequence type	Trufi
Bandwidth	1085 Hz/Px

Sequence - Part 2

Define	Shots
Shots per slice	1
Segments	72
Trufi delta freq.	0 Hz
RF pulse type	Fast
Gradient mode	Fast
Excitation	Slice-sel.
Flip angle mode	Variable
Cine	Off

Sequence - Special

Partition coeff map	Off
ECV map	Off
Synth ECV map	Off
T1 array 0	550 ms
T1 array 1	550 ms
T1 array 2	550 ms
T1 array 3	550 ms
T1 array 4	550 ms
T1 array 5	550 ms
T1 array 6	550 ms
T1 array 7	550 ms
T1 array 8	550 ms
High-contrast acquisition	On
Error map	Off
SPAIR delay	5 ms
T1 sampling scheme	Custom

Sequence - Assistant

Mode	Off
Allowed delay	10 s

\\CMIC\Physics\T2MES\T2MES_low_res\pre FB T1T2 high HR

TA: 1:05 PM: FIX Voxel size: 1.4×1.4×8.0 mmPAT: 2 Rel. SNR: 1.00 : tfi

Properties

Prio recon	Off
Load images to viewer	On
Inline movie	Off
Auto store images	On
Load images to stamp segments	Off
Load images to graphic segments	On
Auto open inline display	Off
Auto close inline display	Off
Start measurement without further preparation	On
Wait for user to start	Off
Start measurements	Single measurement

Routine

Slice group	1
Slices	1
Dist. factor	20 %
Position	Isocenter
Orientation	Transversal
Phase enc. dir.	A >> P
AutoAlign	---
Phase oversampling	0 %
FoV read	360 mm
FoV phase	75.0 %
Slice thickness	8.0 mm
TR	824.00 ms
TE	1.26 ms
Averages	3
Concatenations	1
Filter	Distortion Corr.(2D)
Coil elements	BO1-3;SP2-4

Contrast - Common

TR	824.00 ms
TE	1.26 ms
Magn. preparation	Non-sel. SR perf
T1	350 ms
T2 prep. duration 1	0 ms
T2 prep. duration 2	0 ms
T2 prep. duration 3	0 ms
T2 prep. duration 4	0 ms
T2 prep. duration 5	0 ms
T2 prep. duration 6	0 ms
T2 prep. duration 7	0 ms
T2 prep. duration 8	55 ms
T2 prep. duration 9	55 ms
Flip angle	100 deg
Fat suppr.	None
Wrap-up Magn.	None

Contrast - Dynamic

Averages	3
Averaging mode	Short term
Reconstruction	Magnitude
Measurements	1
Multiple series	Off

Resolution - Common

FoV read	360 mm
FoV phase	75.0 %

Resolution - Common

Slice thickness	8.0 mm
Base resolution	256
Phase resolution	75 %
Phase partial Fourier	7/8
Trajectory	Cartesian
Interpolation	Off

Resolution - iPAT

PAT mode	GRAPPA
Accel. factor PE	2
Ref. lines PE	36
Reference scan mode	GRE/separate

Resolution - Filter Image

Image Filter	Off
Distortion Corr.	On
Mode	2D
Unfiltered images	Off
Prescan Normalize	Off
Normalize	Off
B1 filter	Off

Resolution - Filter Rawdata

Raw filter	Off
Elliptical filter	Off
POCS	Off

Geometry - Common

Slice group	1
Slices	1
Dist. factor	20 %
Position	Isocenter
Orientation	Transversal
Phase enc. dir.	A >> P
FoV read	360 mm
FoV phase	75.0 %
Slice thickness	8.0 mm
TR	824.00 ms
Multi-slice mode	Sequential
Series	Interleaved
Concatenations	1

Geometry - AutoAlign

Slice group	1
Position	Isocenter
Orientation	Transversal
Phase enc. dir.	A >> P
AutoAlign	---
Initial Position	Isocenter
L	0.0 mm
P	0.0 mm
H	0.0 mm
Initial Rotation	0.00 deg
Initial Orientation	Transversal

Geometry - Saturation

Fat suppr.	None
Wrap-up Magn.	None
Special sat.	None

Geometry - Navigator**System - Miscellaneous**

Positioning mode	FIX
Table position	H
Table position	0 mm
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Coil Combine Mode	Adaptive Combine
Save uncombined	Off
Matrix Optimization	Off
Coil Focus	Flat
AutoAlign	---
Coil Select Mode	Default

System - Adjustments

B0 Shim mode	Cardiac
Adjust with body coil	On
Confirm freq. adjustment	Off
Assume Dominant Fat	Off
Assume Silicone	Off
Adjustment Tolerance	Auto

System - Adjust Volume

! Position	Isocenter
! Orientation	Transversal
! Rotation	0.00 deg
! A >> P	150 mm
! R >> L	150 mm
! F >> H	50 mm
Reset	Off

System - Tx/Rx

Frequency 1H	63.673772 MHz
Correction factor	1
Gain	High
Img. Scale Cor.	1.000
Reset	Off
? Ref. amplitude 1H	0.000 V

Physio - Signal1

1st Signal/Mode	ECG/Trigger
Average cycle	947 ± 122 ms
Average cycle	No Signal ms
Captured cycle	947 ± 122 ms
Acquisition window	830 ms
Trigger pulse	1
Trigger delay	0 ms
TR	824.00 ms
Concatenations	1
Segments	72
Phases	1
Adaptive Triggering	Off

Physio - Cardiac

Tagging	None
Magn. preparation	Non-sel. SR perf
T1	350 ms
T2 prep. duration 1	0 ms
T2 prep. duration 2	0 ms
T2 prep. duration 3	0 ms
T2 prep. duration 4	0 ms
T2 prep. duration 5	0 ms

Physio - Cardiac

T2 prep. duration 6	0 ms
T2 prep. duration 7	0 ms
T2 prep. duration 8	55 ms
T2 prep. duration 9	55 ms
Fat suppr.	None
Dark blood	Off
FoV read	360 mm
FoV phase	75.0 %
Phase resolution	75 %
Cine	Off
Trajectory	Cartesian
Dummy heartbeats	0
Motion Correction	Standard

Physio - PACE

Resp. control	Off
Concatenations	1

Inline - Common

Subtract	Off
Measurements	1
StdDev	Off
Motion Correction	Standard
Save original images	On

Inline - Cardiac

Inline Evaluation	T1 map
Magn. preparation	Non-sel. SR perf
Num. of preps	9
Sampling duration 1	1 beats
Sampling duration 2	1 beats
Sampling duration 3	1 beats
Sampling duration 4	1 beats
Sampling duration 5	1 beats
Sampling duration 6	1 beats
Sampling duration 7	1 beats
Sampling duration 8	1 beats
Sampling duration 9	1 beats
T2 prep. duration 1	0 ms
T2 prep. duration 2	0 ms
T2 prep. duration 3	0 ms
T2 prep. duration 4	0 ms
T2 prep. duration 5	0 ms
T2 prep. duration 6	0 ms
T2 prep. duration 7	0 ms
T2 prep. duration 8	55 ms
T2 prep. duration 9	55 ms
Contrasts	1
TE	1.26 ms
TR	824.00 ms
Recovery duration 1	2 beats
Recovery duration 2	2 beats
Recovery duration 3	0 beats
Recovery duration 4	0 beats
Recovery duration 5	0 beats
Recovery duration 6	0 beats
Recovery duration 7	0 beats
Recovery duration 8	3 beats
Recovery duration 9	3 beats
Motion Correction	Standard
Save original images	On

Inline - MIP

MIP-Sag	Off
---------	-----

Inline - MIP

MIP-Cor	Off
MIP-Tra	Off
MIP-Time	Off
Save original images	On

Inline - Composing

Inline Composing	Off
Distortion Corr.	On
Mode	2D
Unfiltered images	Off

Sequence - Part 1

Introduction	Off
Dimension	2D
Reordering	Linear
Asymmetric echo	Weak
Contrasts	1
Optimization	Min. TE
Multi-slice mode	Sequential
Echo spacing	2.9 ms
Sequence type	Trufi
Bandwidth	1085 Hz/Px

Sequence - Part 2

Define	Shots
Shots per slice	1
Segments	72
Trufi delta freq.	0 Hz
RF pulse type	Fast
Gradient mode	Fast
Excitation	Slice-sel.
Flip angle mode	Variable
Cine	Off

Sequence - Special

Partition coeff map	Off
ECV map	Off
Synth ECV map	Off
T1 array 0	350 ms
T1 array 1	350 ms
T1 array 2	350 ms
T1 array 3	350 ms
T1 array 4	350 ms
T1 array 5	350 ms
T1 array 6	350 ms
T1 array 7	350 ms
T1 array 8	350 ms
High-contrast acquisition	On
Error map	Off
SPAIR delay	5 ms
T1 sampling scheme	Fixed TS

Sequence - Assistant

Mode	Off
Allowed delay	10 s

\\CMIC\Physics\T2MES\T2MES_low_res\post FB T1T2

TA: 0:43 PM: FIX Voxel size: 1.4×1.4×8.0 mmPAT: 2 Rel. SNR: 1.00 : tfi

Properties

Prio recon	Off
Load images to viewer	On
Inline movie	Off
Auto store images	On
Load images to stamp segments	Off
Load images to graphic segments	On
Auto open inline display	Off
Auto close inline display	Off
Start measurement without further preparation	On
Wait for user to start	Off
Start measurements	Single measurement

Routine

Slice group	1
Slices	1
Dist. factor	20 %
Position	Isocenter
Orientation	Transversal
Phase enc. dir.	A >> P
AutoAlign	---
Phase oversampling	0 %
FoV read	360 mm
FoV phase	75.0 %
Slice thickness	8.0 mm
TR	904.00 ms
TE	1.26 ms
Averages	3
Concatenations	1
Filter	Distortion Corr.(2D)
Coil elements	BO1-3;SP2-4

Contrast - Common

TR	904.00 ms
TE	1.26 ms
Magn. preparation	Non-sel. SR perf
T1	300 ms
T2 prep. duration 1	0 ms
T2 prep. duration 2	0 ms
T2 prep. duration 3	0 ms
T2 prep. duration 4	0 ms
T2 prep. duration 5	0 ms
T2 prep. duration 6	0 ms
T2 prep. duration 7	0 ms
T2 prep. duration 8	55 ms
T2 prep. duration 9	55 ms
Flip angle	100 deg
Fat suppr.	None
Wrap-up Magn.	None

Contrast - Dynamic

Averages	3
Averaging mode	Short term
Reconstruction	Magnitude
Measurements	1
Multiple series	Off

Resolution - Common

FoV read	360 mm
FoV phase	75.0 %

Resolution - Common

Slice thickness	8.0 mm
Base resolution	256
Phase resolution	75 %
Phase partial Fourier	7/8
Trajectory	Cartesian
Interpolation	Off

Resolution - iPAT

PAT mode	GRAPPA
Accel. factor PE	2
Ref. lines PE	36
Reference scan mode	GRE/separate

Resolution - Filter Image

Image Filter	Off
Distortion Corr.	On
Mode	2D
Unfiltered images	Off
Prescan Normalize	Off
Normalize	Off
B1 filter	Off

Resolution - Filter Rawdata

Raw filter	Off
Elliptical filter	Off
POCS	Off

Geometry - Common

Slice group	1
Slices	1
Dist. factor	20 %
Position	Isocenter
Orientation	Transversal
Phase enc. dir.	A >> P
FoV read	360 mm
FoV phase	75.0 %
Slice thickness	8.0 mm
TR	904.00 ms
Multi-slice mode	Sequential
Series	Interleaved
Concatenations	1

Geometry - AutoAlign

Slice group	1
Position	Isocenter
Orientation	Transversal
Phase enc. dir.	A >> P
AutoAlign	---
Initial Position	Isocenter
L	0.0 mm
P	0.0 mm
H	0.0 mm
Initial Rotation	0.00 deg
Initial Orientation	Transversal

Geometry - Saturation

Fat suppr.	None
Wrap-up Magn.	None
Special sat.	None

Geometry - Navigator**System - Miscellaneous**

Positioning mode	FIX
Table position	H
Table position	0 mm
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Coil Combine Mode	Adaptive Combine
Save uncombined	Off
Matrix Optimization	Off
Coil Focus	Flat
AutoAlign	---
Coil Select Mode	Default

System - Adjustments

B0 Shim mode	Cardiac
Adjust with body coil	On
Confirm freq. adjustment	Off
Assume Dominant Fat	Off
Assume Silicone	Off
Adjustment Tolerance	Auto

System - Adjust Volume

! Position	Isocenter
! Orientation	Transversal
! Rotation	0.00 deg
! A >> P	150 mm
! R >> L	150 mm
! F >> H	50 mm
Reset	Off

System - Tx/Rx

Frequency 1H	63.673772 MHz
Correction factor	1
Gain	High
Img. Scale Cor.	1.000
Reset	Off
? Ref. amplitude 1H	0.000 V

Physio - Signal1

1st Signal/Mode	ECG/Trigger
Average cycle	947 ± 122 ms
Average cycle	No Signal ms
Captured cycle	947 ± 122 ms
Acquisition window	930 ms
Trigger pulse	1
Trigger delay	0 ms
TR	904.00 ms
Concatenations	1
Segments	72
Phases	1
Adaptive Triggering	Off

Physio - Cardiac

Tagging	None
Magn. preparation	Non-sel. SR perf
T1	300 ms
T2 prep. duration 1	0 ms
T2 prep. duration 2	0 ms
T2 prep. duration 3	0 ms
T2 prep. duration 4	0 ms
T2 prep. duration 5	0 ms

Physio - Cardiac

T2 prep. duration 6	0 ms
T2 prep. duration 7	0 ms
T2 prep. duration 8	55 ms
T2 prep. duration 9	55 ms
Fat suppr.	None
Dark blood	Off
FoV read	360 mm
FoV phase	75.0 %
Phase resolution	75 %
Cine	Off
Trajectory	Cartesian
Dummy heartbeats	0
Motion Correction	Standard

Physio - PACE

Resp. control	Off
Concatenations	1

Inline - Common

Subtract	Off
Measurements	1
StdDev	Off
Motion Correction	Standard
Save original images	On

Inline - Cardiac

Inline Evaluation	T1 map
Magn. preparation	Non-sel. SR perf
Num. of preps	9
Sampling duration 1	1 beats
Sampling duration 2	1 beats
Sampling duration 3	1 beats
Sampling duration 4	1 beats
Sampling duration 5	1 beats
Sampling duration 6	1 beats
Sampling duration 7	1 beats
Sampling duration 8	1 beats
Sampling duration 9	1 beats
T2 prep. duration 1	0 ms
T2 prep. duration 2	0 ms
T2 prep. duration 3	0 ms
T2 prep. duration 4	0 ms
T2 prep. duration 5	0 ms
T2 prep. duration 6	0 ms
T2 prep. duration 7	0 ms
T2 prep. duration 8	55 ms
T2 prep. duration 9	55 ms
Contrasts	1
TE	1.26 ms
TR	904.00 ms
Recovery duration 1	0 beats
Recovery duration 2	0 beats
Recovery duration 3	0 beats
Recovery duration 4	0 beats
Recovery duration 5	0 beats
Recovery duration 6	0 beats
Recovery duration 7	0 beats
Recovery duration 8	1 beats
Recovery duration 9	1 beats
Motion Correction	Standard
Save original images	On

Inline - MIP

MIP-Sag	Off
---------	-----

Inline - MIP

MIP-Cor	Off
MIP-Tra	Off
MIP-Time	Off
Save original images	On

Inline - Composing

Inline Composing	Off
Distortion Corr.	On
Mode	2D
Unfiltered images	Off

Sequence - Part 1

Introduction	Off
Dimension	2D
Reordering	Linear
Asymmetric echo	Weak
Contrasts	1
Optimization	Min. TE
Multi-slice mode	Sequential
Echo spacing	2.9 ms
Sequence type	Trufi
Bandwidth	1085 Hz/Px

Sequence - Part 2

Define	Shots
Shots per slice	1
Segments	72
Trufi delta freq.	0 Hz
RF pulse type	Fast
Gradient mode	Fast
Excitation	Slice-sel.
Flip angle mode	Variable
Cine	Off

Sequence - Special

Partition coeff map	Off
ECV map	Off
Synth ECV map	Off
T1 array 0	300 ms
T1 array 1	300 ms
T1 array 2	300 ms
T1 array 3	300 ms
T1 array 4	300 ms
T1 array 5	300 ms
T1 array 6	300 ms
T1 array 7	300 ms
T1 array 8	300 ms
High-contrast acquisition	On
Error map	Off
SPAIR delay	5 ms
T1 sampling scheme	Fixed TS

Sequence - Assistant

Mode	Off
Allowed delay	10 s

\\CMIC\Physics\T2MES\T2MES_low_res\pre_MOLLI_5s(3s)3s_256

TA: 0:11 PM: FIX Voxel size: 1.4×1.4×8.0 mmPAT: 2 Rel. SNR: 1.00 : tfi

Properties

Prio recon	Off
Load images to viewer	On
Inline movie	Off
Auto store images	On
Load images to stamp segments	Off
Load images to graphic segments	On
Auto open inline display	Off
Auto close inline display	Off
Start measurement without further preparation	On
Wait for user to start	Off
Start measurements	Single measurement

Routine

Slice group	1
Slices	1
Dist. factor	20 %
Position	Isocenter
Orientation	Transversal
Phase enc. dir.	A >> P
AutoAlign	---
Phase oversampling	0 %
FoV read	360 mm
FoV phase	75.0 %
Slice thickness	8.0 mm
TR	280.56 ms
TE	1.12 ms
Averages	1
Concatenations	1
Filter	Distortion Corr.(2D)
Coil elements	BO1-3;SP2-4

Contrast - Common

TR	280.56 ms
TE	1.12 ms
Magn. preparation	Non-sel. IR T1map
TI	180 ms
Flip angle	35 deg
Fat suppr.	None
Wrap-up Magn.	None

Contrast - Dynamic

Averages	1
Averaging mode	Short term
Reconstruction	Magnitude
Measurements	1
Multiple series	Off

Resolution - Common

FoV read	360 mm
FoV phase	75.0 %
Slice thickness	8.0 mm
Base resolution	256
Phase resolution	75 %
Phase partial Fourier	7/8
Trajectory	Cartesian
Interpolation	Off

Resolution - iPAT

PAT mode	GRAPPA
----------	--------

Resolution - iPAT

Accel. factor PE	2
Ref. lines PE	36
Reference scan mode	GRE/separate

Resolution - Filter Image

Image Filter	Off
Distortion Corr.	On
Mode	2D
Unfiltered images	Off
Prescan Normalize	Off
Normalize	Off
B1 filter	Off

Resolution - Filter Rawdata

Raw filter	Off
Elliptical filter	Off
POCS	Off

Geometry - Common

Slice group	1
Slices	1
Dist. factor	20 %
Position	Isocenter
Orientation	Transversal
Phase enc. dir.	A >> P
FoV read	360 mm
FoV phase	75.0 %
Slice thickness	8.0 mm
TR	280.56 ms
Multi-slice mode	Sequential
Series	Interleaved
Concatenations	1

Geometry - AutoAlign

Slice group	1
Position	Isocenter
Orientation	Transversal
Phase enc. dir.	A >> P
AutoAlign	---
Initial Position	Isocenter
L	0.0 mm
P	0.0 mm
H	0.0 mm
Initial Rotation	0.00 deg
Initial Orientation	Transversal

Geometry - Saturation

Fat suppr.	None
Wrap-up Magn.	None
Special sat.	None

Geometry - Navigator**System - Miscellaneous**

Positioning mode	FIX
Table position	H
Table position	0 mm
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P

System - Miscellaneous

Transversal	F >> H
Coil Combine Mode	Adaptive Combine
Save uncombined	Off
Matrix Optimization	Off
Coil Focus	Flat
AutoAlign	---
Coil Select Mode	Default

System - Adjustments

B0 Shim mode	Standard
Adjust with body coil	On
Confirm freq. adjustment	Off
Assume Dominant Fat	Off
Assume Silicone	Off
Adjustment Tolerance	Auto

System - Adjust Volume

! Position	Isocenter
! Orientation	Transversal
! Rotation	0.00 deg
! A >> P	100 mm
! R >> L	100 mm
! F >> H	25 mm
Reset	Off

System - Tx/Rx

Frequency 1H	63.673772 MHz
Correction factor	1
Gain	High
Img. Scale Cor.	1.000
Reset	Off
? Ref. amplitude 1H	0.000 V

Physio - Signal1

1st Signal/Mode	ECG/Trigger
Average cycle	947 ± 122 ms
Average cycle	No Signal ms
Captured cycle	947 ± 122 ms
Acquisition window	700 ms
Trigger pulse	1
Trigger delay	418 ms
TR	280.56 ms
Concatenations	1
Segments	72
Phases	1
Adaptive Triggering	Off

Physio - Cardiac

Tagging	None
Magn. preparation	Non-sel. IR T1map
TI	180 ms
Fat suppr.	None
Dark blood	Off
FoV read	360 mm
FoV phase	75.0 %
Phase resolution	75 %
Cine	Off
Trajectory	Cartesian
Dummy heartbeats	0
Motion Correction	Standard

Physio - PACE

Resp. control	Off
Concatenations	1

Sequence - Part 1

Introduction	Off
Dimension	2D
Reordering	Linear
Asymmetric echo	Weak
Contrasts	1
Optimization	Min. TE TR
Multi-slice mode	Sequential
Sequence type	Trufi
Bandwidth	1085 Hz/Px

Sequence - Part 2

Define	Shots
Shots per slice	1
Segments	72
Trufi delta freq.	0 Hz
RF pulse type	Fast
Gradient mode	Fast
Excitation	Slice-sel.
Flip angle mode	Constant
Cine	Off

Sequence - Special

Partition coeff map	Off
ECV map	Off
Synth ECV map	Off
T1 scout mode	Off
Error map	On
Synth PSIR	Off
Periods in seconds	On
16 bit images	Off
T1 sampling scheme	Native

Sequence - Assistant

Mode	Off
Allowed delay	0 s

\\CMIC\Physics\T2MES\T2MES_low_res\post_MOLLI_4s(1s)3s(1s)2s_256

TA: 0:12 PM: FIX Voxel size: 1.4×1.4×8.0 mmPAT: 2 Rel. SNR: 1.00 : tfi

Properties

Prio recon	Off
Load images to viewer	On
Inline movie	Off
Auto store images	On
Load images to stamp segments	Off
Load images to graphic segments	On
Auto open inline display	Off
Auto close inline display	Off
Start measurement without further preparation	On
Wait for user to start	Off
Start measurements	Single measurement

Resolution - iPAT

Accel. factor PE	2
Ref. lines PE	36
Reference scan mode	GRE/separate

Resolution - Filter Image

Image Filter	Off
Distortion Corr.	On
Mode	2D
Unfiltered images	Off
Prescan Normalize	Off
Normalize	Off
B1 filter	Off

Routine

Slice group	1
Slices	1
Dist. factor	20 %
Position	Isocenter
Orientation	Transversal
Phase enc. dir.	A >> P
AutoAlign	---
Phase oversampling	0 %
FoV read	360 mm
FoV phase	75.0 %
Slice thickness	8.0 mm
TR	360.56 ms
TE	1.12 ms
Averages	1
Concatenations	1
Filter	Distortion Corr.(2D)
Coil elements	BO1-3;SP2-4

Resolution - Filter Rawdata

Raw filter	Off
Elliptical filter	Off
POCS	Off

Geometry - Common

Slice group	1
Slices	1
Dist. factor	20 %
Position	Isocenter
Orientation	Transversal
Phase enc. dir.	A >> P
FoV read	360 mm
FoV phase	75.0 %
Slice thickness	8.0 mm
TR	360.56 ms
Multi-slice mode	Sequential
Series	Interleaved
Concatenations	1

Contrast - Common

TR	360.56 ms
TE	1.12 ms
Magn. preparation	Non-sel. IR T1map
T1	260 ms
Flip angle	35 deg
Fat suppr.	None
Wrap-up Magn.	None

Geometry - AutoAlign

Slice group	1
Position	Isocenter
Orientation	Transversal
Phase enc. dir.	A >> P
AutoAlign	---
Initial Position	Isocenter
L	0.0 mm
P	0.0 mm
H	0.0 mm
Initial Rotation	0.00 deg
Initial Orientation	Transversal

Contrast - Dynamic

Averages	1
Averaging mode	Short term
Reconstruction	Magnitude
Measurements	1
Multiple series	Off

Geometry - Saturation

Fat suppr.	None
Wrap-up Magn.	None
Special sat.	None

Resolution - Common

FoV read	360 mm
FoV phase	75.0 %
Slice thickness	8.0 mm
Base resolution	256
Phase resolution	75 %
Phase partial Fourier	7/8
Trajectory	Cartesian
Interpolation	Off

Geometry - Navigator**System - Miscellaneous**

Positioning mode	FIX
Table position	H
Table position	0 mm
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P

Resolution - iPAT

PAT mode	GRAPPA
----------	--------

System - Miscellaneous

Transversal	F >> H
Coil Combine Mode	Adaptive Combine
Save uncombined	Off
Matrix Optimization	Off
Coil Focus	Flat
AutoAlign	---
Coil Select Mode	Default

System - Adjustments

B0 Shim mode	Standard
Adjust with body coil	On
Confirm freq. adjustment	Off
Assume Dominant Fat	Off
Assume Silicone	Off
Adjustment Tolerance	Auto

System - Adjust Volume

! Position	Isocenter
! Orientation	Transversal
! Rotation	0.00 deg
! A >> P	100 mm
! R >> L	100 mm
! F >> H	25 mm
Reset	Off

System - Tx/Rx

Frequency 1H	63.673772 MHz
Correction factor	1
Gain	High
Img. Scale Cor.	1.000
Reset	Off
? Ref. amplitude 1H	0.000 V

Physio - Signal1

1st Signal/Mode	ECG/Trigger
Average cycle	947 ± 122 ms
Average cycle	No Signal ms
Captured cycle	947 ± 122 ms
Acquisition window	831 ms
Trigger pulse	1
Trigger delay	470 ms
TR	360.56 ms
Concatenations	1
Segments	72
Phases	1
Adaptive Triggering	Off

Physio - Cardiac

Tagging	None
Magn. preparation	Non-sel. IR T1map
TI	260 ms
Fat suppr.	None
Dark blood	Off
FoV read	360 mm
FoV phase	75.0 %
Phase resolution	75 %
Cine	Off
Trajectory	Cartesian
Dummy heartbeats	0
Motion Correction	Standard

Physio - PACE

Resp. control	Off
Concatenations	1

Sequence - Part 1

Introduction	Off
Dimension	2D
Reordering	Linear
Asymmetric echo	Weak
Contrasts	1
Optimization	Min. TE TR
Multi-slice mode	Sequential
Sequence type	Trufi
Bandwidth	1085 Hz/Px

Sequence - Part 2

Define	Shots
Shots per slice	1
Segments	72
Trufi delta freq.	0 Hz
RF pulse type	Fast
Gradient mode	Fast
Excitation	Slice-sel.
Flip angle mode	Constant
Cine	Off

Sequence - Special

Partition coeff map	On
ECV map	Off
Synth ECV map	On
T1 scout mode	Off
Error map	On
Synth PSIR	On
Periods in seconds	On
16 bit images	Off
T1 sampling scheme	Post Gd

Sequence - Assistant

Mode	Off
Allowed delay	0 s

\\CMIC\Physics\T2MES\T2MES_low_res\T2map 0-25-55 new SSFP

TA: 0:12 PM: FIX Voxel size: 1.9×1.9×8.0 mmPAT: 2 Rel. SNR: 1.00 : tti

Properties

Prio recon	Off
Load images to viewer	On
Inline movie	Off
Auto store images	On
Load images to stamp segments	Off
Load images to graphic segments	On
Auto open inline display	Off
Auto close inline display	Off
Start measurement without further preparation	On
Wait for user to start	Off
Start measurements	Single measurement

Routine

Slice group	1
Slices	1
Dist. factor	20 %
Position	Isocenter
Orientation	Transversal
Phase enc. dir.	A >> P
AutoAlign	---
Phase oversampling	0 %
FoV read	360 mm
FoV phase	75.0 %
Slice thickness	8.0 mm
TR	207.98 ms
TE	1.06 ms
Averages	1
Concatenations	1
Filter	Distortion Corr.(2D)
Coil elements	BO1-3;SP2-4

Contrast - Common

TR	207.98 ms
TE	1.06 ms
Magn. preparation	T2 prep. adiab.
T2 prep. duration 1	0 ms
T2 prep. duration 2	25 ms
T2 prep. duration 3	55 ms
Flip angle	70 deg
Fat suppr.	None
Wrap-up Magn.	None

Contrast - Dynamic

Averages	1
Averaging mode	Short term
Reconstruction	Magnitude
Measurements	1
Multiple series	Off

Resolution - Common

FoV read	360 mm
FoV phase	75.0 %
Slice thickness	8.0 mm
Base resolution	192
Phase resolution	83 %
Phase partial Fourier	7/8
Trajectory	Cartesian
Interpolation	Off

Resolution - iPAT

PAT mode	GRAPPA
Accel. factor PE	2
Ref. lines PE	32
Reference scan mode	GRE/separate

Resolution - Filter Image

Image Filter	Off
Distortion Corr.	On
Mode	2D
Unfiltered images	Off
Prescan Normalize	Off
Normalize	Off
B1 filter	Off

Resolution - Filter Rawdata

Raw filter	Off
Elliptical filter	Off
POCS	Off

Geometry - Common

Slice group	1
Slices	1
Dist. factor	20 %
Position	Isocenter
Orientation	Transversal
Phase enc. dir.	A >> P
FoV read	360 mm
FoV phase	75.0 %
Slice thickness	8.0 mm
TR	207.98 ms
Multi-slice mode	Sequential
Series	Interleaved
Concatenations	1

Geometry - AutoAlign

Slice group	1
Position	Isocenter
Orientation	Transversal
Phase enc. dir.	A >> P
AutoAlign	---
Initial Position	Isocenter
L	0.0 mm
P	0.0 mm
H	0.0 mm
Initial Rotation	0.00 deg
Initial Orientation	Transversal

Geometry - Saturation

Fat suppr.	None
Wrap-up Magn.	None
Special sat.	None

Geometry - Navigator**System - Miscellaneous**

Positioning mode	FIX
Table position	H
Table position	0 mm
MSMA	S - C - T
Sagittal	R >> L

System - Miscellaneous

Coronal	A >> P
Transversal	F >> H
Coil Combine Mode	Adaptive Combine
Save uncombined	Off
Matrix Optimization	Off
Coil Focus	Flat
AutoAlign	---
Coil Select Mode	Default

System - Adjustments

B0 Shim mode	Standard
Adjust with body coil	On
Confirm freq. adjustment	Off
Assume Dominant Fat	Off
Assume Silicone	Off
Adjustment Tolerance	Auto

System - Adjust Volume

! Position	Isocenter
! Orientation	Transversal
! Rotation	0.00 deg
! A >> P	100 mm
! R >> L	100 mm
! F >> H	25 mm
Reset	Off

System - Tx/Rx

Frequency 1H	63.673772 MHz
Correction factor	1
Gain	High
Img. Scale Cor.	1.000
Reset	Off
? Ref. amplitude 1H	0.000 V

Physio - Signal1

1st Signal/Mode	ECG/Trigger
Average cycle	947 ± 122 ms
Average cycle	No Signal ms
Captured cycle	947 ± 122 ms
Acquisition window	831 ms
Trigger pulse	1
Trigger delay	623 ms
TR	207.98 ms
Concatenations	1
Segments	60
Phases	1
Adaptive Triggering	Off

Physio - Cardiac

Tagging	None
Magn. preparation	T2 prep. adiab.
T2 prep. duration 1	0 ms
T2 prep. duration 2	25 ms
T2 prep. duration 3	55 ms
Fat suppr.	None
Dark blood	Off
FoV read	360 mm
FoV phase	75.0 %
Phase resolution	83 %
Cine	Off
Trajectory	Cartesian
Dummy heartbeats	0
Motion Correction	None

Physio - PACE

Resp. control	Off
Concatenations	1

Sequence - Part 1

Introduction	Off
Dimension	2D
Reordering	Linear
Asymmetric echo	Weak
Contrasts	1
Optimization	Min. TE TR
Multi-slice mode	Sequential
Sequence type	Trufi
Bandwidth	1184 Hz/Px

Sequence - Part 2

Define	Shots
Shots per slice	1
Segments	60
Trufi delta freq.	0 Hz
RF pulse type	Fast
Gradient mode	Fast
Excitation	Slice-sel.
Flip angle mode	Constant
Cine	Off

Sequence - Special

Periods in seconds	On
16 bit images	Off

Sequence - Assistant

Mode	Off
Allowed delay	0 s

\\CMIC\Physics\T2MES\T2MES_low_res\pre FB T1T2

TA: 0:48 PM: FIX Voxel size: 1.4×1.4×8.0 mmPAT: 2 Rel. SNR: 1.00 : tfi

Properties

Prio recon	Off
Load images to viewer	On
Inline movie	Off
Auto store images	On
Load images to stamp segments	Off
Load images to graphic segments	On
Auto open inline display	Off
Auto close inline display	Off
Start measurement without further preparation	On
Wait for user to start	Off
Start measurements	Single measurement

Routine

Slice group	1
Slices	1
Dist. factor	20 %
Position	Isocenter
Orientation	Transversal
Phase enc. dir.	A >> P
AutoAlign	---
Phase oversampling	0 %
FoV read	360 mm
FoV phase	75.0 %
Slice thickness	8.0 mm
TR	824.00 ms
TE	1.26 ms
Averages	3
Concatenations	1
Filter	Distortion Corr.(2D)
Coil elements	BO1-3;SP2-4

Contrast - Common

TR	824.00 ms
TE	1.26 ms
Magn. preparation	Non-sel. SR perf
T1	550 ms
T2 prep. duration 1	0 ms
T2 prep. duration 2	0 ms
T2 prep. duration 3	0 ms
T2 prep. duration 4	0 ms
T2 prep. duration 5	0 ms
T2 prep. duration 6	0 ms
T2 prep. duration 7	0 ms
T2 prep. duration 8	55 ms
T2 prep. duration 9	55 ms
Flip angle	100 deg
Fat suppr.	None
Wrap-up Magn.	None

Contrast - Dynamic

Averages	3
Averaging mode	Short term
Reconstruction	Magnitude
Measurements	1
Multiple series	Off

Resolution - Common

FoV read	360 mm
FoV phase	75.0 %

Resolution - Common

Slice thickness	8.0 mm
Base resolution	256
Phase resolution	75 %
Phase partial Fourier	7/8
Trajectory	Cartesian
Interpolation	Off

Resolution - iPAT

PAT mode	GRAPPA
Accel. factor PE	2
Ref. lines PE	36
Reference scan mode	GRE/separate

Resolution - Filter Image

Image Filter	Off
Distortion Corr.	On
Mode	2D
Unfiltered images	Off
Prescan Normalize	Off
Normalize	Off
B1 filter	Off

Resolution - Filter Rawdata

Raw filter	Off
Elliptical filter	Off
POCS	Off

Geometry - Common

Slice group	1
Slices	1
Dist. factor	20 %
Position	Isocenter
Orientation	Transversal
Phase enc. dir.	A >> P
FoV read	360 mm
FoV phase	75.0 %
Slice thickness	8.0 mm
TR	824.00 ms
Multi-slice mode	Sequential
Series	Interleaved
Concatenations	1

Geometry - AutoAlign

Slice group	1
Position	Isocenter
Orientation	Transversal
Phase enc. dir.	A >> P
AutoAlign	---
Initial Position	Isocenter
L	0.0 mm
P	0.0 mm
H	0.0 mm
Initial Rotation	0.00 deg
Initial Orientation	Transversal

Geometry - Saturation

Fat suppr.	None
Wrap-up Magn.	None
Special sat.	None

Geometry - Navigator**System - Miscellaneous**

Positioning mode	FIX
Table position	H
Table position	0 mm
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Coil Combine Mode	Adaptive Combine
Save uncombined	Off
Matrix Optimization	Off
Coil Focus	Flat
AutoAlign	---
Coil Select Mode	Default

System - Adjustments

B0 Shim mode	Cardiac
Adjust with body coil	On
Confirm freq. adjustment	Off
Assume Dominant Fat	Off
Assume Silicone	Off
Adjustment Tolerance	Auto

System - Adjust Volume

! Position	Isocenter
! Orientation	Transversal
! Rotation	0.00 deg
! A >> P	100 mm
! R >> L	100 mm
! F >> H	50 mm
Reset	Off

System - Tx/Rx

Frequency 1H	63.673772 MHz
Correction factor	1
Gain	High
Img. Scale Cor.	1.000
Reset	Off
? Ref. amplitude 1H	0.000 V

Physio - Signal1

1st Signal/Mode	ECG/Trigger
Average cycle	947 ± 122 ms
Average cycle	No Signal ms
Captured cycle	947 ± 122 ms
Acquisition window	830 ms
Trigger pulse	1
Trigger delay	0 ms
TR	824.00 ms
Concatenations	1
Segments	72
Phases	1
Adaptive Triggering	Off

Physio - Cardiac

Tagging	None
Magn. preparation	Non-sel. SR perf
T1	550 ms
T2 prep. duration 1	0 ms
T2 prep. duration 2	0 ms
T2 prep. duration 3	0 ms
T2 prep. duration 4	0 ms
T2 prep. duration 5	0 ms

Physio - Cardiac

T2 prep. duration 6	0 ms
T2 prep. duration 7	0 ms
T2 prep. duration 8	55 ms
T2 prep. duration 9	55 ms
Fat suppr.	None
Dark blood	Off
FoV read	360 mm
FoV phase	75.0 %
Phase resolution	75 %
Cine	Off
Trajectory	Cartesian
Dummy heartbeats	0
Motion Correction	Standard

Physio - PACE

Resp. control	Off
Concatenations	1

Inline - Common

Subtract	Off
Measurements	1
StdDev	Off
Motion Correction	Standard
Save original images	On

Inline - Cardiac

Inline Evaluation	T1 map
Magn. preparation	Non-sel. SR perf
Num. of preps	9
Sampling duration 1	1 beats
Sampling duration 2	1 beats
Sampling duration 3	1 beats
Sampling duration 4	1 beats
Sampling duration 5	1 beats
Sampling duration 6	1 beats
Sampling duration 7	1 beats
Sampling duration 8	1 beats
Sampling duration 9	1 beats
T2 prep. duration 1	0 ms
T2 prep. duration 2	0 ms
T2 prep. duration 3	0 ms
T2 prep. duration 4	0 ms
T2 prep. duration 5	0 ms
T2 prep. duration 6	0 ms
T2 prep. duration 7	0 ms
T2 prep. duration 8	55 ms
T2 prep. duration 9	55 ms
Contrasts	1
TE	1.26 ms
TR	824.00 ms
Recovery duration 1	1 beats
Recovery duration 2	1 beats
Recovery duration 3	0 beats
Recovery duration 4	0 beats
Recovery duration 5	0 beats
Recovery duration 6	0 beats
Recovery duration 7	0 beats
Recovery duration 8	1 beats
Recovery duration 9	1 beats
Motion Correction	Standard
Save original images	On

Inline - MIP

MIP-Sag	Off
---------	-----

Inline - MIP

MIP-Cor	Off
MIP-Tra	Off
MIP-Time	Off
Save original images	On

Inline - Composing

Inline Composing	Off
Distortion Corr.	On
Mode	2D
Unfiltered images	Off

Sequence - Part 1

Introduction	Off
Dimension	2D
Reordering	Linear
Asymmetric echo	Weak
Contrasts	1
Optimization	Min. TE
Multi-slice mode	Sequential
Echo spacing	2.9 ms
Sequence type	Trufi
Bandwidth	1085 Hz/Px

Sequence - Part 2

Define	Shots
Shots per slice	1
Segments	72
Trufi delta freq.	0 Hz
RF pulse type	Fast
Gradient mode	Fast
Excitation	Slice-sel.
Flip angle mode	Variable
Cine	Off

Sequence - Special

Partition coeff map	Off
ECV map	Off
Synth ECV map	Off
T1 array 0	550 ms
T1 array 1	550 ms
T1 array 2	550 ms
T1 array 3	550 ms
T1 array 4	550 ms
T1 array 5	550 ms
T1 array 6	550 ms
T1 array 7	550 ms
T1 array 8	550 ms
High-contrast acquisition	On
Error map	Off
SPAIR delay	5 ms
T1 sampling scheme	Custom

Sequence - Assistant

Mode	Off
Allowed delay	10 s

\\CMIC\Physics\T2MES\T2MES_low_res\pre FB T1T2 high HR

TA: 1:05 PM: FIX Voxel size: 1.4×1.4×8.0 mmPAT: 2 Rel. SNR: 1.00 : tfi

Properties

Prio recon	Off
Load images to viewer	On
Inline movie	Off
Auto store images	On
Load images to stamp segments	Off
Load images to graphic segments	On
Auto open inline display	Off
Auto close inline display	Off
Start measurement without further preparation	On
Wait for user to start	Off
Start measurements	Single measurement

Routine

Slice group	1
Slices	1
Dist. factor	20 %
Position	Isocenter
Orientation	Transversal
Phase enc. dir.	A >> P
AutoAlign	---
Phase oversampling	0 %
FoV read	360 mm
FoV phase	75.0 %
Slice thickness	8.0 mm
TR	824.00 ms
TE	1.26 ms
Averages	3
Concatenations	1
Filter	Distortion Corr.(2D)
Coil elements	BO1-3;SP2-4

Contrast - Common

TR	824.00 ms
TE	1.26 ms
Magn. preparation	Non-sel. SR perf
T1	350 ms
T2 prep. duration 1	0 ms
T2 prep. duration 2	0 ms
T2 prep. duration 3	0 ms
T2 prep. duration 4	0 ms
T2 prep. duration 5	0 ms
T2 prep. duration 6	0 ms
T2 prep. duration 7	0 ms
T2 prep. duration 8	55 ms
T2 prep. duration 9	55 ms
Flip angle	100 deg
Fat suppr.	None
Wrap-up Magn.	None

Contrast - Dynamic

Averages	3
Averaging mode	Short term
Reconstruction	Magnitude
Measurements	1
Multiple series	Off

Resolution - Common

FoV read	360 mm
FoV phase	75.0 %

Resolution - Common

Slice thickness	8.0 mm
Base resolution	256
Phase resolution	75 %
Phase partial Fourier	7/8
Trajectory	Cartesian
Interpolation	Off

Resolution - iPAT

PAT mode	GRAPPA
Accel. factor PE	2
Ref. lines PE	36
Reference scan mode	GRE/separate

Resolution - Filter Image

Image Filter	Off
Distortion Corr.	On
Mode	2D
Unfiltered images	Off
Prescan Normalize	Off
Normalize	Off
B1 filter	Off

Resolution - Filter Rawdata

Raw filter	Off
Elliptical filter	Off
POCS	Off

Geometry - Common

Slice group	1
Slices	1
Dist. factor	20 %
Position	Isocenter
Orientation	Transversal
Phase enc. dir.	A >> P
FoV read	360 mm
FoV phase	75.0 %
Slice thickness	8.0 mm
TR	824.00 ms
Multi-slice mode	Sequential
Series	Interleaved
Concatenations	1

Geometry - AutoAlign

Slice group	1
Position	Isocenter
Orientation	Transversal
Phase enc. dir.	A >> P
AutoAlign	---
Initial Position	Isocenter
L	0.0 mm
P	0.0 mm
H	0.0 mm
Initial Rotation	0.00 deg
Initial Orientation	Transversal

Geometry - Saturation

Fat suppr.	None
Wrap-up Magn.	None
Special sat.	None

Geometry - Navigator**System - Miscellaneous**

Positioning mode	FIX
Table position	H
Table position	0 mm
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Coil Combine Mode	Adaptive Combine
Save uncombined	Off
Matrix Optimization	Off
Coil Focus	Flat
AutoAlign	---
Coil Select Mode	Default

System - Adjustments

B0 Shim mode	Cardiac
Adjust with body coil	On
Confirm freq. adjustment	Off
Assume Dominant Fat	Off
Assume Silicone	Off
Adjustment Tolerance	Auto

System - Adjust Volume

! Position	Isocenter
! Orientation	Transversal
! Rotation	0.00 deg
! A >> P	150 mm
! R >> L	150 mm
! F >> H	50 mm
Reset	Off

System - Tx/Rx

Frequency 1H	63.673772 MHz
Correction factor	1
Gain	High
Img. Scale Cor.	1.000
Reset	Off
? Ref. amplitude 1H	0.000 V

Physio - Signal1

1st Signal/Mode	ECG/Trigger
Average cycle	947 ± 122 ms
Average cycle	No Signal ms
Captured cycle	947 ± 122 ms
Acquisition window	830 ms
Trigger pulse	1
Trigger delay	0 ms
TR	824.00 ms
Concatenations	1
Segments	72
Phases	1
Adaptive Triggering	Off

Physio - Cardiac

Tagging	None
Magn. preparation	Non-sel. SR perf
T1	350 ms
T2 prep. duration 1	0 ms
T2 prep. duration 2	0 ms
T2 prep. duration 3	0 ms
T2 prep. duration 4	0 ms
T2 prep. duration 5	0 ms

Physio - Cardiac

T2 prep. duration 6	0 ms
T2 prep. duration 7	0 ms
T2 prep. duration 8	55 ms
T2 prep. duration 9	55 ms
Fat suppr.	None
Dark blood	Off
FoV read	360 mm
FoV phase	75.0 %
Phase resolution	75 %
Cine	Off
Trajectory	Cartesian
Dummy heartbeats	0
Motion Correction	Standard

Physio - PACE

Resp. control	Off
Concatenations	1

Inline - Common

Subtract	Off
Measurements	1
StdDev	Off
Motion Correction	Standard
Save original images	On

Inline - Cardiac

Inline Evaluation	T1 map
Magn. preparation	Non-sel. SR perf
Num. of preps	9
Sampling duration 1	1 beats
Sampling duration 2	1 beats
Sampling duration 3	1 beats
Sampling duration 4	1 beats
Sampling duration 5	1 beats
Sampling duration 6	1 beats
Sampling duration 7	1 beats
Sampling duration 8	1 beats
Sampling duration 9	1 beats
T2 prep. duration 1	0 ms
T2 prep. duration 2	0 ms
T2 prep. duration 3	0 ms
T2 prep. duration 4	0 ms
T2 prep. duration 5	0 ms
T2 prep. duration 6	0 ms
T2 prep. duration 7	0 ms
T2 prep. duration 8	55 ms
T2 prep. duration 9	55 ms
Contrasts	1
TE	1.26 ms
TR	824.00 ms
Recovery duration 1	2 beats
Recovery duration 2	2 beats
Recovery duration 3	0 beats
Recovery duration 4	0 beats
Recovery duration 5	0 beats
Recovery duration 6	0 beats
Recovery duration 7	0 beats
Recovery duration 8	3 beats
Recovery duration 9	3 beats
Motion Correction	Standard
Save original images	On

Inline - MIP

MIP-Sag	Off
---------	-----

Inline - MIP

MIP-Cor	Off
MIP-Tra	Off
MIP-Time	Off
Save original images	On

Inline - Composing

Inline Composing	Off
Distortion Corr.	On
Mode	2D
Unfiltered images	Off

Sequence - Part 1

Introduction	Off
Dimension	2D
Reordering	Linear
Asymmetric echo	Weak
Contrasts	1
Optimization	Min. TE
Multi-slice mode	Sequential
Echo spacing	2.9 ms
Sequence type	Trufi
Bandwidth	1085 Hz/Px

Sequence - Part 2

Define	Shots
Shots per slice	1
Segments	72
Trufi delta freq.	0 Hz
RF pulse type	Fast
Gradient mode	Fast
Excitation	Slice-sel.
Flip angle mode	Variable
Cine	Off

Sequence - Special

Partition coeff map	Off
ECV map	Off
Synth ECV map	Off
T1 array 0	350 ms
T1 array 1	350 ms
T1 array 2	350 ms
T1 array 3	350 ms
T1 array 4	350 ms
T1 array 5	350 ms
T1 array 6	350 ms
T1 array 7	350 ms
T1 array 8	350 ms
High-contrast acquisition	On
Error map	Off
SPAIR delay	5 ms
T1 sampling scheme	Fixed TS

Sequence - Assistant

Mode	Off
Allowed delay	10 s

\\CMIC\Physics\T2MES\T2MES_low_res\post FB T1T2

TA: 0:43 PM: FIX Voxel size: 1.4×1.4×8.0 mmPAT: 2 Rel. SNR: 1.00 : tfi

Properties

Prio recon	Off
Load images to viewer	On
Inline movie	Off
Auto store images	On
Load images to stamp segments	Off
Load images to graphic segments	On
Auto open inline display	Off
Auto close inline display	Off
Start measurement without further preparation	On
Wait for user to start	Off
Start measurements	Single measurement

Routine

Slice group	1
Slices	1
Dist. factor	20 %
Position	Isocenter
Orientation	Transversal
Phase enc. dir.	A >> P
AutoAlign	---
Phase oversampling	0 %
FoV read	360 mm
FoV phase	75.0 %
Slice thickness	8.0 mm
TR	904.00 ms
TE	1.26 ms
Averages	3
Concatenations	1
Filter	Distortion Corr.(2D)
Coil elements	BO1-3;SP2-4

Contrast - Common

TR	904.00 ms
TE	1.26 ms
Magn. preparation	Non-sel. SR perf
T1	300 ms
T2 prep. duration 1	0 ms
T2 prep. duration 2	0 ms
T2 prep. duration 3	0 ms
T2 prep. duration 4	0 ms
T2 prep. duration 5	0 ms
T2 prep. duration 6	0 ms
T2 prep. duration 7	0 ms
T2 prep. duration 8	55 ms
T2 prep. duration 9	55 ms
Flip angle	100 deg
Fat suppr.	None
Wrap-up Magn.	None

Contrast - Dynamic

Averages	3
Averaging mode	Short term
Reconstruction	Magnitude
Measurements	1
Multiple series	Off

Resolution - Common

FoV read	360 mm
FoV phase	75.0 %

Resolution - Common

Slice thickness	8.0 mm
Base resolution	256
Phase resolution	75 %
Phase partial Fourier	7/8
Trajectory	Cartesian
Interpolation	Off

Resolution - iPAT

PAT mode	GRAPPA
Accel. factor PE	2
Ref. lines PE	36
Reference scan mode	GRE/separate

Resolution - Filter Image

Image Filter	Off
Distortion Corr.	On
Mode	2D
Unfiltered images	Off
Prescan Normalize	Off
Normalize	Off
B1 filter	Off

Resolution - Filter Rawdata

Raw filter	Off
Elliptical filter	Off
POCS	Off

Geometry - Common

Slice group	1
Slices	1
Dist. factor	20 %
Position	Isocenter
Orientation	Transversal
Phase enc. dir.	A >> P
FoV read	360 mm
FoV phase	75.0 %
Slice thickness	8.0 mm
TR	904.00 ms
Multi-slice mode	Sequential
Series	Interleaved
Concatenations	1

Geometry - AutoAlign

Slice group	1
Position	Isocenter
Orientation	Transversal
Phase enc. dir.	A >> P
AutoAlign	---
Initial Position	Isocenter
L	0.0 mm
P	0.0 mm
H	0.0 mm
Initial Rotation	0.00 deg
Initial Orientation	Transversal

Geometry - Saturation

Fat suppr.	None
Wrap-up Magn.	None
Special sat.	None

Geometry - Navigator**System - Miscellaneous**

Positioning mode	FIX
Table position	H
Table position	0 mm
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Coil Combine Mode	Adaptive Combine
Save uncombined	Off
Matrix Optimization	Off
Coil Focus	Flat
AutoAlign	---
Coil Select Mode	Default

System - Adjustments

B0 Shim mode	Cardiac
Adjust with body coil	On
Confirm freq. adjustment	Off
Assume Dominant Fat	Off
Assume Silicone	Off
Adjustment Tolerance	Auto

System - Adjust Volume

! Position	Isocenter
! Orientation	Transversal
! Rotation	0.00 deg
! A >> P	150 mm
! R >> L	150 mm
! F >> H	50 mm
Reset	Off

System - Tx/Rx

Frequency 1H	63.673772 MHz
Correction factor	1
Gain	High
Img. Scale Cor.	1.000
Reset	Off
? Ref. amplitude 1H	0.000 V

Physio - Signal1

1st Signal/Mode	ECG/Trigger
Average cycle	947 ± 122 ms
Average cycle	No Signal ms
Captured cycle	947 ± 122 ms
Acquisition window	930 ms
Trigger pulse	1
Trigger delay	0 ms
TR	904.00 ms
Concatenations	1
Segments	72
Phases	1
Adaptive Triggering	Off

Physio - Cardiac

Tagging	None
Magn. preparation	Non-sel. SR perf
T1	300 ms
T2 prep. duration 1	0 ms
T2 prep. duration 2	0 ms
T2 prep. duration 3	0 ms
T2 prep. duration 4	0 ms
T2 prep. duration 5	0 ms

Physio - Cardiac

T2 prep. duration 6	0 ms
T2 prep. duration 7	0 ms
T2 prep. duration 8	55 ms
T2 prep. duration 9	55 ms
Fat suppr.	None
Dark blood	Off
FoV read	360 mm
FoV phase	75.0 %
Phase resolution	75 %
Cine	Off
Trajectory	Cartesian
Dummy heartbeats	0
Motion Correction	Standard

Physio - PACE

Resp. control	Off
Concatenations	1

Inline - Common

Subtract	Off
Measurements	1
StdDev	Off
Motion Correction	Standard
Save original images	On

Inline - Cardiac

Inline Evaluation	T1 map
Magn. preparation	Non-sel. SR perf
Num. of preps	9
Sampling duration 1	1 beats
Sampling duration 2	1 beats
Sampling duration 3	1 beats
Sampling duration 4	1 beats
Sampling duration 5	1 beats
Sampling duration 6	1 beats
Sampling duration 7	1 beats
Sampling duration 8	1 beats
Sampling duration 9	1 beats
T2 prep. duration 1	0 ms
T2 prep. duration 2	0 ms
T2 prep. duration 3	0 ms
T2 prep. duration 4	0 ms
T2 prep. duration 5	0 ms
T2 prep. duration 6	0 ms
T2 prep. duration 7	0 ms
T2 prep. duration 8	55 ms
T2 prep. duration 9	55 ms
Contrasts	1
TE	1.26 ms
TR	904.00 ms
Recovery duration 1	0 beats
Recovery duration 2	0 beats
Recovery duration 3	0 beats
Recovery duration 4	0 beats
Recovery duration 5	0 beats
Recovery duration 6	0 beats
Recovery duration 7	0 beats
Recovery duration 8	1 beats
Recovery duration 9	1 beats
Motion Correction	Standard
Save original images	On

Inline - MIP

MIP-Sag	Off
---------	-----

Inline - MIP

MIP-Cor	Off
MIP-Tra	Off
MIP-Time	Off
Save original images	On

Inline - Composing

Inline Composing	Off
Distortion Corr.	On
Mode	2D
Unfiltered images	Off

Sequence - Part 1

Introduction	Off
Dimension	2D
Reordering	Linear
Asymmetric echo	Weak
Contrasts	1
Optimization	Min. TE
Multi-slice mode	Sequential
Echo spacing	2.9 ms
Sequence type	Trufi
Bandwidth	1085 Hz/Px

Sequence - Part 2

Define	Shots
Shots per slice	1
Segments	72
Trufi delta freq.	0 Hz
RF pulse type	Fast
Gradient mode	Fast
Excitation	Slice-sel.
Flip angle mode	Variable
Cine	Off

Sequence - Special

Partition coeff map	Off
ECV map	Off
Synth ECV map	Off
T1 array 0	300 ms
T1 array 1	300 ms
T1 array 2	300 ms
T1 array 3	300 ms
T1 array 4	300 ms
T1 array 5	300 ms
T1 array 6	300 ms
T1 array 7	300 ms
T1 array 8	300 ms
High-contrast acquisition	On
Error map	Off
SPAIR delay	5 ms
T1 sampling scheme	Fixed TS

Sequence - Assistant

Mode	Off
Allowed delay	10 s

\\CMIC\Physics\T2MES\T2MES_low_res\pre_MOLLI_5s(3s)3s_256

TA: 0:11 PM: FIX Voxel size: 1.4×1.4×8.0 mmPAT: 2 Rel. SNR: 1.00 : tfi

Properties

Prio recon	Off
Load images to viewer	On
Inline movie	Off
Auto store images	On
Load images to stamp segments	Off
Load images to graphic segments	On
Auto open inline display	Off
Auto close inline display	Off
Start measurement without further preparation	On
Wait for user to start	Off
Start measurements	Single measurement

Routine

Slice group	1
Slices	1
Dist. factor	20 %
Position	Isocenter
Orientation	Transversal
Phase enc. dir.	A >> P
AutoAlign	---
Phase oversampling	0 %
FoV read	360 mm
FoV phase	75.0 %
Slice thickness	8.0 mm
TR	280.56 ms
TE	1.12 ms
Averages	1
Concatenations	1
Filter	Distortion Corr.(2D)
Coil elements	BO1-3;SP2-4

Contrast - Common

TR	280.56 ms
TE	1.12 ms
Magn. preparation	Non-sel. IR T1map
T1	180 ms
Flip angle	35 deg
Fat suppr.	None
Wrap-up Magn.	None

Contrast - Dynamic

Averages	1
Averaging mode	Short term
Reconstruction	Magnitude
Measurements	1
Multiple series	Off

Resolution - Common

FoV read	360 mm
FoV phase	75.0 %
Slice thickness	8.0 mm
Base resolution	256
Phase resolution	75 %
Phase partial Fourier	7/8
Trajectory	Cartesian
Interpolation	Off

Resolution - iPAT

PAT mode	GRAPPA
----------	--------

Resolution - iPAT

Accel. factor PE	2
Ref. lines PE	36
Reference scan mode	GRE/separate

Resolution - Filter Image

Image Filter	Off
Distortion Corr.	On
Mode	2D
Unfiltered images	Off
Prescan Normalize	Off
Normalize	Off
B1 filter	Off

Resolution - Filter Rawdata

Raw filter	Off
Elliptical filter	Off
POCS	Off

Geometry - Common

Slice group	1
Slices	1
Dist. factor	20 %
Position	Isocenter
Orientation	Transversal
Phase enc. dir.	A >> P
FoV read	360 mm
FoV phase	75.0 %
Slice thickness	8.0 mm
TR	280.56 ms
Multi-slice mode	Sequential
Series	Interleaved
Concatenations	1

Geometry - AutoAlign

Slice group	1
Position	Isocenter
Orientation	Transversal
Phase enc. dir.	A >> P
AutoAlign	---
Initial Position	Isocenter
L	0.0 mm
P	0.0 mm
H	0.0 mm
Initial Rotation	0.00 deg
Initial Orientation	Transversal

Geometry - Saturation

Fat suppr.	None
Wrap-up Magn.	None
Special sat.	None

Geometry - Navigator**System - Miscellaneous**

Positioning mode	FIX
Table position	H
Table position	0 mm
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P

System - Miscellaneous

Transversal	F >> H
Coil Combine Mode	Adaptive Combine
Save uncombined	Off
Matrix Optimization	Off
Coil Focus	Flat
AutoAlign	---
Coil Select Mode	Default

System - Adjustments

B0 Shim mode	Standard
Adjust with body coil	On
Confirm freq. adjustment	Off
Assume Dominant Fat	Off
Assume Silicone	Off
Adjustment Tolerance	Auto

System - Adjust Volume

! Position	Isocenter
! Orientation	Transversal
! Rotation	0.00 deg
! A >> P	100 mm
! R >> L	100 mm
! F >> H	25 mm
Reset	Off

System - Tx/Rx

Frequency 1H	63.673772 MHz
Correction factor	1
Gain	High
Img. Scale Cor.	1.000
Reset	Off
? Ref. amplitude 1H	0.000 V

Physio - Signal1

1st Signal/Mode	ECG/Trigger
Average cycle	947 ± 122 ms
Average cycle	No Signal ms
Captured cycle	947 ± 122 ms
Acquisition window	700 ms
Trigger pulse	1
Trigger delay	418 ms
TR	280.56 ms
Concatenations	1
Segments	72
Phases	1
Adaptive Triggering	Off

Physio - Cardiac

Tagging	None
Magn. preparation	Non-sel. IR T1map
TI	180 ms
Fat suppr.	None
Dark blood	Off
FoV read	360 mm
FoV phase	75.0 %
Phase resolution	75 %
Cine	Off
Trajectory	Cartesian
Dummy heartbeats	0
Motion Correction	Standard

Physio - PACE

Resp. control	Off
Concatenations	1

Sequence - Part 1

Introduction	Off
Dimension	2D
Reordering	Linear
Asymmetric echo	Weak
Contrasts	1
Optimization	Min. TE TR
Multi-slice mode	Sequential
Sequence type	Trufi
Bandwidth	1085 Hz/Px

Sequence - Part 2

Define	Shots
Shots per slice	1
Segments	72
Trufi delta freq.	0 Hz
RF pulse type	Fast
Gradient mode	Fast
Excitation	Slice-sel.
Flip angle mode	Constant
Cine	Off

Sequence - Special

Partition coeff map	Off
ECV map	Off
Synth ECV map	Off
T1 scout mode	Off
Error map	On
Synth PSIR	Off
Periods in seconds	On
16 bit images	Off
T1 sampling scheme	Native

Sequence - Assistant

Mode	Off
Allowed delay	0 s

\\CMIC\Physics\T2MES\T2MES_low_res\post_MOLLI_4s(1s)3s(1s)2s_256

TA: 0:12 PM: FIX Voxel size: 1.4×1.4×8.0 mmPAT: 2 Rel. SNR: 1.00 : tfi

Properties

Prio recon	Off
Load images to viewer	On
Inline movie	Off
Auto store images	On
Load images to stamp segments	Off
Load images to graphic segments	On
Auto open inline display	Off
Auto close inline display	Off
Start measurement without further preparation	On
Wait for user to start	Off
Start measurements	Single measurement

Resolution - iPAT

Accel. factor PE	2
Ref. lines PE	36
Reference scan mode	GRE/separate

Resolution - Filter Image

Image Filter	Off
Distortion Corr.	On
Mode	2D
Unfiltered images	Off
Prescan Normalize	Off
Normalize	Off
B1 filter	Off

Routine

Slice group	1
Slices	1
Dist. factor	20 %
Position	Isocenter
Orientation	Transversal
Phase enc. dir.	A >> P
AutoAlign	---
Phase oversampling	0 %
FoV read	360 mm
FoV phase	75.0 %
Slice thickness	8.0 mm
TR	360.56 ms
TE	1.12 ms
Averages	1
Concatenations	1
Filter	Distortion Corr.(2D)
Coil elements	BO1-3;SP2-4

Resolution - Filter Rawdata

Raw filter	Off
Elliptical filter	Off
POCS	Off

Geometry - Common

Slice group	1
Slices	1
Dist. factor	20 %
Position	Isocenter
Orientation	Transversal
Phase enc. dir.	A >> P
FoV read	360 mm
FoV phase	75.0 %
Slice thickness	8.0 mm
TR	360.56 ms
Multi-slice mode	Sequential
Series	Interleaved
Concatenations	1

Contrast - Common

TR	360.56 ms
TE	1.12 ms
Magn. preparation	Non-sel. IR T1map
T1	260 ms
Flip angle	35 deg
Fat suppr.	None
Wrap-up Magn.	None

Geometry - AutoAlign

Slice group	1
Position	Isocenter
Orientation	Transversal
Phase enc. dir.	A >> P
AutoAlign	---
Initial Position	Isocenter
L	0.0 mm
P	0.0 mm
H	0.0 mm
Initial Rotation	0.00 deg
Initial Orientation	Transversal

Contrast - Dynamic

Averages	1
Averaging mode	Short term
Reconstruction	Magnitude
Measurements	1
Multiple series	Off

Resolution - Common

FoV read	360 mm
FoV phase	75.0 %
Slice thickness	8.0 mm
Base resolution	256
Phase resolution	75 %
Phase partial Fourier	7/8
Trajectory	Cartesian
Interpolation	Off

Geometry - Saturation

Fat suppr.	None
Wrap-up Magn.	None
Special sat.	None

Geometry - Navigator**System - Miscellaneous**

Positioning mode	FIX
Table position	H
Table position	0 mm
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P

Resolution - iPAT

PAT mode	GRAPPA
----------	--------

System - Miscellaneous

Transversal	F >> H
Coil Combine Mode	Adaptive Combine
Save uncombined	Off
Matrix Optimization	Off
Coil Focus	Flat
AutoAlign	---
Coil Select Mode	Default

System - Adjustments

B0 Shim mode	Standard
Adjust with body coil	On
Confirm freq. adjustment	Off
Assume Dominant Fat	Off
Assume Silicone	Off
Adjustment Tolerance	Auto

System - Adjust Volume

! Position	Isocenter
! Orientation	Transversal
! Rotation	0.00 deg
! A >> P	100 mm
! R >> L	100 mm
! F >> H	25 mm
Reset	Off

System - Tx/Rx

Frequency 1H	63.673772 MHz
Correction factor	1
Gain	High
Img. Scale Cor.	1.000
Reset	Off
? Ref. amplitude 1H	0.000 V

Physio - Signal1

1st Signal/Mode	ECG/Trigger
Average cycle	947 ± 122 ms
Average cycle	No Signal ms
Captured cycle	947 ± 122 ms
Acquisition window	831 ms
Trigger pulse	1
Trigger delay	470 ms
TR	360.56 ms
Concatenations	1
Segments	72
Phases	1
Adaptive Triggering	Off

Physio - Cardiac

Tagging	None
Magn. preparation	Non-sel. IR T1map
TI	260 ms
Fat suppr.	None
Dark blood	Off
FoV read	360 mm
FoV phase	75.0 %
Phase resolution	75 %
Cine	Off
Trajectory	Cartesian
Dummy heartbeats	0
Motion Correction	Standard

Physio - PACE

Resp. control	Off
Concatenations	1

Sequence - Part 1

Introduction	Off
Dimension	2D
Reordering	Linear
Asymmetric echo	Weak
Contrasts	1
Optimization	Min. TE TR
Multi-slice mode	Sequential
Sequence type	Trufi
Bandwidth	1085 Hz/Px

Sequence - Part 2

Define	Shots
Shots per slice	1
Segments	72
Trufi delta freq.	0 Hz
RF pulse type	Fast
Gradient mode	Fast
Excitation	Slice-sel.
Flip angle mode	Constant
Cine	Off

Sequence - Special

Partition coeff map	On
ECV map	Off
Synth ECV map	On
T1 scout mode	Off
Error map	On
Synth PSIR	On
Periods in seconds	On
16 bit images	Off
T1 sampling scheme	Post Gd

Sequence - Assistant

Mode	Off
Allowed delay	0 s

\\CMIC\Physics\T2MES\T2MES_low_res\T2map 0-25-55 new SSFP

TA: 0:12 PM: FIX Voxel size: 1.9×1.9×8.0 mmPAT: 2 Rel. SNR: 1.00 : tfi

Properties

Prio recon	Off
Load images to viewer	On
Inline movie	Off
Auto store images	On
Load images to stamp segments	Off
Load images to graphic segments	On
Auto open inline display	Off
Auto close inline display	Off
Start measurement without further preparation	On
Wait for user to start	Off
Start measurements	Single measurement

Routine

Slice group	1
Slices	1
Dist. factor	20 %
Position	Isocenter
Orientation	Transversal
Phase enc. dir.	A >> P
AutoAlign	---
Phase oversampling	0 %
FoV read	360 mm
FoV phase	75.0 %
Slice thickness	8.0 mm
TR	207.98 ms
TE	1.06 ms
Averages	1
Concatenations	1
Filter	Distortion Corr.(2D)
Coil elements	BO1-3;SP2-4

Contrast - Common

TR	207.98 ms
TE	1.06 ms
Magn. preparation	T2 prep. adiab.
T2 prep. duration 1	0 ms
T2 prep. duration 2	25 ms
T2 prep. duration 3	55 ms
Flip angle	70 deg
Fat suppr.	None
Wrap-up Magn.	None

Contrast - Dynamic

Averages	1
Averaging mode	Short term
Reconstruction	Magnitude
Measurements	1
Multiple series	Off

Resolution - Common

FoV read	360 mm
FoV phase	75.0 %
Slice thickness	8.0 mm
Base resolution	192
Phase resolution	83 %
Phase partial Fourier	7/8
Trajectory	Cartesian
Interpolation	Off

Resolution - iPAT

PAT mode	GRAPPA
Accel. factor PE	2
Ref. lines PE	32
Reference scan mode	GRE/separate

Resolution - Filter Image

Image Filter	Off
Distortion Corr.	On
Mode	2D
Unfiltered images	Off
Prescan Normalize	Off
Normalize	Off
B1 filter	Off

Resolution - Filter Rawdata

Raw filter	Off
Elliptical filter	Off
POCS	Off

Geometry - Common

Slice group	1
Slices	1
Dist. factor	20 %
Position	Isocenter
Orientation	Transversal
Phase enc. dir.	A >> P
FoV read	360 mm
FoV phase	75.0 %
Slice thickness	8.0 mm
TR	207.98 ms
Multi-slice mode	Sequential
Series	Interleaved
Concatenations	1

Geometry - AutoAlign

Slice group	1
Position	Isocenter
Orientation	Transversal
Phase enc. dir.	A >> P
AutoAlign	---
Initial Position	Isocenter
L	0.0 mm
P	0.0 mm
H	0.0 mm
Initial Rotation	0.00 deg
Initial Orientation	Transversal

Geometry - Saturation

Fat suppr.	None
Wrap-up Magn.	None
Special sat.	None

Geometry - Navigator**System - Miscellaneous**

Positioning mode	FIX
Table position	H
Table position	0 mm
MSMA	S - C - T
Sagittal	R >> L

System - Miscellaneous

Coronal	A >> P
Transversal	F >> H
Coil Combine Mode	Adaptive Combine
Save uncombined	Off
Matrix Optimization	Off
Coil Focus	Flat
AutoAlign	---
Coil Select Mode	Default

System - Adjustments

B0 Shim mode	Standard
Adjust with body coil	On
Confirm freq. adjustment	Off
Assume Dominant Fat	Off
Assume Silicone	Off
Adjustment Tolerance	Auto

System - Adjust Volume

! Position	Isocenter
! Orientation	Transversal
! Rotation	0.00 deg
! A >> P	100 mm
! R >> L	100 mm
! F >> H	25 mm
Reset	Off

System - Tx/Rx

Frequency 1H	63.673772 MHz
Correction factor	1
Gain	High
Img. Scale Cor.	1.000
Reset	Off
? Ref. amplitude 1H	0.000 V

Physio - Signal1

1st Signal/Mode	ECG/Trigger
Average cycle	947 ± 122 ms
Average cycle	No Signal ms
Captured cycle	947 ± 122 ms
Acquisition window	831 ms
Trigger pulse	1
Trigger delay	623 ms
TR	207.98 ms
Concatenations	1
Segments	60
Phases	1
Adaptive Triggering	Off

Physio - Cardiac

Tagging	None
Magn. preparation	T2 prep. adiab.
T2 prep. duration 1	0 ms
T2 prep. duration 2	25 ms
T2 prep. duration 3	55 ms
Fat suppr.	None
Dark blood	Off
FoV read	360 mm
FoV phase	75.0 %
Phase resolution	83 %
Cine	Off
Trajectory	Cartesian
Dummy heartbeats	0
Motion Correction	None

Physio - PACE

Resp. control	Off
Concatenations	1

Sequence - Part 1

Introduction	Off
Dimension	2D
Reordering	Linear
Asymmetric echo	Weak
Contrasts	1
Optimization	Min. TE TR
Multi-slice mode	Sequential
Sequence type	Trufi
Bandwidth	1184 Hz/Px

Sequence - Part 2

Define	Shots
Shots per slice	1
Segments	60
Trufi delta freq.	0 Hz
RF pulse type	Fast
Gradient mode	Fast
Excitation	Slice-sel.
Flip angle mode	Constant
Cine	Off

Sequence - Special

Periods in seconds	On
16 bit images	Off

Sequence - Assistant

Mode	Off
Allowed delay	0 s

\\CMIC\Physics\T2MES\T2MES_low_res\pre FB T1T2

TA: 0:48 PM: FIX Voxel size: 1.4×1.4×8.0 mmPAT: 2 Rel. SNR: 1.00 : tfi

Properties

Prio recon	Off
Load images to viewer	On
Inline movie	Off
Auto store images	On
Load images to stamp segments	Off
Load images to graphic segments	On
Auto open inline display	Off
Auto close inline display	Off
Start measurement without further preparation	On
Wait for user to start	Off
Start measurements	Single measurement

Routine

Slice group	1
Slices	1
Dist. factor	20 %
Position	Isocenter
Orientation	Transversal
Phase enc. dir.	A >> P
AutoAlign	---
Phase oversampling	0 %
FoV read	360 mm
FoV phase	75.0 %
Slice thickness	8.0 mm
TR	824.00 ms
TE	1.26 ms
Averages	3
Concatenations	1
Filter	Distortion Corr.(2D)
Coil elements	BO1-3;SP2-4

Contrast - Common

TR	824.00 ms
TE	1.26 ms
Magn. preparation	Non-sel. SR perf
T1	550 ms
T2 prep. duration 1	0 ms
T2 prep. duration 2	0 ms
T2 prep. duration 3	0 ms
T2 prep. duration 4	0 ms
T2 prep. duration 5	0 ms
T2 prep. duration 6	0 ms
T2 prep. duration 7	0 ms
T2 prep. duration 8	55 ms
T2 prep. duration 9	55 ms
Flip angle	100 deg
Fat suppr.	None
Wrap-up Magn.	None

Contrast - Dynamic

Averages	3
Averaging mode	Short term
Reconstruction	Magnitude
Measurements	1
Multiple series	Off

Resolution - Common

FoV read	360 mm
FoV phase	75.0 %

Resolution - Common

Slice thickness	8.0 mm
Base resolution	256
Phase resolution	75 %
Phase partial Fourier	7/8
Trajectory	Cartesian
Interpolation	Off

Resolution - iPAT

PAT mode	GRAPPA
Accel. factor PE	2
Ref. lines PE	36
Reference scan mode	GRE/separate

Resolution - Filter Image

Image Filter	Off
Distortion Corr.	On
Mode	2D
Unfiltered images	Off
Prescan Normalize	Off
Normalize	Off
B1 filter	Off

Resolution - Filter Rawdata

Raw filter	Off
Elliptical filter	Off
POCS	Off

Geometry - Common

Slice group	1
Slices	1
Dist. factor	20 %
Position	Isocenter
Orientation	Transversal
Phase enc. dir.	A >> P
FoV read	360 mm
FoV phase	75.0 %
Slice thickness	8.0 mm
TR	824.00 ms
Multi-slice mode	Sequential
Series	Interleaved
Concatenations	1

Geometry - AutoAlign

Slice group	1
Position	Isocenter
Orientation	Transversal
Phase enc. dir.	A >> P
AutoAlign	---
Initial Position	Isocenter
L	0.0 mm
P	0.0 mm
H	0.0 mm
Initial Rotation	0.00 deg
Initial Orientation	Transversal

Geometry - Saturation

Fat suppr.	None
Wrap-up Magn.	None
Special sat.	None

Geometry - Navigator**System - Miscellaneous**

Positioning mode	FIX
Table position	H
Table position	0 mm
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Coil Combine Mode	Adaptive Combine
Save uncombined	Off
Matrix Optimization	Off
Coil Focus	Flat
AutoAlign	---
Coil Select Mode	Default

System - Adjustments

B0 Shim mode	Cardiac
Adjust with body coil	On
Confirm freq. adjustment	Off
Assume Dominant Fat	Off
Assume Silicone	Off
Adjustment Tolerance	Auto

System - Adjust Volume

! Position	Isocenter
! Orientation	Transversal
! Rotation	0.00 deg
! A >> P	100 mm
! R >> L	100 mm
! F >> H	50 mm
Reset	Off

System - Tx/Rx

Frequency 1H	63.673772 MHz
Correction factor	1
Gain	High
Img. Scale Cor.	1.000
Reset	Off
? Ref. amplitude 1H	0.000 V

Physio - Signal1

1st Signal/Mode	ECG/Trigger
Average cycle	947 ± 122 ms
Average cycle	No Signal ms
Captured cycle	947 ± 122 ms
Acquisition window	830 ms
Trigger pulse	1
Trigger delay	0 ms
TR	824.00 ms
Concatenations	1
Segments	72
Phases	1
Adaptive Triggering	Off

Physio - Cardiac

Tagging	None
Magn. preparation	Non-sel. SR perf
T1	550 ms
T2 prep. duration 1	0 ms
T2 prep. duration 2	0 ms
T2 prep. duration 3	0 ms
T2 prep. duration 4	0 ms
T2 prep. duration 5	0 ms

Physio - Cardiac

T2 prep. duration 6	0 ms
T2 prep. duration 7	0 ms
T2 prep. duration 8	55 ms
T2 prep. duration 9	55 ms
Fat suppr.	None
Dark blood	Off
FoV read	360 mm
FoV phase	75.0 %
Phase resolution	75 %
Cine	Off
Trajectory	Cartesian
Dummy heartbeats	0
Motion Correction	Standard

Physio - PACE

Resp. control	Off
Concatenations	1

Inline - Common

Subtract	Off
Measurements	1
StdDev	Off
Motion Correction	Standard
Save original images	On

Inline - Cardiac

Inline Evaluation	T1 map
Magn. preparation	Non-sel. SR perf
Num. of preps	9
Sampling duration 1	1 beats
Sampling duration 2	1 beats
Sampling duration 3	1 beats
Sampling duration 4	1 beats
Sampling duration 5	1 beats
Sampling duration 6	1 beats
Sampling duration 7	1 beats
Sampling duration 8	1 beats
Sampling duration 9	1 beats
T2 prep. duration 1	0 ms
T2 prep. duration 2	0 ms
T2 prep. duration 3	0 ms
T2 prep. duration 4	0 ms
T2 prep. duration 5	0 ms
T2 prep. duration 6	0 ms
T2 prep. duration 7	0 ms
T2 prep. duration 8	55 ms
T2 prep. duration 9	55 ms
Contrasts	1
TE	1.26 ms
TR	824.00 ms
Recovery duration 1	1 beats
Recovery duration 2	1 beats
Recovery duration 3	0 beats
Recovery duration 4	0 beats
Recovery duration 5	0 beats
Recovery duration 6	0 beats
Recovery duration 7	0 beats
Recovery duration 8	1 beats
Recovery duration 9	1 beats
Motion Correction	Standard
Save original images	On

Inline - MIP

MIP-Sag	Off
---------	-----

Inline - MIP

MIP-Cor	Off
MIP-Tra	Off
MIP-Time	Off
Save original images	On

Inline - Composing

Inline Composing	Off
Distortion Corr.	On
Mode	2D
Unfiltered images	Off

Sequence - Part 1

Introduction	Off
Dimension	2D
Reordering	Linear
Asymmetric echo	Weak
Contrasts	1
Optimization	Min. TE
Multi-slice mode	Sequential
Echo spacing	2.9 ms
Sequence type	Trufi
Bandwidth	1085 Hz/Px

Sequence - Part 2

Define	Shots
Shots per slice	1
Segments	72
Trufi delta freq.	0 Hz
RF pulse type	Fast
Gradient mode	Fast
Excitation	Slice-sel.
Flip angle mode	Variable
Cine	Off

Sequence - Special

Partition coeff map	Off
ECV map	Off
Synth ECV map	Off
T1 array 0	550 ms
T1 array 1	550 ms
T1 array 2	550 ms
T1 array 3	550 ms
T1 array 4	550 ms
T1 array 5	550 ms
T1 array 6	550 ms
T1 array 7	550 ms
T1 array 8	550 ms
High-contrast acquisition	On
Error map	Off
SPAIR delay	5 ms
T1 sampling scheme	Custom

Sequence - Assistant

Mode	Off
Allowed delay	10 s

\\CMIC\Physics\T2MES\T2MES_low_res\pre FB T1T2 high HR

TA: 1:05 PM: FIX Voxel size: 1.4×1.4×8.0 mmPAT: 2 Rel. SNR: 1.00 : tfi

Properties

Prio recon	Off
Load images to viewer	On
Inline movie	Off
Auto store images	On
Load images to stamp segments	Off
Load images to graphic segments	On
Auto open inline display	Off
Auto close inline display	Off
Start measurement without further preparation	On
Wait for user to start	Off
Start measurements	Single measurement

Routine

Slice group	1
Slices	1
Dist. factor	20 %
Position	Isocenter
Orientation	Transversal
Phase enc. dir.	A >> P
AutoAlign	---
Phase oversampling	0 %
FoV read	360 mm
FoV phase	75.0 %
Slice thickness	8.0 mm
TR	824.00 ms
TE	1.26 ms
Averages	3
Concatenations	1
Filter	Distortion Corr.(2D)
Coil elements	BO1-3;SP2-4

Contrast - Common

TR	824.00 ms
TE	1.26 ms
Magn. preparation	Non-sel. SR perf
T1	350 ms
T2 prep. duration 1	0 ms
T2 prep. duration 2	0 ms
T2 prep. duration 3	0 ms
T2 prep. duration 4	0 ms
T2 prep. duration 5	0 ms
T2 prep. duration 6	0 ms
T2 prep. duration 7	0 ms
T2 prep. duration 8	55 ms
T2 prep. duration 9	55 ms
Flip angle	100 deg
Fat suppr.	None
Wrap-up Magn.	None

Contrast - Dynamic

Averages	3
Averaging mode	Short term
Reconstruction	Magnitude
Measurements	1
Multiple series	Off

Resolution - Common

FoV read	360 mm
FoV phase	75.0 %

Resolution - Common

Slice thickness	8.0 mm
Base resolution	256
Phase resolution	75 %
Phase partial Fourier	7/8
Trajectory	Cartesian
Interpolation	Off

Resolution - iPAT

PAT mode	GRAPPA
Accel. factor PE	2
Ref. lines PE	36
Reference scan mode	GRE/separate

Resolution - Filter Image

Image Filter	Off
Distortion Corr.	On
Mode	2D
Unfiltered images	Off
Prescan Normalize	Off
Normalize	Off
B1 filter	Off

Resolution - Filter Rawdata

Raw filter	Off
Elliptical filter	Off
POCS	Off

Geometry - Common

Slice group	1
Slices	1
Dist. factor	20 %
Position	Isocenter
Orientation	Transversal
Phase enc. dir.	A >> P
FoV read	360 mm
FoV phase	75.0 %
Slice thickness	8.0 mm
TR	824.00 ms
Multi-slice mode	Sequential
Series	Interleaved
Concatenations	1

Geometry - AutoAlign

Slice group	1
Position	Isocenter
Orientation	Transversal
Phase enc. dir.	A >> P
AutoAlign	---
Initial Position	Isocenter
L	0.0 mm
P	0.0 mm
H	0.0 mm
Initial Rotation	0.00 deg
Initial Orientation	Transversal

Geometry - Saturation

Fat suppr.	None
Wrap-up Magn.	None
Special sat.	None

Geometry - Navigator**System - Miscellaneous**

Positioning mode	FIX
Table position	H
Table position	0 mm
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Coil Combine Mode	Adaptive Combine
Save uncombined	Off
Matrix Optimization	Off
Coil Focus	Flat
AutoAlign	---
Coil Select Mode	Default

System - Adjustments

B0 Shim mode	Cardiac
Adjust with body coil	On
Confirm freq. adjustment	Off
Assume Dominant Fat	Off
Assume Silicone	Off
Adjustment Tolerance	Auto

System - Adjust Volume

! Position	Isocenter
! Orientation	Transversal
! Rotation	0.00 deg
! A >> P	150 mm
! R >> L	150 mm
! F >> H	50 mm
Reset	Off

System - Tx/Rx

Frequency 1H	63.673772 MHz
Correction factor	1
Gain	High
Img. Scale Cor.	1.000
Reset	Off
? Ref. amplitude 1H	0.000 V

Physio - Signal1

1st Signal/Mode	ECG/Trigger
Average cycle	947 ± 122 ms
Average cycle	No Signal ms
Captured cycle	947 ± 122 ms
Acquisition window	830 ms
Trigger pulse	1
Trigger delay	0 ms
TR	824.00 ms
Concatenations	1
Segments	72
Phases	1
Adaptive Triggering	Off

Physio - Cardiac

Tagging	None
Magn. preparation	Non-sel. SR perf
T1	350 ms
T2 prep. duration 1	0 ms
T2 prep. duration 2	0 ms
T2 prep. duration 3	0 ms
T2 prep. duration 4	0 ms
T2 prep. duration 5	0 ms

Physio - Cardiac

T2 prep. duration 6	0 ms
T2 prep. duration 7	0 ms
T2 prep. duration 8	55 ms
T2 prep. duration 9	55 ms
Fat suppr.	None
Dark blood	Off
FoV read	360 mm
FoV phase	75.0 %
Phase resolution	75 %
Cine	Off
Trajectory	Cartesian
Dummy heartbeats	0
Motion Correction	Standard

Physio - PACE

Resp. control	Off
Concatenations	1

Inline - Common

Subtract	Off
Measurements	1
StdDev	Off
Motion Correction	Standard
Save original images	On

Inline - Cardiac

Inline Evaluation	T1 map
Magn. preparation	Non-sel. SR perf
Num. of preps	9
Sampling duration 1	1 beats
Sampling duration 2	1 beats
Sampling duration 3	1 beats
Sampling duration 4	1 beats
Sampling duration 5	1 beats
Sampling duration 6	1 beats
Sampling duration 7	1 beats
Sampling duration 8	1 beats
Sampling duration 9	1 beats
T2 prep. duration 1	0 ms
T2 prep. duration 2	0 ms
T2 prep. duration 3	0 ms
T2 prep. duration 4	0 ms
T2 prep. duration 5	0 ms
T2 prep. duration 6	0 ms
T2 prep. duration 7	0 ms
T2 prep. duration 8	55 ms
T2 prep. duration 9	55 ms
Contrasts	1
TE	1.26 ms
TR	824.00 ms
Recovery duration 1	2 beats
Recovery duration 2	2 beats
Recovery duration 3	0 beats
Recovery duration 4	0 beats
Recovery duration 5	0 beats
Recovery duration 6	0 beats
Recovery duration 7	0 beats
Recovery duration 8	3 beats
Recovery duration 9	3 beats
Motion Correction	Standard
Save original images	On

Inline - MIP

MIP-Sag	Off
---------	-----

Inline - MIP

MIP-Cor	Off
MIP-Tra	Off
MIP-Time	Off
Save original images	On

Inline - Composing

Inline Composing	Off
Distortion Corr.	On
Mode	2D
Unfiltered images	Off

Sequence - Part 1

Introduction	Off
Dimension	2D
Reordering	Linear
Asymmetric echo	Weak
Contrasts	1
Optimization	Min. TE
Multi-slice mode	Sequential
Echo spacing	2.9 ms
Sequence type	Trufi
Bandwidth	1085 Hz/Px

Sequence - Part 2

Define	Shots
Shots per slice	1
Segments	72
Trufi delta freq.	0 Hz
RF pulse type	Fast
Gradient mode	Fast
Excitation	Slice-sel.
Flip angle mode	Variable
Cine	Off

Sequence - Special

Partition coeff map	Off
ECV map	Off
Synth ECV map	Off
T1 array 0	350 ms
T1 array 1	350 ms
T1 array 2	350 ms
T1 array 3	350 ms
T1 array 4	350 ms
T1 array 5	350 ms
T1 array 6	350 ms
T1 array 7	350 ms
T1 array 8	350 ms
High-contrast acquisition	On
Error map	Off
SPAIR delay	5 ms
T1 sampling scheme	Fixed TS

Sequence - Assistant

Mode	Off
Allowed delay	10 s

\\CMIC\Physics\T2MES\T2MES_low_res\post FB T1T2

TA: 0:43 PM: FIX Voxel size: 1.4×1.4×8.0 mmPAT: 2 Rel. SNR: 1.00 : tfi

Properties

Prio recon	Off
Load images to viewer	On
Inline movie	Off
Auto store images	On
Load images to stamp segments	Off
Load images to graphic segments	On
Auto open inline display	Off
Auto close inline display	Off
Start measurement without further preparation	On
Wait for user to start	Off
Start measurements	Single measurement

Routine

Slice group	1
Slices	1
Dist. factor	20 %
Position	Isocenter
Orientation	Transversal
Phase enc. dir.	A >> P
AutoAlign	---
Phase oversampling	0 %
FoV read	360 mm
FoV phase	75.0 %
Slice thickness	8.0 mm
TR	904.00 ms
TE	1.26 ms
Averages	3
Concatenations	1
Filter	Distortion Corr.(2D)
Coil elements	BO1-3;SP2-4

Contrast - Common

TR	904.00 ms
TE	1.26 ms
Magn. preparation	Non-sel. SR perf
T1	300 ms
T2 prep. duration 1	0 ms
T2 prep. duration 2	0 ms
T2 prep. duration 3	0 ms
T2 prep. duration 4	0 ms
T2 prep. duration 5	0 ms
T2 prep. duration 6	0 ms
T2 prep. duration 7	0 ms
T2 prep. duration 8	55 ms
T2 prep. duration 9	55 ms
Flip angle	100 deg
Fat suppr.	None
Wrap-up Magn.	None

Contrast - Dynamic

Averages	3
Averaging mode	Short term
Reconstruction	Magnitude
Measurements	1
Multiple series	Off

Resolution - Common

FoV read	360 mm
FoV phase	75.0 %

Resolution - Common

Slice thickness	8.0 mm
Base resolution	256
Phase resolution	75 %
Phase partial Fourier	7/8
Trajectory	Cartesian
Interpolation	Off

Resolution - iPAT

PAT mode	GRAPPA
Accel. factor PE	2
Ref. lines PE	36
Reference scan mode	GRE/separate

Resolution - Filter Image

Image Filter	Off
Distortion Corr.	On
Mode	2D
Unfiltered images	Off
Prescan Normalize	Off
Normalize	Off
B1 filter	Off

Resolution - Filter Rawdata

Raw filter	Off
Elliptical filter	Off
POCS	Off

Geometry - Common

Slice group	1
Slices	1
Dist. factor	20 %
Position	Isocenter
Orientation	Transversal
Phase enc. dir.	A >> P
FoV read	360 mm
FoV phase	75.0 %
Slice thickness	8.0 mm
TR	904.00 ms
Multi-slice mode	Sequential
Series	Interleaved
Concatenations	1

Geometry - AutoAlign

Slice group	1
Position	Isocenter
Orientation	Transversal
Phase enc. dir.	A >> P
AutoAlign	---
Initial Position	Isocenter
L	0.0 mm
P	0.0 mm
H	0.0 mm
Initial Rotation	0.00 deg
Initial Orientation	Transversal

Geometry - Saturation

Fat suppr.	None
Wrap-up Magn.	None
Special sat.	None

Geometry - Navigator**System - Miscellaneous**

Positioning mode	FIX
Table position	H
Table position	0 mm
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Coil Combine Mode	Adaptive Combine
Save uncombined	Off
Matrix Optimization	Off
Coil Focus	Flat
AutoAlign	---
Coil Select Mode	Default

System - Adjustments

B0 Shim mode	Cardiac
Adjust with body coil	On
Confirm freq. adjustment	Off
Assume Dominant Fat	Off
Assume Silicone	Off
Adjustment Tolerance	Auto

System - Adjust Volume

! Position	Isocenter
! Orientation	Transversal
! Rotation	0.00 deg
! A >> P	150 mm
! R >> L	150 mm
! F >> H	50 mm
Reset	Off

System - Tx/Rx

Frequency 1H	63.673772 MHz
Correction factor	1
Gain	High
Img. Scale Cor.	1.000
Reset	Off
? Ref. amplitude 1H	0.000 V

Physio - Signal1

1st Signal/Mode	ECG/Trigger
Average cycle	947 ± 122 ms
Average cycle	No Signal ms
Captured cycle	947 ± 122 ms
Acquisition window	930 ms
Trigger pulse	1
Trigger delay	0 ms
TR	904.00 ms
Concatenations	1
Segments	72
Phases	1
Adaptive Triggering	Off

Physio - Cardiac

Tagging	None
Magn. preparation	Non-sel. SR perf
T1	300 ms
T2 prep. duration 1	0 ms
T2 prep. duration 2	0 ms
T2 prep. duration 3	0 ms
T2 prep. duration 4	0 ms
T2 prep. duration 5	0 ms

Physio - Cardiac

T2 prep. duration 6	0 ms
T2 prep. duration 7	0 ms
T2 prep. duration 8	55 ms
T2 prep. duration 9	55 ms
Fat suppr.	None
Dark blood	Off
FoV read	360 mm
FoV phase	75.0 %
Phase resolution	75 %
Cine	Off
Trajectory	Cartesian
Dummy heartbeats	0
Motion Correction	Standard

Physio - PACE

Resp. control	Off
Concatenations	1

Inline - Common

Subtract	Off
Measurements	1
StdDev	Off
Motion Correction	Standard
Save original images	On

Inline - Cardiac

Inline Evaluation	T1 map
Magn. preparation	Non-sel. SR perf
Num. of preps	9
Sampling duration 1	1 beats
Sampling duration 2	1 beats
Sampling duration 3	1 beats
Sampling duration 4	1 beats
Sampling duration 5	1 beats
Sampling duration 6	1 beats
Sampling duration 7	1 beats
Sampling duration 8	1 beats
Sampling duration 9	1 beats
T2 prep. duration 1	0 ms
T2 prep. duration 2	0 ms
T2 prep. duration 3	0 ms
T2 prep. duration 4	0 ms
T2 prep. duration 5	0 ms
T2 prep. duration 6	0 ms
T2 prep. duration 7	0 ms
T2 prep. duration 8	55 ms
T2 prep. duration 9	55 ms
Contrasts	1
TE	1.26 ms
TR	904.00 ms
Recovery duration 1	0 beats
Recovery duration 2	0 beats
Recovery duration 3	0 beats
Recovery duration 4	0 beats
Recovery duration 5	0 beats
Recovery duration 6	0 beats
Recovery duration 7	0 beats
Recovery duration 8	1 beats
Recovery duration 9	1 beats
Motion Correction	Standard
Save original images	On

Inline - MIP

MIP-Sag	Off
---------	-----

Inline - MIP

MIP-Cor	Off
MIP-Tra	Off
MIP-Time	Off
Save original images	On

Inline - Composing

Inline Composing	Off
Distortion Corr.	On
Mode	2D
Unfiltered images	Off

Sequence - Part 1

Introduction	Off
Dimension	2D
Reordering	Linear
Asymmetric echo	Weak
Contrasts	1
Optimization	Min. TE
Multi-slice mode	Sequential
Echo spacing	2.9 ms
Sequence type	Trufi
Bandwidth	1085 Hz/Px

Sequence - Part 2

Define	Shots
Shots per slice	1
Segments	72
Trufi delta freq.	0 Hz
RF pulse type	Fast
Gradient mode	Fast
Excitation	Slice-sel.
Flip angle mode	Variable
Cine	Off

Sequence - Special

Partition coeff map	Off
ECV map	Off
Synth ECV map	Off
T1 array 0	300 ms
T1 array 1	300 ms
T1 array 2	300 ms
T1 array 3	300 ms
T1 array 4	300 ms
T1 array 5	300 ms
T1 array 6	300 ms
T1 array 7	300 ms
T1 array 8	300 ms
High-contrast acquisition	On
Error map	Off
SPAIR delay	5 ms
T1 sampling scheme	Fixed TS

Sequence - Assistant

Mode	Off
Allowed delay	10 s

\\CMIC\Physics\T2MES\T2MES_low_res\BH_FW_4e_GT

TA: 8.5 s PM: FIX Voxel size: 1.4×1.4×8.0 mmPAT: Off Rel. SNR: 1.00 : fl

Properties

Prio recon	Off
Load images to viewer	Off
Inline movie	Off
Auto store images	On
Load images to stamp segments	Off
Load images to graphic segments	On
Auto open inline display	Off
Auto close inline display	Off
Start measurement without further preparation	On
Wait for user to start	Off
Start measurements	Single measurement

Routine

Slice group	1
Slices	1
Dist. factor	20 %
Position	Isocenter
Orientation	Transversal
Phase enc. dir.	A >> P
AutoAlign	---
Phase oversampling	0 %
FoV read	360 mm
FoV phase	75.0 %
Slice thickness	8.0 mm
TR	208.00 ms
TE 1	1.56 ms
TE 2	3.88 ms
TE 3	6.2 ms
TE 4	8.52 ms
Averages	1
Concatenations	1
Filter	Distortion Corr.(2D)
Coil elements	BO1-3;SP2-4

Contrast - Common

TR	208.00 ms
TE 1	1.56 ms
TE 2	3.88 ms
TE 3	6.2 ms
TE 4	8.52 ms
Magn. preparation	None
Flip angle	12 deg
Fat suppr.	None
Wrap-up Magn.	None

Contrast - Dynamic

Averages	1
Averaging mode	Short term
Reconstruction	Magnitude
Measurements	1
Multiple series	Off

Resolution - Common

FoV read	360 mm
FoV phase	75.0 %
Slice thickness	8.0 mm
Base resolution	256
Phase resolution	75 %
Phase partial Fourier	Off

Resolution - Common

Trajectory	Cartesian
Interpolation	Off

Resolution - iPAT

PAT mode	None
----------	------

Resolution - Filter Image

Image Filter	Off
Distortion Corr.	On
Mode	2D
Unfiltered images	Off
Prescan Normalize	Off
Normalize	Off
B1 filter	Off

Resolution - Filter Rawdata

Raw filter	Off
Elliptical filter	Off
POCS	Off

Geometry - Common

Slice group	1
Slices	1
Dist. factor	20 %
Position	Isocenter
Orientation	Transversal
Phase enc. dir.	A >> P
FoV read	360 mm
FoV phase	75.0 %
Slice thickness	8.0 mm
TR	208.00 ms
Multi-slice mode	Sequential
Series	Interleaved
Concatenations	1

Geometry - AutoAlign

Slice group	1
Position	Isocenter
Orientation	Transversal
Phase enc. dir.	A >> P
AutoAlign	---
Initial Position	Isocenter
L	0.0 mm
P	0.0 mm
H	0.0 mm
Initial Rotation	0.00 deg
Initial Orientation	Transversal

Geometry - Saturation

Fat suppr.	None
Wrap-up Magn.	None
Special sat.	None

Geometry - Navigator**System - Miscellaneous**

Positioning mode	FIX
Table position	H
Table position	0 mm

System - Miscellaneous

MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Coil Combine Mode	Adaptive Combine
Save uncombined	Off
Matrix Optimization	Off
Coil Focus	Flat
AutoAlign	---
Coil Select Mode	Default

System - Adjustments

B0 Shim mode	Cardiac
Adjust with body coil	On
Confirm freq. adjustment	Off
Assume Dominant Fat	Off
Assume Silicone	Off
Adjustment Tolerance	Auto

System - Adjust Volume

! Position	Isocenter
! Orientation	Transversal
! Rotation	0.00 deg
! A >> P	100 mm
! R >> L	100 mm
! F >> H	25 mm
Reset	Off

System - Tx/Rx

Frequency 1H	63.673772 MHz
Correction factor	1
Gain	High
Img. Scale Cor.	1.000
Reset	Off
? Ref. amplitude 1H	0.000 V

Physio - Signal1

1st Signal/Mode	ECG/Trigger
Average cycle	947 ± 122 ms
Average cycle	No Signal ms
Captured cycle	947 ± 122 ms
Acquisition window	830 ms
Trigger pulse	1
Trigger delay	622 ms
TR	208.00 ms
Concatenations	1
Segments	20
Phases	1
Adaptive Triggering	Off

Physio - Cardiac

Tagging	None
Magn. preparation	None
Fat suppr.	None
Dark blood	Off
FoV read	360 mm
FoV phase	75.0 %
Phase resolution	75 %
Cine	Off
Trajectory	Cartesian
Dummy heartbeats	1

Physio - PACE

Resp. control	Off
---------------	-----

Physio - PACE

Concatenations	1
----------------	---

Inline - Common

Subtract	Off
Measurements	1
StdDev	Off
Save original images	On

Inline - Cardiac

Inline Evaluation	Off
Magn. preparation	None
Contrasts	4
TE 1	1.56 ms
TE 2	3.88 ms
TE 3	6.2 ms
TE 4	8.52 ms
TR	208.00 ms
Save original images	On

Inline - MIP

MIP-Sag	Off
MIP-Cor	Off
MIP-Tra	Off
MIP-Time	Off
Save original images	On

Inline - Composing

Inline Composing	Off
Distortion Corr.	On
Mode	2D
Unfiltered images	Off

Sequence - Part 1

Introduction	Off
Dimension	2D
Reordering	Linear
Asymmetric echo	Off
Contrasts	4
Flow comp. 1	No
Readout mode	Monopolar
Optimization	Min. TE TR
Multi-slice mode	Sequential
Echo spacing	10.4 ms
Sequence type	Gre
Bandwidth 1	977 Hz/Px
Bandwidth 2	977 Hz/Px
Bandwidth 3	977 Hz/Px
Bandwidth 4	977 Hz/Px

Sequence - Part 2

Define	Segments
Segments	20
RF pulse type	Fast
Gradient mode	Fast
Excitation	Slice-sel.
Flip angle mode	Constant
RF spoiling	On
Phase Enc. Rewinder	On
Cine	Off

Sequence - Special

FatWater Separation	On
Multi-echo Images	On

Sequence - Special

In-Opp Phase Images	Off
Frequency Map	On
T2* Map	Off
Motion Correction	Off
MoCo Averaging Mode	Complex MoCo
MoCo Images Only?	Off
Frequency Offset	0 Hz
No. of Interleaves	0

Sequence - Assistant

Mode	Off
Allowed delay	0 s

\\CMIC\Physics\T2MES\T2MES_low_res\se te = 10

TA: 24:10 PM: FIX Voxel size: 1.9×1.9×8.0 mmRel. SNR: 1.00 : se

Properties

Prio recon	Off
Load images to viewer	Off
Inline movie	Off
Auto store images	On
Load images to stamp segments	Off
Load images to graphic segments	On
Auto open inline display	Off
Auto close inline display	Off
Start measurement without further preparation	On
Wait for user to start	Off
Start measurements	Single measurement

Routine

Slice group	1
Slices	1
Dist. factor	50 %
Position	Isocenter
Orientation	Transversal
Phase enc. dir.	A >> P
AutoAlign	---
Phase oversampling	0 %
FoV read	360 mm
FoV phase	75.0 %
Slice thickness	8.0 mm
TR	10000.0 ms
TE	10.0 ms
Averages	1
Concatenations	1
Filter	Distortion Corr.(2D)
Coil elements	BO1-3;SP2-4

Contrast - Common

TR	10000.0 ms
TE	10.0 ms
MTC	Off
Magn. preparation	None
Flip angle	90 deg
Fat suppr.	None
Water suppr.	None
Blood suppr.	Off

Contrast - Dynamic

Averages	1
Averaging mode	Short term
Reconstruction	Magnitude
Measurements	1
Multiple series	Each measurement

Resolution - Common

FoV read	360 mm
FoV phase	75.0 %
Slice thickness	8.0 mm
Base resolution	192
Phase resolution	100 %
Phase partial Fourier	Off
Interpolation	Off

Resolution - Filter Image

Image Filter	Off
--------------	-----

Resolution - Filter Image

Distortion Corr.	On
Mode	2D
Unfiltered images	Off
Prescan Normalize	Off
Normalize	Off
B1 filter	Off

Resolution - Filter Rawdata

Raw filter	Off
Elliptical filter	Off

Geometry - Common

Slice group	1
Slices	1
Dist. factor	50 %
Position	Isocenter
Orientation	Transversal
Phase enc. dir.	A >> P
FoV read	360 mm
FoV phase	75.0 %
Slice thickness	8.0 mm
TR	10000.0 ms
Multi-slice mode	Interleaved
Series	Interleaved
Concatenations	1

Geometry - AutoAlign

Slice group	1
Position	Isocenter
Orientation	Transversal
Phase enc. dir.	A >> P
AutoAlign	---
Initial Position	Isocenter
L	0.0 mm
P	0.0 mm
H	0.0 mm
Initial Rotation	0.00 deg
Initial Orientation	Transversal

Geometry - Saturation

Fat suppr.	None
Water suppr.	None
Special sat.	None

System - Miscellaneous

Positioning mode	FIX
Table position	H
Table position	0 mm
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Coil Combine Mode	Adaptive Combine
Save uncombined	Off
Matrix Optimization	Off
Coil Focus	Flat
AutoAlign	---
Coil Select Mode	Default

System - Adjustments

B0 Shim mode	Tune up
Adjust with body coil	On
Confirm freq. adjustment	Off
Assume Dominant Fat	Off
Assume Silicone	Off
Adjustment Tolerance	Auto

System - Adjust Volume

Position	Isocenter
Orientation	Transversal
Rotation	0.00 deg
A >> P	263 mm
R >> L	350 mm
F >> H	350 mm
Reset	Off

System - Tx/Rx

Frequency 1H	63.673772 MHz
Correction factor	1
Gain	High
Img. Scale Cor.	1.000
Reset	Off
? Ref. amplitude 1H	0.000 V

Physio - Signal1

1st Signal/Mode	None
TR	10000.0 ms
Concatenations	1

Physio - Cardiac

Magn. preparation	None
Fat suppr.	None
Dark blood	Off
FoV read	360 mm
FoV phase	75.0 %
Phase resolution	100 %

Inline - Common

Subtract	Off
Measurements	1
StdDev	Off
Save original images	On

Inline - MIP

MIP-Sag	Off
MIP-Cor	Off
MIP-Tra	Off
MIP-Time	Off
Save original images	On

Inline - Composing

Inline Composing	Off
Distortion Corr.	On
Mode	2D
Unfiltered images	Off

Sequence - Part 1

Introduction	Off
Asymmetric echo	Off
Contrasts	1
Multi-slice mode	Interleaved
Bandwidth	501 Hz/Px

Sequence - Part 2

Acoustic noise reduction	None
RF pulse type	Normal
Gradient mode	Fast

Sequence - Assistant

Mode	Off
Allowed delay	0 s

\\CMIC\Physics\T2MES\T2MES_low_res\se te = 20

TA: 24:10 PM: FIX Voxel size: 1.9×1.9×8.0 mmRel. SNR: 1.00 : se

Properties

Prio recon	Off
Load images to viewer	Off
Inline movie	Off
Auto store images	On
Load images to stamp segments	Off
Load images to graphic segments	On
Auto open inline display	Off
Auto close inline display	Off
Start measurement without further preparation	On
Wait for user to start	Off
Start measurements	Single measurement

Routine

Slice group	1
Slices	1
Dist. factor	50 %
Position	Isocenter
Orientation	Transversal
Phase enc. dir.	A >> P
AutoAlign	---
Phase oversampling	0 %
FoV read	360 mm
FoV phase	75.0 %
Slice thickness	8.0 mm
TR	10000.0 ms
TE	20.0 ms
Averages	1
Concatenations	1
Filter	Distortion Corr.(2D)
Coil elements	BO1-3;SP2-4

Contrast - Common

TR	10000.0 ms
TE	20.0 ms
MTC	Off
Magn. preparation	None
Flip angle	90 deg
Fat suppr.	None
Water suppr.	None
Blood suppr.	Off

Contrast - Dynamic

Averages	1
Averaging mode	Short term
Reconstruction	Magnitude
Measurements	1
Multiple series	Each measurement

Resolution - Common

FoV read	360 mm
FoV phase	75.0 %
Slice thickness	8.0 mm
Base resolution	192
Phase resolution	100 %
Phase partial Fourier	Off
Interpolation	Off

Resolution - Filter Image

Image Filter	Off
--------------	-----

Resolution - Filter Image

Distortion Corr.	On
Mode	2D
Unfiltered images	Off
Prescan Normalize	Off
Normalize	Off
B1 filter	Off

Resolution - Filter Rawdata

Raw filter	Off
Elliptical filter	Off

Geometry - Common

Slice group	1
Slices	1
Dist. factor	50 %
Position	Isocenter
Orientation	Transversal
Phase enc. dir.	A >> P
FoV read	360 mm
FoV phase	75.0 %
Slice thickness	8.0 mm
TR	10000.0 ms
Multi-slice mode	Interleaved
Series	Interleaved
Concatenations	1

Geometry - AutoAlign

Slice group	1
Position	Isocenter
Orientation	Transversal
Phase enc. dir.	A >> P
AutoAlign	---
Initial Position	Isocenter
L	0.0 mm
P	0.0 mm
H	0.0 mm
Initial Rotation	0.00 deg
Initial Orientation	Transversal

Geometry - Saturation

Fat suppr.	None
Water suppr.	None
Special sat.	None

System - Miscellaneous

Positioning mode	FIX
Table position	H
Table position	0 mm
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Coil Combine Mode	Adaptive Combine
Save uncombined	Off
Matrix Optimization	Off
Coil Focus	Flat
AutoAlign	---
Coil Select Mode	Default

System - Adjustments

B0 Shim mode	Standard
Adjust with body coil	On
Confirm freq. adjustment	Off
Assume Dominant Fat	Off
Assume Silicone	Off
Adjustment Tolerance	Auto

System - Adjust Volume

Position	Isocenter
Orientation	Transversal
Rotation	0.00 deg
A >> P	270 mm
R >> L	360 mm
F >> H	8 mm
Reset	Off

System - Tx/Rx

Frequency 1H	63.673772 MHz
Correction factor	1
Gain	High
Img. Scale Cor.	1.000
Reset	Off
? Ref. amplitude 1H	0.000 V

Physio - Signal1

1st Signal/Mode	None
TR	10000.0 ms
Concatenations	1

Physio - Cardiac

Magn. preparation	None
Fat suppr.	None
Dark blood	Off
FoV read	360 mm
FoV phase	75.0 %
Phase resolution	100 %

Inline - Common

Subtract	Off
Measurements	1
StdDev	Off
Save original images	On

Inline - MIP

MIP-Sag	Off
MIP-Cor	Off
MIP-Tra	Off
MIP-Time	Off
Save original images	On

Inline - Composing

Inline Composing	Off
Distortion Corr.	On
Mode	2D
Unfiltered images	Off

Sequence - Part 1

Introduction	Off
Asymmetric echo	Off
Contrasts	1
Multi-slice mode	Interleaved
Bandwidth	501 Hz/Px

Sequence - Part 2

Acoustic noise reduction	None
RF pulse type	Normal
Gradient mode	Fast

Sequence - Assistant

Mode	Off
Allowed delay	0 s

\\CMIC\Physics\T2MES\T2MES_low_res\se te = 30

TA: 24:10 PM: FIX Voxel size: 1.9×1.9×8.0 mmRel. SNR: 1.00 : se

Properties

Prio recon	Off
Load images to viewer	Off
Inline movie	Off
Auto store images	On
Load images to stamp segments	Off
Load images to graphic segments	On
Auto open inline display	Off
Auto close inline display	Off
Start measurement without further preparation	On
Wait for user to start	Off
Start measurements	Single measurement

Routine

Slice group	1
Slices	1
Dist. factor	50 %
Position	Isocenter
Orientation	Transversal
Phase enc. dir.	A >> P
AutoAlign	---
Phase oversampling	0 %
FoV read	360 mm
FoV phase	75.0 %
Slice thickness	8.0 mm
TR	10000.0 ms
TE	30.0 ms
Averages	1
Concatenations	1
Filter	Distortion Corr.(2D)
Coil elements	BO1-3;SP2-4

Contrast - Common

TR	10000.0 ms
TE	30.0 ms
MTC	Off
Magn. preparation	None
Flip angle	90 deg
Fat suppr.	None
Water suppr.	None
Blood suppr.	Off

Contrast - Dynamic

Averages	1
Averaging mode	Short term
Reconstruction	Magnitude
Measurements	1
Multiple series	Each measurement

Resolution - Common

FoV read	360 mm
FoV phase	75.0 %
Slice thickness	8.0 mm
Base resolution	192
Phase resolution	100 %
Phase partial Fourier	Off
Interpolation	Off

Resolution - Filter Image

Image Filter	Off
--------------	-----

Resolution - Filter Image

Distortion Corr.	On
Mode	2D
Unfiltered images	Off
Prescan Normalize	Off
Normalize	Off
B1 filter	Off

Resolution - Filter Rawdata

Raw filter	Off
Elliptical filter	Off

Geometry - Common

Slice group	1
Slices	1
Dist. factor	50 %
Position	Isocenter
Orientation	Transversal
Phase enc. dir.	A >> P
FoV read	360 mm
FoV phase	75.0 %
Slice thickness	8.0 mm
TR	10000.0 ms
Multi-slice mode	Interleaved
Series	Interleaved
Concatenations	1

Geometry - AutoAlign

Slice group	1
Position	Isocenter
Orientation	Transversal
Phase enc. dir.	A >> P
AutoAlign	---
Initial Position	Isocenter
L	0.0 mm
P	0.0 mm
H	0.0 mm
Initial Rotation	0.00 deg
Initial Orientation	Transversal

Geometry - Saturation

Fat suppr.	None
Water suppr.	None
Special sat.	None

System - Miscellaneous

Positioning mode	FIX
Table position	H
Table position	0 mm
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Coil Combine Mode	Adaptive Combine
Save uncombined	Off
Matrix Optimization	Off
Coil Focus	Flat
AutoAlign	---
Coil Select Mode	Default

System - Adjustments

B0 Shim mode	Standard
Adjust with body coil	On
Confirm freq. adjustment	Off
Assume Dominant Fat	Off
Assume Silicone	Off
Adjustment Tolerance	Auto

System - Adjust Volume

Position	Isocenter
Orientation	Transversal
Rotation	0.00 deg
A >> P	270 mm
R >> L	360 mm
F >> H	8 mm
Reset	Off

System - Tx/Rx

Frequency 1H	63.673772 MHz
Correction factor	1
Gain	High
Img. Scale Cor.	1.000
Reset	Off
? Ref. amplitude 1H	0.000 V

Physio - Signal1

1st Signal/Mode	None
TR	10000.0 ms
Concatenations	1

Physio - Cardiac

Magn. preparation	None
Fat suppr.	None
Dark blood	Off
FoV read	360 mm
FoV phase	75.0 %
Phase resolution	100 %

Inline - Common

Subtract	Off
Measurements	1
StdDev	Off
Save original images	On

Inline - MIP

MIP-Sag	Off
MIP-Cor	Off
MIP-Tra	Off
MIP-Time	Off
Save original images	On

Inline - Composing

Inline Composing	Off
Distortion Corr.	On
Mode	2D
Unfiltered images	Off

Sequence - Part 1

Introduction	Off
Asymmetric echo	Off
Contrasts	1
Multi-slice mode	Interleaved
Bandwidth	501 Hz/Px

Sequence - Part 2

Acoustic noise reduction	None
RF pulse type	Normal
Gradient mode	Fast

Sequence - Assistant

Mode	Off
Allowed delay	0 s

\\CMIC\Physics\T2MES\T2MES_low_res\se te = 40

TA: 24:10 PM: FIX Voxel size: 1.9×1.9×8.0 mmRel. SNR: 1.00 : se

Properties

Prio recon	Off
Load images to viewer	Off
Inline movie	Off
Auto store images	On
Load images to stamp segments	Off
Load images to graphic segments	On
Auto open inline display	Off
Auto close inline display	Off
Start measurement without further preparation	On
Wait for user to start	Off
Start measurements	Single measurement

Routine

Slice group	1
Slices	1
Dist. factor	50 %
Position	Isocenter
Orientation	Transversal
Phase enc. dir.	A >> P
AutoAlign	---
Phase oversampling	0 %
FoV read	360 mm
FoV phase	75.0 %
Slice thickness	8.0 mm
TR	10000.0 ms
TE	40.0 ms
Averages	1
Concatenations	1
Filter	Distortion Corr.(2D)
Coil elements	BO1-3;SP2-4

Contrast - Common

TR	10000.0 ms
TE	40.0 ms
MTC	Off
Magn. preparation	None
Flip angle	90 deg
Fat suppr.	None
Water suppr.	None
Blood suppr.	Off

Contrast - Dynamic

Averages	1
Averaging mode	Short term
Reconstruction	Magnitude
Measurements	1
Multiple series	Each measurement

Resolution - Common

FoV read	360 mm
FoV phase	75.0 %
Slice thickness	8.0 mm
Base resolution	192
Phase resolution	100 %
Phase partial Fourier	Off
Interpolation	Off

Resolution - Filter Image

Image Filter	Off
--------------	-----

Resolution - Filter Image

Distortion Corr.	On
Mode	2D
Unfiltered images	Off
Prescan Normalize	Off
Normalize	Off
B1 filter	Off

Resolution - Filter Rawdata

Raw filter	Off
Elliptical filter	Off

Geometry - Common

Slice group	1
Slices	1
Dist. factor	50 %
Position	Isocenter
Orientation	Transversal
Phase enc. dir.	A >> P
FoV read	360 mm
FoV phase	75.0 %
Slice thickness	8.0 mm
TR	10000.0 ms
Multi-slice mode	Interleaved
Series	Interleaved
Concatenations	1

Geometry - AutoAlign

Slice group	1
Position	Isocenter
Orientation	Transversal
Phase enc. dir.	A >> P
AutoAlign	---
Initial Position	Isocenter
L	0.0 mm
P	0.0 mm
H	0.0 mm
Initial Rotation	0.00 deg
Initial Orientation	Transversal

Geometry - Saturation

Fat suppr.	None
Water suppr.	None
Special sat.	None

System - Miscellaneous

Positioning mode	FIX
Table position	H
Table position	0 mm
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Coil Combine Mode	Adaptive Combine
Save uncombined	Off
Matrix Optimization	Off
Coil Focus	Flat
AutoAlign	---
Coil Select Mode	Default

System - Adjustments

B0 Shim mode	Tune up
Adjust with body coil	On
Confirm freq. adjustment	Off
Assume Dominant Fat	Off
Assume Silicone	Off
Adjustment Tolerance	Auto

System - Adjust Volume

Position	Isocenter
Orientation	Transversal
Rotation	0.00 deg
A >> P	263 mm
R >> L	350 mm
F >> H	350 mm
Reset	Off

System - Tx/Rx

Frequency 1H	63.673772 MHz
Correction factor	1
Gain	High
Img. Scale Cor.	1.000
Reset	Off
? Ref. amplitude 1H	0.000 V

Physio - Signal1

1st Signal/Mode	None
TR	10000.0 ms
Concatenations	1

Physio - Cardiac

Magn. preparation	None
Fat suppr.	None
Dark blood	Off
FoV read	360 mm
FoV phase	75.0 %
Phase resolution	100 %

Inline - Common

Subtract	Off
Measurements	1
StdDev	Off
Save original images	On

Inline - MIP

MIP-Sag	Off
MIP-Cor	Off
MIP-Tra	Off
MIP-Time	Off
Save original images	On

Inline - Composing

Inline Composing	Off
Distortion Corr.	On
Mode	2D
Unfiltered images	Off

Sequence - Part 1

Introduction	Off
Asymmetric echo	Off
Contrasts	1
Multi-slice mode	Interleaved
Bandwidth	501 Hz/Px

Sequence - Part 2

Acoustic noise reduction	None
RF pulse type	Normal
Gradient mode	Fast

Sequence - Assistant

Mode	Off
Allowed delay	0 s

\\CMIC\Physics\T2MES\T2MES_low_res\se te = 50

TA: 24:10 PM: FIX Voxel size: 1.9×1.9×8.0 mmRel. SNR: 1.00 : se

Properties

Prio recon	Off
Load images to viewer	Off
Inline movie	Off
Auto store images	On
Load images to stamp segments	Off
Load images to graphic segments	On
Auto open inline display	Off
Auto close inline display	Off
Start measurement without further preparation	On
Wait for user to start	Off
Start measurements	Single measurement

Routine

Slice group	1
Slices	1
Dist. factor	50 %
Position	Isocenter
Orientation	Transversal
Phase enc. dir.	A >> P
AutoAlign	---
Phase oversampling	0 %
FoV read	360 mm
FoV phase	75.0 %
Slice thickness	8.0 mm
TR	10000.0 ms
TE	50.0 ms
Averages	1
Concatenations	1
Filter	Distortion Corr.(2D)
Coil elements	BO1-3;SP2-4

Contrast - Common

TR	10000.0 ms
TE	50.0 ms
MTC	Off
Magn. preparation	None
Flip angle	90 deg
Fat suppr.	None
Water suppr.	None
Blood suppr.	Off

Contrast - Dynamic

Averages	1
Averaging mode	Short term
Reconstruction	Magnitude
Measurements	1
Multiple series	Each measurement

Resolution - Common

FoV read	360 mm
FoV phase	75.0 %
Slice thickness	8.0 mm
Base resolution	192
Phase resolution	100 %
Phase partial Fourier	Off
Interpolation	Off

Resolution - Filter Image

Image Filter	Off
--------------	-----

Resolution - Filter Image

Distortion Corr.	On
Mode	2D
Unfiltered images	Off
Prescan Normalize	Off
Normalize	Off
B1 filter	Off

Resolution - Filter Rawdata

Raw filter	Off
Elliptical filter	Off

Geometry - Common

Slice group	1
Slices	1
Dist. factor	50 %
Position	Isocenter
Orientation	Transversal
Phase enc. dir.	A >> P
FoV read	360 mm
FoV phase	75.0 %
Slice thickness	8.0 mm
TR	10000.0 ms
Multi-slice mode	Interleaved
Series	Interleaved
Concatenations	1

Geometry - AutoAlign

Slice group	1
Position	Isocenter
Orientation	Transversal
Phase enc. dir.	A >> P
AutoAlign	---
Initial Position	Isocenter
L	0.0 mm
P	0.0 mm
H	0.0 mm
Initial Rotation	0.00 deg
Initial Orientation	Transversal

Geometry - Saturation

Fat suppr.	None
Water suppr.	None
Special sat.	None

System - Miscellaneous

Positioning mode	FIX
Table position	H
Table position	0 mm
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Coil Combine Mode	Adaptive Combine
Save uncombined	Off
Matrix Optimization	Off
Coil Focus	Flat
AutoAlign	---
Coil Select Mode	Default

System - Adjustments

B0 Shim mode	Tune up
Adjust with body coil	On
Confirm freq. adjustment	Off
Assume Dominant Fat	Off
Assume Silicone	Off
Adjustment Tolerance	Auto

System - Adjust Volume

Position	Isocenter
Orientation	Transversal
Rotation	0.00 deg
A >> P	263 mm
R >> L	350 mm
F >> H	350 mm
Reset	Off

System - Tx/Rx

Frequency 1H	63.673772 MHz
Correction factor	1
Gain	High
Img. Scale Cor.	1.000
Reset	Off
? Ref. amplitude 1H	0.000 V

Physio - Signal1

1st Signal/Mode	None
TR	10000.0 ms
Concatenations	1

Physio - Cardiac

Magn. preparation	None
Fat suppr.	None
Dark blood	Off
FoV read	360 mm
FoV phase	75.0 %
Phase resolution	100 %

Inline - Common

Subtract	Off
Measurements	1
StdDev	Off
Save original images	On

Inline - MIP

MIP-Sag	Off
MIP-Cor	Off
MIP-Tra	Off
MIP-Time	Off
Save original images	On

Inline - Composing

Inline Composing	Off
Distortion Corr.	On
Mode	2D
Unfiltered images	Off

Sequence - Part 1

Introduction	Off
Asymmetric echo	Off
Contrasts	1
Multi-slice mode	Interleaved
Bandwidth	501 Hz/Px

Sequence - Part 2

Acoustic noise reduction	None
RF pulse type	Normal
Gradient mode	Fast

Sequence - Assistant

Mode	Off
Allowed delay	0 s

\\CMIC\Physics\T2MES\T2MES_low_res\se te = 60

TA: 24:10 PM: FIX Voxel size: 1.9×1.9×8.0 mmRel. SNR: 1.00 : se

Properties

Prio recon	Off
Load images to viewer	Off
Inline movie	Off
Auto store images	On
Load images to stamp segments	Off
Load images to graphic segments	On
Auto open inline display	Off
Auto close inline display	Off
Start measurement without further preparation	On
Wait for user to start	Off
Start measurements	Single measurement

Routine

Slice group	1
Slices	1
Dist. factor	50 %
Position	Isocenter
Orientation	Transversal
Phase enc. dir.	A >> P
AutoAlign	---
Phase oversampling	0 %
FoV read	360 mm
FoV phase	75.0 %
Slice thickness	8.0 mm
TR	10000.0 ms
TE	60.0 ms
Averages	1
Concatenations	1
Filter	Distortion Corr.(2D)
Coil elements	BO1-3;SP2-4

Contrast - Common

TR	10000.0 ms
TE	60.0 ms
MTC	Off
Magn. preparation	None
Flip angle	90 deg
Fat suppr.	None
Water suppr.	None
Blood suppr.	Off

Contrast - Dynamic

Averages	1
Averaging mode	Short term
Reconstruction	Magnitude
Measurements	1
Multiple series	Each measurement

Resolution - Common

FoV read	360 mm
FoV phase	75.0 %
Slice thickness	8.0 mm
Base resolution	192
Phase resolution	100 %
Phase partial Fourier	Off
Interpolation	Off

Resolution - Filter Image

Image Filter	Off
--------------	-----

Resolution - Filter Image

Distortion Corr.	On
Mode	2D
Unfiltered images	Off
Prescan Normalize	Off
Normalize	Off
B1 filter	Off

Resolution - Filter Rawdata

Raw filter	Off
Elliptical filter	Off

Geometry - Common

Slice group	1
Slices	1
Dist. factor	50 %
Position	Isocenter
Orientation	Transversal
Phase enc. dir.	A >> P
FoV read	360 mm
FoV phase	75.0 %
Slice thickness	8.0 mm
TR	10000.0 ms
Multi-slice mode	Interleaved
Series	Interleaved
Concatenations	1

Geometry - AutoAlign

Slice group	1
Position	Isocenter
Orientation	Transversal
Phase enc. dir.	A >> P
AutoAlign	---
Initial Position	Isocenter
L	0.0 mm
P	0.0 mm
H	0.0 mm
Initial Rotation	0.00 deg
Initial Orientation	Transversal

Geometry - Saturation

Fat suppr.	None
Water suppr.	None
Special sat.	None

System - Miscellaneous

Positioning mode	FIX
Table position	H
Table position	0 mm
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Coil Combine Mode	Adaptive Combine
Save uncombined	Off
Matrix Optimization	Off
Coil Focus	Flat
AutoAlign	---
Coil Select Mode	Default

System - Adjustments

B0 Shim mode	Tune up
Adjust with body coil	On
Confirm freq. adjustment	Off
Assume Dominant Fat	Off
Assume Silicone	Off
Adjustment Tolerance	Auto

System - Adjust Volume

Position	Isocenter
Orientation	Transversal
Rotation	0.00 deg
A >> P	263 mm
R >> L	350 mm
F >> H	350 mm
Reset	Off

System - Tx/Rx

Frequency 1H	63.673772 MHz
Correction factor	1
Gain	High
Img. Scale Cor.	1.000
Reset	Off
? Ref. amplitude 1H	0.000 V

Physio - Signal1

1st Signal/Mode	None
TR	10000.0 ms
Concatenations	1

Physio - Cardiac

Magn. preparation	None
Fat suppr.	None
Dark blood	Off
FoV read	360 mm
FoV phase	75.0 %
Phase resolution	100 %

Inline - Common

Subtract	Off
Measurements	1
StdDev	Off
Save original images	On

Inline - MIP

MIP-Sag	Off
MIP-Cor	Off
MIP-Tra	Off
MIP-Time	Off
Save original images	On

Inline - Composing

Inline Composing	Off
Distortion Corr.	On
Mode	2D
Unfiltered images	Off

Sequence - Part 1

Introduction	Off
Asymmetric echo	Off
Contrasts	1
Multi-slice mode	Interleaved
Bandwidth	501 Hz/Px

Sequence - Part 2

Acoustic noise reduction	None
RF pulse type	Normal
Gradient mode	Fast

Sequence - Assistant

Mode	Off
Allowed delay	0 s

\\CMIC\Physics\T2MES\T2MES_low_res\se te = 80

TA: 24:10 PM: FIX Voxel size: 1.9×1.9×8.0 mmRel. SNR: 1.00 : se

Properties

Prio recon	Off
Load images to viewer	Off
Inline movie	Off
Auto store images	On
Load images to stamp segments	Off
Load images to graphic segments	On
Auto open inline display	Off
Auto close inline display	Off
Start measurement without further preparation	On
Wait for user to start	Off
Start measurements	Single measurement

Routine

Slice group	1
Slices	1
Dist. factor	50 %
Position	Isocenter
Orientation	Transversal
Phase enc. dir.	A >> P
AutoAlign	---
Phase oversampling	0 %
FoV read	360 mm
FoV phase	75.0 %
Slice thickness	8.0 mm
TR	10000.0 ms
TE	80.0 ms
Averages	1
Concatenations	1
Filter	Distortion Corr.(2D)
Coil elements	BO1-3;SP2-4

Contrast - Common

TR	10000.0 ms
TE	80.0 ms
MTC	Off
Magn. preparation	None
Flip angle	90 deg
Fat suppr.	None
Water suppr.	None
Blood suppr.	Off

Contrast - Dynamic

Averages	1
Averaging mode	Short term
Reconstruction	Magnitude
Measurements	1
Multiple series	Each measurement

Resolution - Common

FoV read	360 mm
FoV phase	75.0 %
Slice thickness	8.0 mm
Base resolution	192
Phase resolution	100 %
Phase partial Fourier	Off
Interpolation	Off

Resolution - Filter Image

Image Filter	Off
--------------	-----

Resolution - Filter Image

Distortion Corr.	On
Mode	2D
Unfiltered images	Off
Prescan Normalize	Off
Normalize	Off
B1 filter	Off

Resolution - Filter Rawdata

Raw filter	Off
Elliptical filter	Off

Geometry - Common

Slice group	1
Slices	1
Dist. factor	50 %
Position	Isocenter
Orientation	Transversal
Phase enc. dir.	A >> P
FoV read	360 mm
FoV phase	75.0 %
Slice thickness	8.0 mm
TR	10000.0 ms
Multi-slice mode	Interleaved
Series	Interleaved
Concatenations	1

Geometry - AutoAlign

Slice group	1
Position	Isocenter
Orientation	Transversal
Phase enc. dir.	A >> P
AutoAlign	---
Initial Position	Isocenter
L	0.0 mm
P	0.0 mm
H	0.0 mm
Initial Rotation	0.00 deg
Initial Orientation	Transversal

Geometry - Saturation

Fat suppr.	None
Water suppr.	None
Special sat.	None

System - Miscellaneous

Positioning mode	FIX
Table position	H
Table position	0 mm
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Coil Combine Mode	Adaptive Combine
Save uncombined	Off
Matrix Optimization	Off
Coil Focus	Flat
AutoAlign	---
Coil Select Mode	Default

System - Adjustments

B0 Shim mode	Tune up
Adjust with body coil	On
Confirm freq. adjustment	Off
Assume Dominant Fat	Off
Assume Silicone	Off
Adjustment Tolerance	Auto

System - Adjust Volume

Position	Isocenter
Orientation	Transversal
Rotation	0.00 deg
A >> P	263 mm
R >> L	350 mm
F >> H	350 mm
Reset	Off

System - Tx/Rx

Frequency 1H	63.673772 MHz
Correction factor	1
Gain	High
Img. Scale Cor.	1.000
Reset	Off
? Ref. amplitude 1H	0.000 V

Physio - Signal1

1st Signal/Mode	None
TR	10000.0 ms
Concatenations	1

Physio - Cardiac

Magn. preparation	None
Fat suppr.	None
Dark blood	Off
FoV read	360 mm
FoV phase	75.0 %
Phase resolution	100 %

Inline - Common

Subtract	Off
Measurements	1
StdDev	Off
Save original images	On

Inline - MIP

MIP-Sag	Off
MIP-Cor	Off
MIP-Tra	Off
MIP-Time	Off
Save original images	On

Inline - Composing

Inline Composing	Off
Distortion Corr.	On
Mode	2D
Unfiltered images	Off

Sequence - Part 1

Introduction	Off
Asymmetric echo	Off
Contrasts	1
Multi-slice mode	Interleaved
Bandwidth	501 Hz/Px

Sequence - Part 2

Acoustic noise reduction	None
RF pulse type	Normal
Gradient mode	Fast

Sequence - Assistant

Mode	Off
Allowed delay	0 s

\\CMIC\Physics\T2MES\T2MES_low_res\se te = 100

TA: 24:10 PM: FIX Voxel size: 1.9×1.9×8.0 mmRel. SNR: 1.00 : se

Properties

Prio recon	Off
Load images to viewer	Off
Inline movie	Off
Auto store images	On
Load images to stamp segments	Off
Load images to graphic segments	On
Auto open inline display	Off
Auto close inline display	Off
Start measurement without further preparation	On
Wait for user to start	Off
Start measurements	Single measurement

Routine

Slice group	1
Slices	1
Dist. factor	50 %
Position	Isocenter
Orientation	Transversal
Phase enc. dir.	A >> P
AutoAlign	---
Phase oversampling	0 %
FoV read	360 mm
FoV phase	75.0 %
Slice thickness	8.0 mm
TR	10000.0 ms
TE	100.0 ms
Averages	1
Concatenations	1
Filter	Distortion Corr.(2D)
Coil elements	BO1-3;SP2-4

Contrast - Common

TR	10000.0 ms
TE	100.0 ms
MTC	Off
Magn. preparation	None
Flip angle	90 deg
Fat suppr.	None
Water suppr.	None
Blood suppr.	Off

Contrast - Dynamic

Averages	1
Averaging mode	Short term
Reconstruction	Magnitude
Measurements	1
Multiple series	Each measurement

Resolution - Common

FoV read	360 mm
FoV phase	75.0 %
Slice thickness	8.0 mm
Base resolution	192
Phase resolution	100 %
Phase partial Fourier	Off
Interpolation	Off

Resolution - Filter Image

Image Filter	Off
--------------	-----

Resolution - Filter Image

Distortion Corr.	On
Mode	2D
Unfiltered images	Off
Prescan Normalize	Off
Normalize	Off
B1 filter	Off

Resolution - Filter Rawdata

Raw filter	Off
Elliptical filter	Off

Geometry - Common

Slice group	1
Slices	1
Dist. factor	50 %
Position	Isocenter
Orientation	Transversal
Phase enc. dir.	A >> P
FoV read	360 mm
FoV phase	75.0 %
Slice thickness	8.0 mm
TR	10000.0 ms
Multi-slice mode	Interleaved
Series	Interleaved
Concatenations	1

Geometry - AutoAlign

Slice group	1
Position	Isocenter
Orientation	Transversal
Phase enc. dir.	A >> P
AutoAlign	---
Initial Position	Isocenter
L	0.0 mm
P	0.0 mm
H	0.0 mm
Initial Rotation	0.00 deg
Initial Orientation	Transversal

Geometry - Saturation

Fat suppr.	None
Water suppr.	None
Special sat.	None

System - Miscellaneous

Positioning mode	FIX
Table position	H
Table position	0 mm
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Coil Combine Mode	Adaptive Combine
Save uncombined	Off
Matrix Optimization	Off
Coil Focus	Flat
AutoAlign	---
Coil Select Mode	Default

System - Adjustments

B0 Shim mode	Tune up
Adjust with body coil	On
Confirm freq. adjustment	Off
Assume Dominant Fat	Off
Assume Silicone	Off
Adjustment Tolerance	Auto

System - Adjust Volume

Position	Isocenter
Orientation	Transversal
Rotation	0.00 deg
A >> P	263 mm
R >> L	350 mm
F >> H	350 mm
Reset	Off

System - Tx/Rx

Frequency 1H	63.673772 MHz
Correction factor	1
Gain	High
Img. Scale Cor.	1.000
Reset	Off
? Ref. amplitude 1H	0.000 V

Physio - Signal1

1st Signal/Mode	None
TR	10000.0 ms
Concatenations	1

Physio - Cardiac

Magn. preparation	None
Fat suppr.	None
Dark blood	Off
FoV read	360 mm
FoV phase	75.0 %
Phase resolution	100 %

Inline - Common

Subtract	Off
Measurements	1
StdDev	Off
Save original images	On

Inline - MIP

MIP-Sag	Off
MIP-Cor	Off
MIP-Tra	Off
MIP-Time	Off
Save original images	On

Inline - Composing

Inline Composing	Off
Distortion Corr.	On
Mode	2D
Unfiltered images	Off

Sequence - Part 1

Introduction	Off
Asymmetric echo	Off
Contrasts	1
Multi-slice mode	Interleaved
Bandwidth	501 Hz/Px

Sequence - Part 2

Acoustic noise reduction	None
RF pulse type	Normal
Gradient mode	Fast

Sequence - Assistant

Mode	Off
Allowed delay	0 s

\\CMIC\Physics\T2MES\T2MES_low_res\se te = 125

TA: 24:10 PM: FIX Voxel size: 1.9×1.9×8.0 mmRel. SNR: 1.00 : se

Properties

Prio recon	Off
Load images to viewer	Off
Inline movie	Off
Auto store images	On
Load images to stamp segments	Off
Load images to graphic segments	On
Auto open inline display	Off
Auto close inline display	Off
Start measurement without further preparation	On
Wait for user to start	Off
Start measurements	Single measurement

Routine

Slice group	1
Slices	1
Dist. factor	50 %
Position	Isocenter
Orientation	Transversal
Phase enc. dir.	A >> P
AutoAlign	---
Phase oversampling	0 %
FoV read	360 mm
FoV phase	75.0 %
Slice thickness	8.0 mm
TR	10000.0 ms
TE	125.0 ms
Averages	1
Concatenations	1
Filter	Distortion Corr.(2D)
Coil elements	BO1-3;SP2-4

Contrast - Common

TR	10000.0 ms
TE	125.0 ms
MTC	Off
Magn. preparation	None
Flip angle	90 deg
Fat suppr.	None
Water suppr.	None
Blood suppr.	Off

Contrast - Dynamic

Averages	1
Averaging mode	Short term
Reconstruction	Magnitude
Measurements	1
Multiple series	Each measurement

Resolution - Common

FoV read	360 mm
FoV phase	75.0 %
Slice thickness	8.0 mm
Base resolution	192
Phase resolution	100 %
Phase partial Fourier	Off
Interpolation	Off

Resolution - Filter Image

Image Filter	Off
--------------	-----

Resolution - Filter Image

Distortion Corr.	On
Mode	2D
Unfiltered images	Off
Prescan Normalize	Off
Normalize	Off
B1 filter	Off

Resolution - Filter Rawdata

Raw filter	Off
Elliptical filter	Off

Geometry - Common

Slice group	1
Slices	1
Dist. factor	50 %
Position	Isocenter
Orientation	Transversal
Phase enc. dir.	A >> P
FoV read	360 mm
FoV phase	75.0 %
Slice thickness	8.0 mm
TR	10000.0 ms
Multi-slice mode	Interleaved
Series	Interleaved
Concatenations	1

Geometry - AutoAlign

Slice group	1
Position	Isocenter
Orientation	Transversal
Phase enc. dir.	A >> P
AutoAlign	---
Initial Position	Isocenter
L	0.0 mm
P	0.0 mm
H	0.0 mm
Initial Rotation	0.00 deg
Initial Orientation	Transversal

Geometry - Saturation

Fat suppr.	None
Water suppr.	None
Special sat.	None

System - Miscellaneous

Positioning mode	FIX
Table position	H
Table position	0 mm
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Coil Combine Mode	Adaptive Combine
Save uncombined	Off
Matrix Optimization	Off
Coil Focus	Flat
AutoAlign	---
Coil Select Mode	Default

System - Adjustments

B0 Shim mode	Tune up
Adjust with body coil	On
Confirm freq. adjustment	Off
Assume Dominant Fat	Off
Assume Silicone	Off
Adjustment Tolerance	Auto

System - Adjust Volume

Position	Isocenter
Orientation	Transversal
Rotation	0.00 deg
A >> P	263 mm
R >> L	350 mm
F >> H	350 mm
Reset	Off

System - Tx/Rx

Frequency 1H	63.673772 MHz
Correction factor	1
Gain	High
Img. Scale Cor.	1.000
Reset	Off
? Ref. amplitude 1H	0.000 V

Physio - Signal1

1st Signal/Mode	None
TR	10000.0 ms
Concatenations	1

Physio - Cardiac

Magn. preparation	None
Fat suppr.	None
Dark blood	Off
FoV read	360 mm
FoV phase	75.0 %
Phase resolution	100 %

Inline - Common

Subtract	Off
Measurements	1
StdDev	Off
Save original images	On

Inline - MIP

MIP-Sag	Off
MIP-Cor	Off
MIP-Tra	Off
MIP-Time	Off
Save original images	On

Inline - Composing

Inline Composing	Off
Distortion Corr.	On
Mode	2D
Unfiltered images	Off

Sequence - Part 1

Introduction	Off
Asymmetric echo	Off
Contrasts	1
Multi-slice mode	Interleaved
Bandwidth	501 Hz/Px

Sequence - Part 2

Acoustic noise reduction	None
RF pulse type	Normal
Gradient mode	Fast

Sequence - Assistant

Mode	Off
Allowed delay	0 s

\\CMIC\Physics\T2MES\T2MES_low_res\se te = 150

TA: 24:10 PM: FIX Voxel size: 1.9×1.9×8.0 mmRel. SNR: 1.00 : se

Properties

Prio recon	Off
Load images to viewer	Off
Inline movie	Off
Auto store images	On
Load images to stamp segments	Off
Load images to graphic segments	On
Auto open inline display	Off
Auto close inline display	Off
Start measurement without further preparation	On
Wait for user to start	Off
Start measurements	Single measurement

Routine

Slice group	1
Slices	1
Dist. factor	50 %
Position	Isocenter
Orientation	Transversal
Phase enc. dir.	A >> P
AutoAlign	---
Phase oversampling	0 %
FoV read	360 mm
FoV phase	75.0 %
Slice thickness	8.0 mm
TR	10000.0 ms
TE	150.0 ms
Averages	1
Concatenations	1
Filter	Distortion Corr.(2D)
Coil elements	BO1-3;SP2-4

Contrast - Common

TR	10000.0 ms
TE	150.0 ms
MTC	Off
Magn. preparation	None
Flip angle	90 deg
Fat suppr.	None
Water suppr.	None
Blood suppr.	Off

Contrast - Dynamic

Averages	1
Averaging mode	Short term
Reconstruction	Magnitude
Measurements	1
Multiple series	Each measurement

Resolution - Common

FoV read	360 mm
FoV phase	75.0 %
Slice thickness	8.0 mm
Base resolution	192
Phase resolution	100 %
Phase partial Fourier	Off
Interpolation	Off

Resolution - Filter Image

Image Filter	Off
--------------	-----

Resolution - Filter Image

Distortion Corr.	On
Mode	2D
Unfiltered images	Off
Prescan Normalize	Off
Normalize	Off
B1 filter	Off

Resolution - Filter Rawdata

Raw filter	Off
Elliptical filter	Off

Geometry - Common

Slice group	1
Slices	1
Dist. factor	50 %
Position	Isocenter
Orientation	Transversal
Phase enc. dir.	A >> P
FoV read	360 mm
FoV phase	75.0 %
Slice thickness	8.0 mm
TR	10000.0 ms
Multi-slice mode	Interleaved
Series	Interleaved
Concatenations	1

Geometry - AutoAlign

Slice group	1
Position	Isocenter
Orientation	Transversal
Phase enc. dir.	A >> P
AutoAlign	---
Initial Position	Isocenter
L	0.0 mm
P	0.0 mm
H	0.0 mm
Initial Rotation	0.00 deg
Initial Orientation	Transversal

Geometry - Saturation

Fat suppr.	None
Water suppr.	None
Special sat.	None

System - Miscellaneous

Positioning mode	FIX
Table position	H
Table position	0 mm
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Coil Combine Mode	Adaptive Combine
Save uncombined	Off
Matrix Optimization	Off
Coil Focus	Flat
AutoAlign	---
Coil Select Mode	Default

System - Adjustments

B0 Shim mode	Tune up
Adjust with body coil	On
Confirm freq. adjustment	Off
Assume Dominant Fat	Off
Assume Silicone	Off
Adjustment Tolerance	Auto

System - Adjust Volume

Position	Isocenter
Orientation	Transversal
Rotation	0.00 deg
A >> P	263 mm
R >> L	350 mm
F >> H	350 mm
Reset	Off

System - Tx/Rx

Frequency 1H	63.673772 MHz
Correction factor	1
Gain	High
Img. Scale Cor.	1.000
Reset	Off
? Ref. amplitude 1H	0.000 V

Physio - Signal1

1st Signal/Mode	None
TR	10000.0 ms
Concatenations	1

Physio - Cardiac

Magn. preparation	None
Fat suppr.	None
Dark blood	Off
FoV read	360 mm
FoV phase	75.0 %
Phase resolution	100 %

Inline - Common

Subtract	Off
Measurements	1
StdDev	Off
Save original images	On

Inline - MIP

MIP-Sag	Off
MIP-Cor	Off
MIP-Tra	Off
MIP-Time	Off
Save original images	On

Inline - Composing

Inline Composing	Off
Distortion Corr.	On
Mode	2D
Unfiltered images	Off

Sequence - Part 1

Introduction	Off
Asymmetric echo	Off
Contrasts	1
Multi-slice mode	Interleaved
Bandwidth	501 Hz/Px

Sequence - Part 2

Acoustic noise reduction	None
RF pulse type	Normal
Gradient mode	Fast

Sequence - Assistant

Mode	Off
Allowed delay	0 s

\\CMIC\Physics\T2MES\T2MES_low_res\IR GRE TI=20

TA: 24:00 PM: FIX Voxel size: 1.9×1.9×8.0 mmRel. SNR: 1.00 : fl

Properties

Prio recon	Off
Load images to viewer	On
Inline movie	Off
Auto store images	On
Load images to stamp segments	Off
Load images to graphic segments	Off
Auto open inline display	Off
Auto close inline display	Off
Start measurement without further preparation	On
Wait for user to start	Off
Start measurements	Single measurement

Routine

Slice group	1
Slices	1
Position	Isocenter
Orientation	Transversal
Phase enc. dir.	A >> P
AutoAlign	---
FoV read	360 mm
FoV phase	75.0 %
Slice thickness	8.0 mm
TR	10000.00 ms
TE	10.00 ms
Filter	Distortion Corr.(2D)
Coil elements	BO1-3;SP2-4

Contrast - Common

TR	10000.00 ms
TE	10.00 ms
Magn. preparation	Non-sel. IR
TI	20 ms
Flip angle	90 deg

Resolution - Common

FoV read	360 mm
FoV phase	75.0 %
Slice thickness	8.0 mm
Base resolution	192
Phase resolution	100 %

Resolution - Filter Image

Image Filter	Off
Distortion Corr.	On
Mode	2D
Unfiltered images	Off
Prescan Normalize	Off
Normalize	Off
B1 filter	Off

Resolution - Filter Rawdata

Raw filter	Off
Elliptical filter	Off

Geometry - Common

Slice group	1
Slices	1
Position	Isocenter
Orientation	Transversal

Geometry - Common

Phase enc. dir.	A >> P
FoV read	360 mm
FoV phase	75.0 %
Slice thickness	8.0 mm
TR	10000.00 ms

Geometry - AutoAlign

Slice group	1
Position	Isocenter
Orientation	Transversal
Phase enc. dir.	A >> P
AutoAlign	---
Initial Position	Isocenter
L	0.0 mm
P	0.0 mm
H	0.0 mm
Initial Rotation	0.00 deg
Initial Orientation	Transversal

System - Miscellaneous

Positioning mode	FIX
Table position	H
Table position	0 mm
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Coil Combine Mode	Adaptive Combine
Save uncombined	Off
Matrix Optimization	Off
Coil Focus	Flat
AutoAlign	---
Coil Select Mode	Default

System - Adjustments

B0 Shim mode	Cardiac
Adjust with body coil	On
Confirm freq. adjustment	Off
Assume Dominant Fat	Off
Assume Silicone	Off
Adjustment Tolerance	Auto

System - Adjust Volume

! Position	Isocenter
! Orientation	Transversal
! Rotation	0.00 deg
! A >> P	100 mm
! R >> L	100 mm
! F >> H	50 mm
Reset	Off

System - Tx/Rx

Frequency 1H	63.673772 MHz
Correction factor	1
Gain	High
Img. Scale Cor.	1.000
Reset	Off
? Ref. amplitude 1H	0.000 V

Sequence - Part 1

Bandwidth	500 Hz/Px
-----------	-----------

Sequence - Special

Adiab IR Pulse	Tan/Tanh
RF Phase Spoiling	ON

Sequence - Assistant

Mode	Off
------	-----

\\CMIC\Physics\T2MES\T2MES_low_res\IR GRE TI=50

TA: 24:00 PM: FIX Voxel size: 1.9×1.9×8.0 mmRel. SNR: 1.00 : fl

Properties

Prio recon	Off
Load images to viewer	On
Inline movie	Off
Auto store images	On
Load images to stamp segments	Off
Load images to graphic segments	Off
Auto open inline display	Off
Auto close inline display	Off
Start measurement without further preparation	On
Wait for user to start	Off
Start measurements	Single measurement

Routine

Slice group	1
Slices	1
Position	Isocenter
Orientation	Transversal
Phase enc. dir.	A >> P
AutoAlign	---
FoV read	360 mm
FoV phase	75.0 %
Slice thickness	8.0 mm
TR	10000.00 ms
TE	10.00 ms
Filter	Distortion Corr.(2D)
Coil elements	BO1-3;SP2-4

Contrast - Common

TR	10000.00 ms
TE	10.00 ms
Magn. preparation	Non-sel. IR
TI	50 ms
Flip angle	90 deg

Resolution - Common

FoV read	360 mm
FoV phase	75.0 %
Slice thickness	8.0 mm
Base resolution	192
Phase resolution	100 %

Resolution - Filter Image

Image Filter	Off
Distortion Corr.	On
Mode	2D
Unfiltered images	Off
Prescan Normalize	Off
Normalize	Off
B1 filter	Off

Resolution - Filter Rawdata

Raw filter	Off
Elliptical filter	Off

Geometry - Common

Slice group	1
Slices	1
Position	Isocenter
Orientation	Transversal

Geometry - Common

Phase enc. dir.	A >> P
FoV read	360 mm
FoV phase	75.0 %
Slice thickness	8.0 mm
TR	10000.00 ms

Geometry - AutoAlign

Slice group	1
Position	Isocenter
Orientation	Transversal
Phase enc. dir.	A >> P
AutoAlign	---
Initial Position	Isocenter
L	0.0 mm
P	0.0 mm
H	0.0 mm
Initial Rotation	0.00 deg
Initial Orientation	Transversal

System - Miscellaneous

Positioning mode	FIX
Table position	H
Table position	0 mm
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Coil Combine Mode	Adaptive Combine
Save uncombined	Off
Matrix Optimization	Off
Coil Focus	Flat
AutoAlign	---
Coil Select Mode	Default

System - Adjustments

B0 Shim mode	Cardiac
Adjust with body coil	On
Confirm freq. adjustment	Off
Assume Dominant Fat	Off
Assume Silicone	Off
Adjustment Tolerance	Auto

System - Adjust Volume

! Position	Isocenter
! Orientation	Transversal
! Rotation	0.00 deg
! A >> P	100 mm
! R >> L	100 mm
! F >> H	50 mm
Reset	Off

System - Tx/Rx

Frequency 1H	63.673772 MHz
Correction factor	1
Gain	High
Img. Scale Cor.	1.000
Reset	Off
? Ref. amplitude 1H	0.000 V

Sequence - Part 1

Bandwidth	500 Hz/Px
-----------	-----------

Sequence - Special

Adiab IR Pulse	Tan/Tanh
RF Phase Spoiling	ON

Sequence - Assistant

Mode	Off
------	-----

\\CMIC\Physics\T2MES\T2MES_low_res\IR GRE TI=100

TA: 24:00 PM: FIX Voxel size: 1.9×1.9×8.0 mmRel. SNR: 1.00 : fl

Properties

Prio recon	Off
Load images to viewer	On
Inline movie	Off
Auto store images	On
Load images to stamp segments	Off
Load images to graphic segments	Off
Auto open inline display	Off
Auto close inline display	Off
Start measurement without further preparation	On
Wait for user to start	Off
Start measurements	Single measurement

Routine

Slice group	1
Slices	1
Position	Isocenter
Orientation	Transversal
Phase enc. dir.	A >> P
AutoAlign	---
FoV read	360 mm
FoV phase	75.0 %
Slice thickness	8.0 mm
TR	10000.00 ms
TE	10.00 ms
Filter	Distortion Corr.(2D)
Coil elements	BO1-3;SP2-4

Contrast - Common

TR	10000.00 ms
TE	10.00 ms
Magn. preparation	Non-sel. IR
TI	100 ms
Flip angle	90 deg

Resolution - Common

FoV read	360 mm
FoV phase	75.0 %
Slice thickness	8.0 mm
Base resolution	192
Phase resolution	100 %

Resolution - Filter Image

Image Filter	Off
Distortion Corr.	On
Mode	2D
Unfiltered images	Off
Prescan Normalize	Off
Normalize	Off
B1 filter	Off

Resolution - Filter Rawdata

Raw filter	Off
Elliptical filter	Off

Geometry - Common

Slice group	1
Slices	1
Position	Isocenter
Orientation	Transversal

Geometry - Common

Phase enc. dir.	A >> P
FoV read	360 mm
FoV phase	75.0 %
Slice thickness	8.0 mm
TR	10000.00 ms

Geometry - AutoAlign

Slice group	1
Position	Isocenter
Orientation	Transversal
Phase enc. dir.	A >> P
AutoAlign	---
Initial Position	Isocenter
L	0.0 mm
P	0.0 mm
H	0.0 mm
Initial Rotation	0.00 deg
Initial Orientation	Transversal

System - Miscellaneous

Positioning mode	FIX
Table position	H
Table position	0 mm
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Coil Combine Mode	Adaptive Combine
Save uncombined	Off
Matrix Optimization	Off
Coil Focus	Flat
AutoAlign	---
Coil Select Mode	Default

System - Adjustments

B0 Shim mode	Cardiac
Adjust with body coil	On
Confirm freq. adjustment	Off
Assume Dominant Fat	Off
Assume Silicone	Off
Adjustment Tolerance	Auto

System - Adjust Volume

! Position	Isocenter
! Orientation	Transversal
! Rotation	0.00 deg
! A >> P	100 mm
! R >> L	100 mm
! F >> H	50 mm
Reset	Off

System - Tx/Rx

Frequency 1H	63.673772 MHz
Correction factor	1
Gain	High
Img. Scale Cor.	1.000
Reset	Off
? Ref. amplitude 1H	0.000 V

Sequence - Part 1

Bandwidth	500 Hz/Px
-----------	-----------

Sequence - Special

Adiab IR Pulse	Tan/Tanh
RF Phase Spoiling	ON

Sequence - Assistant

Mode	Off
------	-----

\\CMIC\Physics\T2MES\T2MES_low_res\IR GRE TI=200

TA: 24:00 PM: FIX Voxel size: 1.9×1.9×8.0 mmRel. SNR: 1.00 : fl

Properties

Prio recon	Off
Load images to viewer	On
Inline movie	Off
Auto store images	On
Load images to stamp segments	Off
Load images to graphic segments	Off
Auto open inline display	Off
Auto close inline display	Off
Start measurement without further preparation	On
Wait for user to start	Off
Start measurements	Single measurement

Routine

Slice group	1
Slices	1
Position	Isocenter
Orientation	Transversal
Phase enc. dir.	A >> P
AutoAlign	---
FoV read	360 mm
FoV phase	75.0 %
Slice thickness	8.0 mm
TR	10000.00 ms
TE	10.00 ms
Filter	Distortion Corr.(2D)
Coil elements	BO1-3;SP2-4

Contrast - Common

TR	10000.00 ms
TE	10.00 ms
Magn. preparation	Non-sel. IR
TI	200 ms
Flip angle	90 deg

Resolution - Common

FoV read	360 mm
FoV phase	75.0 %
Slice thickness	8.0 mm
Base resolution	192
Phase resolution	100 %

Resolution - Filter Image

Image Filter	Off
Distortion Corr.	On
Mode	2D
Unfiltered images	Off
Prescan Normalize	Off
Normalize	Off
B1 filter	Off

Resolution - Filter Rawdata

Raw filter	Off
Elliptical filter	Off

Geometry - Common

Slice group	1
Slices	1
Position	Isocenter
Orientation	Transversal

Geometry - Common

Phase enc. dir.	A >> P
FoV read	360 mm
FoV phase	75.0 %
Slice thickness	8.0 mm
TR	10000.00 ms

Geometry - AutoAlign

Slice group	1
Position	Isocenter
Orientation	Transversal
Phase enc. dir.	A >> P
AutoAlign	---
Initial Position	Isocenter
L	0.0 mm
P	0.0 mm
H	0.0 mm
Initial Rotation	0.00 deg
Initial Orientation	Transversal

System - Miscellaneous

Positioning mode	FIX
Table position	H
Table position	0 mm
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Coil Combine Mode	Adaptive Combine
Save uncombined	Off
Matrix Optimization	Off
Coil Focus	Flat
AutoAlign	---
Coil Select Mode	Default

System - Adjustments

B0 Shim mode	Cardiac
Adjust with body coil	On
Confirm freq. adjustment	Off
Assume Dominant Fat	Off
Assume Silicone	Off
Adjustment Tolerance	Auto

System - Adjust Volume

! Position	Isocenter
! Orientation	Transversal
! Rotation	0.00 deg
! A >> P	100 mm
! R >> L	100 mm
! F >> H	50 mm
Reset	Off

System - Tx/Rx

Frequency 1H	63.673772 MHz
Correction factor	1
Gain	High
Img. Scale Cor.	1.000
Reset	Off
? Ref. amplitude 1H	0.000 V

Sequence - Part 1

Bandwidth	500 Hz/Px
-----------	-----------

Sequence - Special

Adiab IR Pulse	Tan/Tanh
RF Phase Spoiling	ON

Sequence - Assistant

Mode	Off
------	-----

\\CMIC\Physics\T2MES\T2MES_low_res\IR GRE TI=400

TA: 24:00 PM: FIX Voxel size: 1.9×1.9×8.0 mmRel. SNR: 1.00 : fl

Properties

Prio recon	Off
Load images to viewer	On
Inline movie	Off
Auto store images	On
Load images to stamp segments	Off
Load images to graphic segments	Off
Auto open inline display	Off
Auto close inline display	Off
Start measurement without further preparation	On
Wait for user to start	Off
Start measurements	Single measurement

Routine

Slice group	1
Slices	1
Position	Isocenter
Orientation	Transversal
Phase enc. dir.	A >> P
AutoAlign	---
FoV read	360 mm
FoV phase	75.0 %
Slice thickness	8.0 mm
TR	10000.00 ms
TE	10.00 ms
Filter	Distortion Corr.(2D)
Coil elements	BO1-3;SP2-4

Contrast - Common

TR	10000.00 ms
TE	10.00 ms
Magn. preparation	Non-sel. IR
TI	400 ms
Flip angle	90 deg

Resolution - Common

FoV read	360 mm
FoV phase	75.0 %
Slice thickness	8.0 mm
Base resolution	192
Phase resolution	100 %

Resolution - Filter Image

Image Filter	Off
Distortion Corr.	On
Mode	2D
Unfiltered images	Off
Prescan Normalize	Off
Normalize	Off
B1 filter	Off

Resolution - Filter Rawdata

Raw filter	Off
Elliptical filter	Off

Geometry - Common

Slice group	1
Slices	1
Position	Isocenter
Orientation	Transversal

Geometry - Common

Phase enc. dir.	A >> P
FoV read	360 mm
FoV phase	75.0 %
Slice thickness	8.0 mm
TR	10000.00 ms

Geometry - AutoAlign

Slice group	1
Position	Isocenter
Orientation	Transversal
Phase enc. dir.	A >> P
AutoAlign	---
Initial Position	Isocenter
L	0.0 mm
P	0.0 mm
H	0.0 mm
Initial Rotation	0.00 deg
Initial Orientation	Transversal

System - Miscellaneous

Positioning mode	FIX
Table position	H
Table position	0 mm
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Coil Combine Mode	Adaptive Combine
Save uncombined	Off
Matrix Optimization	Off
Coil Focus	Flat
AutoAlign	---
Coil Select Mode	Default

System - Adjustments

B0 Shim mode	Cardiac
Adjust with body coil	On
Confirm freq. adjustment	Off
Assume Dominant Fat	Off
Assume Silicone	Off
Adjustment Tolerance	Auto

System - Adjust Volume

! Position	Isocenter
! Orientation	Transversal
! Rotation	0.00 deg
! A >> P	100 mm
! R >> L	100 mm
! F >> H	50 mm
Reset	Off

System - Tx/Rx

Frequency 1H	63.673772 MHz
Correction factor	1
Gain	High
Img. Scale Cor.	1.000
Reset	Off
? Ref. amplitude 1H	0.000 V

Sequence - Part 1

Bandwidth	500 Hz/Px
-----------	-----------

Sequence - Special

Adiab IR Pulse	Tan/Tanh
RF Phase Spoiling	ON

Sequence - Assistant

Mode	Off
------	-----

\\CMIC\Physics\T2MES\T2MES_low_res\IR GRE TI=600

TA: 24:00 PM: FIX Voxel size: 1.9×1.9×8.0 mmRel. SNR: 1.00 : fl

Properties

Prio recon	Off
Load images to viewer	On
Inline movie	Off
Auto store images	On
Load images to stamp segments	Off
Load images to graphic segments	Off
Auto open inline display	Off
Auto close inline display	Off
Start measurement without further preparation	On
Wait for user to start	Off
Start measurements	Single measurement

Routine

Slice group	1
Slices	1
Position	Isocenter
Orientation	Transversal
Phase enc. dir.	A >> P
AutoAlign	---
FoV read	360 mm
FoV phase	75.0 %
Slice thickness	8.0 mm
TR	10000.00 ms
TE	10.00 ms
Filter	Distortion Corr.(2D)
Coil elements	BO1-3;SP2-4

Contrast - Common

TR	10000.00 ms
TE	10.00 ms
Magn. preparation	Non-sel. IR
TI	600 ms
Flip angle	90 deg

Resolution - Common

FoV read	360 mm
FoV phase	75.0 %
Slice thickness	8.0 mm
Base resolution	192
Phase resolution	100 %

Resolution - Filter Image

Image Filter	Off
Distortion Corr.	On
Mode	2D
Unfiltered images	Off
Prescan Normalize	Off
Normalize	Off
B1 filter	Off

Resolution - Filter Rawdata

Raw filter	Off
Elliptical filter	Off

Geometry - Common

Slice group	1
Slices	1
Position	Isocenter
Orientation	Transversal

Geometry - Common

Phase enc. dir.	A >> P
FoV read	360 mm
FoV phase	75.0 %
Slice thickness	8.0 mm
TR	10000.00 ms

Geometry - AutoAlign

Slice group	1
Position	Isocenter
Orientation	Transversal
Phase enc. dir.	A >> P
AutoAlign	---
Initial Position	Isocenter
L	0.0 mm
P	0.0 mm
H	0.0 mm
Initial Rotation	0.00 deg
Initial Orientation	Transversal

System - Miscellaneous

Positioning mode	FIX
Table position	H
Table position	0 mm
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Coil Combine Mode	Adaptive Combine
Save uncombined	Off
Matrix Optimization	Off
Coil Focus	Flat
AutoAlign	---
Coil Select Mode	Default

System - Adjustments

B0 Shim mode	Cardiac
Adjust with body coil	On
Confirm freq. adjustment	Off
Assume Dominant Fat	Off
Assume Silicone	Off
Adjustment Tolerance	Auto

System - Adjust Volume

! Position	Isocenter
! Orientation	Transversal
! Rotation	0.00 deg
! A >> P	100 mm
! R >> L	100 mm
! F >> H	50 mm
Reset	Off

System - Tx/Rx

Frequency 1H	63.673772 MHz
Correction factor	1
Gain	High
Img. Scale Cor.	1.000
Reset	Off
? Ref. amplitude 1H	0.000 V

Sequence - Part 1

Bandwidth	500 Hz/Px
-----------	-----------

Sequence - Special

Adiab IR Pulse	Tan/Tanh
RF Phase Spoiling	ON

Sequence - Assistant

Mode	Off
------	-----

\\CMIC\Physics\T2MES\T2MES_low_res\IR GRE TI=800

TA: 24:00 PM: FIX Voxel size: 1.9×1.9×8.0 mmRel. SNR: 1.00 : fl

Properties

Prio recon	Off
Load images to viewer	On
Inline movie	Off
Auto store images	On
Load images to stamp segments	Off
Load images to graphic segments	Off
Auto open inline display	Off
Auto close inline display	Off
Start measurement without further preparation	On
Wait for user to start	Off
Start measurements	Single measurement

Routine

Slice group	1
Slices	1
Position	Isocenter
Orientation	Transversal
Phase enc. dir.	A >> P
AutoAlign	---
FoV read	360 mm
FoV phase	75.0 %
Slice thickness	8.0 mm
TR	10000.00 ms
TE	10.00 ms
Filter	Distortion Corr.(2D)
Coil elements	BO1-3;SP2-4

Contrast - Common

TR	10000.00 ms
TE	10.00 ms
Magn. preparation	Non-sel. IR
TI	800 ms
Flip angle	90 deg

Resolution - Common

FoV read	360 mm
FoV phase	75.0 %
Slice thickness	8.0 mm
Base resolution	192
Phase resolution	100 %

Resolution - Filter Image

Image Filter	Off
Distortion Corr.	On
Mode	2D
Unfiltered images	Off
Prescan Normalize	Off
Normalize	Off
B1 filter	Off

Resolution - Filter Rawdata

Raw filter	Off
Elliptical filter	Off

Geometry - Common

Slice group	1
Slices	1
Position	Isocenter
Orientation	Transversal

Geometry - Common

Phase enc. dir.	A >> P
FoV read	360 mm
FoV phase	75.0 %
Slice thickness	8.0 mm
TR	10000.00 ms

Geometry - AutoAlign

Slice group	1
Position	Isocenter
Orientation	Transversal
Phase enc. dir.	A >> P
AutoAlign	---
Initial Position	Isocenter
L	0.0 mm
P	0.0 mm
H	0.0 mm
Initial Rotation	0.00 deg
Initial Orientation	Transversal

System - Miscellaneous

Positioning mode	FIX
Table position	H
Table position	0 mm
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Coil Combine Mode	Adaptive Combine
Save uncombined	Off
Matrix Optimization	Off
Coil Focus	Flat
AutoAlign	---
Coil Select Mode	Default

System - Adjustments

B0 Shim mode	Cardiac
Adjust with body coil	On
Confirm freq. adjustment	Off
Assume Dominant Fat	Off
Assume Silicone	Off
Adjustment Tolerance	Auto

System - Adjust Volume

! Position	Isocenter
! Orientation	Transversal
! Rotation	0.00 deg
! A >> P	100 mm
! R >> L	100 mm
! F >> H	50 mm
Reset	Off

System - Tx/Rx

Frequency 1H	63.673772 MHz
Correction factor	1
Gain	High
Img. Scale Cor.	1.000
Reset	Off
? Ref. amplitude 1H	0.000 V

Sequence - Part 1

Bandwidth	500 Hz/Px
-----------	-----------

Sequence - Special

Adiab IR Pulse	Tan/Tanh
RF Phase Spoiling	ON

Sequence - Assistant

Mode	Off
------	-----

\\CMIC\Physics\T2MES\T2MES_low_res\IR GRE TI=1000

TA: 24:00 PM: FIX Voxel size: 1.9×1.9×8.0 mmRel. SNR: 1.00 : fl

Properties

Prio recon	Off
Load images to viewer	On
Inline movie	Off
Auto store images	On
Load images to stamp segments	Off
Load images to graphic segments	Off
Auto open inline display	Off
Auto close inline display	Off
Start measurement without further preparation	On
Wait for user to start	Off
Start measurements	Single measurement

Routine

Slice group	1
Slices	1
Position	Isocenter
Orientation	Transversal
Phase enc. dir.	A >> P
AutoAlign	---
FoV read	360 mm
FoV phase	75.0 %
Slice thickness	8.0 mm
TR	10000.00 ms
TE	10.00 ms
Filter	Distortion Corr.(2D)
Coil elements	BO1-3;SP2-4

Contrast - Common

TR	10000.00 ms
TE	10.00 ms
Magn. preparation	Non-sel. IR
TI	1000 ms
Flip angle	90 deg

Resolution - Common

FoV read	360 mm
FoV phase	75.0 %
Slice thickness	8.0 mm
Base resolution	192
Phase resolution	100 %

Resolution - Filter Image

Image Filter	Off
Distortion Corr.	On
Mode	2D
Unfiltered images	Off
Prescan Normalize	Off
Normalize	Off
B1 filter	Off

Resolution - Filter Rawdata

Raw filter	Off
Elliptical filter	Off

Geometry - Common

Slice group	1
Slices	1
Position	Isocenter
Orientation	Transversal

Geometry - Common

Phase enc. dir.	A >> P
FoV read	360 mm
FoV phase	75.0 %
Slice thickness	8.0 mm
TR	10000.00 ms

Geometry - AutoAlign

Slice group	1
Position	Isocenter
Orientation	Transversal
Phase enc. dir.	A >> P
AutoAlign	---
Initial Position	Isocenter
L	0.0 mm
P	0.0 mm
H	0.0 mm
Initial Rotation	0.00 deg
Initial Orientation	Transversal

System - Miscellaneous

Positioning mode	FIX
Table position	H
Table position	0 mm
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Coil Combine Mode	Adaptive Combine
Save uncombined	Off
Matrix Optimization	Off
Coil Focus	Flat
AutoAlign	---
Coil Select Mode	Default

System - Adjustments

B0 Shim mode	Cardiac
Adjust with body coil	On
Confirm freq. adjustment	Off
Assume Dominant Fat	Off
Assume Silicone	Off
Adjustment Tolerance	Auto

System - Adjust Volume

! Position	Isocenter
! Orientation	Transversal
! Rotation	0.00 deg
! A >> P	100 mm
! R >> L	100 mm
! F >> H	50 mm
Reset	Off

System - Tx/Rx

Frequency 1H	63.673772 MHz
Correction factor	1
Gain	High
Img. Scale Cor.	1.000
Reset	Off
? Ref. amplitude 1H	0.000 V

Sequence - Part 1

Bandwidth	500 Hz/Px
-----------	-----------

Sequence - Special

Adiab IR Pulse	Tan/Tanh
RF Phase Spoiling	ON

Sequence - Assistant

Mode	Off
------	-----

\\CMIC\Physics\T2MES\T2MES_low_res\IR GRE TI=1300

TA: 24:00 PM: FIX Voxel size: 1.9×1.9×8.0 mmRel. SNR: 1.00 : fl

Properties

Prio recon	Off
Load images to viewer	On
Inline movie	Off
Auto store images	On
Load images to stamp segments	Off
Load images to graphic segments	Off
Auto open inline display	Off
Auto close inline display	Off
Start measurement without further preparation	On
Wait for user to start	Off
Start measurements	Single measurement

Routine

Slice group	1
Slices	1
Position	Isocenter
Orientation	Transversal
Phase enc. dir.	A >> P
AutoAlign	---
FoV read	360 mm
FoV phase	75.0 %
Slice thickness	8.0 mm
TR	10000.00 ms
TE	10.00 ms
Filter	Distortion Corr.(2D)
Coil elements	BO1-3;SP2-4

Contrast - Common

TR	10000.00 ms
TE	10.00 ms
Magn. preparation	Non-sel. IR
TI	1300 ms
Flip angle	90 deg

Resolution - Common

FoV read	360 mm
FoV phase	75.0 %
Slice thickness	8.0 mm
Base resolution	192
Phase resolution	100 %

Resolution - Filter Image

Image Filter	Off
Distortion Corr.	On
Mode	2D
Unfiltered images	Off
Prescan Normalize	Off
Normalize	Off
B1 filter	Off

Resolution - Filter Rawdata

Raw filter	Off
Elliptical filter	Off

Geometry - Common

Slice group	1
Slices	1
Position	Isocenter
Orientation	Transversal

Geometry - Common

Phase enc. dir.	A >> P
FoV read	360 mm
FoV phase	75.0 %
Slice thickness	8.0 mm
TR	10000.00 ms

Geometry - AutoAlign

Slice group	1
Position	Isocenter
Orientation	Transversal
Phase enc. dir.	A >> P
AutoAlign	---
Initial Position	Isocenter
L	0.0 mm
P	0.0 mm
H	0.0 mm
Initial Rotation	0.00 deg
Initial Orientation	Transversal

System - Miscellaneous

Positioning mode	FIX
Table position	H
Table position	0 mm
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Coil Combine Mode	Adaptive Combine
Save uncombined	Off
Matrix Optimization	Off
Coil Focus	Flat
AutoAlign	---
Coil Select Mode	Default

System - Adjustments

B0 Shim mode	Cardiac
Adjust with body coil	On
Confirm freq. adjustment	Off
Assume Dominant Fat	Off
Assume Silicone	Off
Adjustment Tolerance	Auto

System - Adjust Volume

! Position	Isocenter
! Orientation	Transversal
! Rotation	0.00 deg
! A >> P	100 mm
! R >> L	100 mm
! F >> H	50 mm
Reset	Off

System - Tx/Rx

Frequency 1H	63.673772 MHz
Correction factor	1
Gain	High
Img. Scale Cor.	1.000
Reset	Off
? Ref. amplitude 1H	0.000 V

Sequence - Part 1

Bandwidth	500 Hz/Px
-----------	-----------

Sequence - Special

Adiab IR Pulse	Tan/Tanh
RF Phase Spoiling	ON

Sequence - Assistant

Mode	Off
------	-----

\\CMIC\Physics\T2MES\T2MES_low_res\IR GRE TI=1700

TA: 24:00 PM: FIX Voxel size: 1.9×1.9×8.0 mmRel. SNR: 1.00 : fl

Properties

Prio recon	Off
Load images to viewer	On
Inline movie	Off
Auto store images	On
Load images to stamp segments	Off
Load images to graphic segments	Off
Auto open inline display	Off
Auto close inline display	Off
Start measurement without further preparation	On
Wait for user to start	Off
Start measurements	Single measurement

Routine

Slice group	1
Slices	1
Position	Isocenter
Orientation	Transversal
Phase enc. dir.	A >> P
AutoAlign	---
FoV read	360 mm
FoV phase	75.0 %
Slice thickness	8.0 mm
TR	10000.00 ms
TE	10.00 ms
Filter	Distortion Corr.(2D)
Coil elements	BO1-3;SP2-4

Contrast - Common

TR	10000.00 ms
TE	10.00 ms
Magn. preparation	Non-sel. IR
TI	1700 ms
Flip angle	90 deg

Resolution - Common

FoV read	360 mm
FoV phase	75.0 %
Slice thickness	8.0 mm
Base resolution	192
Phase resolution	100 %

Resolution - Filter Image

Image Filter	Off
Distortion Corr.	On
Mode	2D
Unfiltered images	Off
Prescan Normalize	Off
Normalize	Off
B1 filter	Off

Resolution - Filter Rawdata

Raw filter	Off
Elliptical filter	Off

Geometry - Common

Slice group	1
Slices	1
Position	Isocenter
Orientation	Transversal

Geometry - Common

Phase enc. dir.	A >> P
FoV read	360 mm
FoV phase	75.0 %
Slice thickness	8.0 mm
TR	10000.00 ms

Geometry - AutoAlign

Slice group	1
Position	Isocenter
Orientation	Transversal
Phase enc. dir.	A >> P
AutoAlign	---
Initial Position	Isocenter
L	0.0 mm
P	0.0 mm
H	0.0 mm
Initial Rotation	0.00 deg
Initial Orientation	Transversal

System - Miscellaneous

Positioning mode	FIX
Table position	H
Table position	0 mm
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Coil Combine Mode	Adaptive Combine
Save uncombined	Off
Matrix Optimization	Off
Coil Focus	Flat
AutoAlign	---
Coil Select Mode	Default

System - Adjustments

B0 Shim mode	Cardiac
Adjust with body coil	On
Confirm freq. adjustment	Off
Assume Dominant Fat	Off
Assume Silicone	Off
Adjustment Tolerance	Auto

System - Adjust Volume

! Position	Isocenter
! Orientation	Transversal
! Rotation	0.00 deg
! A >> P	100 mm
! R >> L	100 mm
! F >> H	50 mm
Reset	Off

System - Tx/Rx

Frequency 1H	63.673772 MHz
Correction factor	1
Gain	High
Img. Scale Cor.	1.000
Reset	Off
? Ref. amplitude 1H	0.000 V

Sequence - Part 1

Bandwidth	500 Hz/Px
-----------	-----------

Sequence - Special

Adiab IR Pulse	Tan/Tanh
RF Phase Spoiling	ON

Sequence - Assistant

Mode	Off
------	-----

\\CMIC\Physics\T2MES\T2MES_low_res\IR GRE TI=2100

TA: 24:00 PM: FIX Voxel size: 1.9×1.9×8.0 mmRel. SNR: 1.00 : fl

Properties

Prio recon	Off
Load images to viewer	On
Inline movie	Off
Auto store images	On
Load images to stamp segments	Off
Load images to graphic segments	Off
Auto open inline display	Off
Auto close inline display	Off
Start measurement without further preparation	On
Wait for user to start	Off
Start measurements	Single measurement

Routine

Slice group	1
Slices	1
Position	Isocenter
Orientation	Transversal
Phase enc. dir.	A >> P
AutoAlign	---
FoV read	360 mm
FoV phase	75.0 %
Slice thickness	8.0 mm
TR	10000.00 ms
TE	10.00 ms
Filter	Distortion Corr.(2D)
Coil elements	BO1-3;SP2-4

Contrast - Common

TR	10000.00 ms
TE	10.00 ms
Magn. preparation	Non-sel. IR
TI	2100 ms
Flip angle	90 deg

Resolution - Common

FoV read	360 mm
FoV phase	75.0 %
Slice thickness	8.0 mm
Base resolution	192
Phase resolution	100 %

Resolution - Filter Image

Image Filter	Off
Distortion Corr.	On
Mode	2D
Unfiltered images	Off
Prescan Normalize	Off
Normalize	Off
B1 filter	Off

Resolution - Filter Rawdata

Raw filter	Off
Elliptical filter	Off

Geometry - Common

Slice group	1
Slices	1
Position	Isocenter
Orientation	Transversal

Geometry - Common

Phase enc. dir.	A >> P
FoV read	360 mm
FoV phase	75.0 %
Slice thickness	8.0 mm
TR	10000.00 ms

Geometry - AutoAlign

Slice group	1
Position	Isocenter
Orientation	Transversal
Phase enc. dir.	A >> P
AutoAlign	---
Initial Position	Isocenter
L	0.0 mm
P	0.0 mm
H	0.0 mm
Initial Rotation	0.00 deg
Initial Orientation	Transversal

System - Miscellaneous

Positioning mode	FIX
Table position	H
Table position	0 mm
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Coil Combine Mode	Adaptive Combine
Save uncombined	Off
Matrix Optimization	Off
Coil Focus	Flat
AutoAlign	---
Coil Select Mode	Default

System - Adjustments

B0 Shim mode	Cardiac
Adjust with body coil	On
Confirm freq. adjustment	Off
Assume Dominant Fat	Off
Assume Silicone	Off
Adjustment Tolerance	Auto

System - Adjust Volume

! Position	Isocenter
! Orientation	Transversal
! Rotation	0.00 deg
! A >> P	100 mm
! R >> L	100 mm
! F >> H	50 mm
Reset	Off

System - Tx/Rx

Frequency 1H	63.673772 MHz
Correction factor	1
Gain	High
Img. Scale Cor.	1.000
Reset	Off
? Ref. amplitude 1H	0.000 V

Sequence - Part 1

Bandwidth	500 Hz/Px
-----------	-----------

Sequence - Special

Adiab IR Pulse	Tan/Tanh
RF Phase Spoiling	ON

Sequence - Assistant

Mode	Off
------	-----

\\CMIC\Physics\T2MES\T2MES_low_res\BH_FW_4e2i_GT

TA: 8.5 s PM: FIX Voxel size: 1.4×1.4×8.0 mmPAT: Off Rel. SNR: 1.00 : fl

Properties

Prio recon	Off
Load images to viewer	Off
Inline movie	Off
Auto store images	On
Load images to stamp segments	Off
Load images to graphic segments	On
Auto open inline display	Off
Auto close inline display	Off
Start measurement without further preparation	On
Wait for user to start	Off
Start measurements	Single measurement

Routine

Slice group	1
Slices	1
Dist. factor	20 %
Position	Isocenter
Orientation	Transversal
Phase enc. dir.	A >> P
AutoAlign	---
Phase oversampling	0 %
FoV read	360 mm
FoV phase	75.0 %
Slice thickness	8.0 mm
TR	287.60 ms
TE 1	1.56 ms
TE 2	2.72 ms
TE 3	3.88 ms
TE 4	5.04 ms
Averages	1
Concatenations	1
Filter	Distortion Corr.(2D)
Coil elements	BO1-3;SP2-4

Contrast - Common

TR	287.60 ms
TE 1	1.56 ms
TE 2	2.72 ms
TE 3	3.88 ms
TE 4	5.04 ms
Magn. preparation	None
Flip angle	12 deg
Fat suppr.	None
Wrap-up Magn.	None

Contrast - Dynamic

Averages	1
Averaging mode	Short term
Reconstruction	Magnitude
Measurements	1
Multiple series	Off

Resolution - Common

FoV read	360 mm
FoV phase	75.0 %
Slice thickness	8.0 mm
Base resolution	256
Phase resolution	75 %
Phase partial Fourier	Off

Resolution - Common

Trajectory	Cartesian
Interpolation	Off

Resolution - iPAT

PAT mode	None
----------	------

Resolution - Filter Image

Image Filter	Off
Distortion Corr.	On
Mode	2D
Unfiltered images	Off
Prescan Normalize	Off
Normalize	Off
B1 filter	Off

Resolution - Filter Rawdata

Raw filter	Off
Elliptical filter	Off
POCS	Off

Geometry - Common

Slice group	1
Slices	1
Dist. factor	20 %
Position	Isocenter
Orientation	Transversal
Phase enc. dir.	A >> P
FoV read	360 mm
FoV phase	75.0 %
Slice thickness	8.0 mm
TR	287.60 ms
Multi-slice mode	Sequential
Series	Interleaved
Concatenations	1

Geometry - AutoAlign

Slice group	1
Position	Isocenter
Orientation	Transversal
Phase enc. dir.	A >> P
AutoAlign	---
Initial Position	Isocenter
L	0.0 mm
P	0.0 mm
H	0.0 mm
Initial Rotation	0.00 deg
Initial Orientation	Transversal

Geometry - Saturation

Fat suppr.	None
Wrap-up Magn.	None
Special sat.	None

Geometry - Navigator**System - Miscellaneous**

Positioning mode	FIX
Table position	H
Table position	0 mm

System - Miscellaneous

MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Coil Combine Mode	Adaptive Combine
Save uncombined	Off
Matrix Optimization	Off
Coil Focus	Flat
AutoAlign	---
Coil Select Mode	Default

System - Adjustments

B0 Shim mode	Cardiac
Adjust with body coil	On
Confirm freq. adjustment	Off
Assume Dominant Fat	Off
Assume Silicone	Off
Adjustment Tolerance	Auto

System - Adjust Volume

! Position	Isocenter
! Orientation	Transversal
! Rotation	0.00 deg
! A >> P	100 mm
! R >> L	100 mm
! F >> H	50 mm
Reset	Off

System - Tx/Rx

Frequency 1H	63.673772 MHz
Correction factor	1
Gain	High
Img. Scale Cor.	1.000
Reset	Off
? Ref. amplitude 1H	0.000 V

Physio - Signal1

1st Signal/Mode	ECG/Trigger
Average cycle	947 ± 122 ms
Average cycle	No Signal ms
Captured cycle	947 ± 122 ms
Acquisition window	831 ms
Trigger pulse	1
Trigger delay	543 ms
TR	287.60 ms
Concatenations	1
Segments	20
Phases	1
Adaptive Triggering	Off

Physio - Cardiac

Tagging	None
Magn. preparation	None
Fat suppr.	None
Dark blood	Off
FoV read	360 mm
FoV phase	75.0 %
Phase resolution	75 %
Cine	Off
Trajectory	Cartesian
Dummy heartbeats	1

Physio - PACE

Resp. control	Off
---------------	-----

Physio - PACE

Concatenations	1
----------------	---

Inline - Common

Subtract	Off
Measurements	1
StdDev	Off
Save original images	On

Inline - Cardiac

Inline Evaluation	Off
Magn. preparation	None
Contrasts	4
TE 1	1.56 ms
TE 2	2.72 ms
TE 3	3.88 ms
TE 4	5.04 ms
TR	287.60 ms
Save original images	On

Inline - MIP

MIP-Sag	Off
MIP-Cor	Off
MIP-Tra	Off
MIP-Time	Off
Save original images	On

Inline - Composing

Inline Composing	Off
Distortion Corr.	On
Mode	2D
Unfiltered images	Off

Sequence - Part 1

Introduction	Off
Dimension	2D
Reordering	Linear
Asymmetric echo	Off
Contrasts	4
Flow comp. 1	No
Readout mode	Monopolar
Optimization	None
Multi-slice mode	Sequential
Echo spacing	14.4 ms
Sequence type	Gre
Bandwidth 1	977 Hz/Px
Bandwidth 2	977 Hz/Px
Bandwidth 3	977 Hz/Px
Bandwidth 4	977 Hz/Px

Sequence - Part 2

Define	Segments
Segments	20
RF pulse type	Fast
Gradient mode	Fast
Excitation	Slice-sel.
Flip angle mode	Constant
RF spoiling	On
Phase Enc. Rewinder	On
Cine	Off

Sequence - Special

FatWater Separation	On
Multi-echo Images	Off

Sequence - Special

In-Opp Phase Images	Off
Frequency Map	Off
T2* Map	Off
Motion Correction	Off
MoCo Averaging Mode	Complex MoCo
MoCo Images Only?	Off
Frequency Offset	0 Hz
No. of Interleaves	2

Sequence - Assistant

Mode	Off
Allowed delay	0 s

\\CMIC\Physics\T2MES\T2MES_low_res\highres_3D_of_tubes_x5

TA: 4:49 PM: FIX Voxel size: 1.0×1.0×1.0 mmPAT: 2 Rel. SNR: 1.00 : tfl

Properties

Prio recon	Off
Load images to viewer	On
Inline movie	Off
Auto store images	On
Load images to stamp segments	On
Load images to graphic segments	On
Auto open inline display	Off
Auto close inline display	Off
Start measurement without further preparation	On
Wait for user to start	Off
Start measurements	Single measurement

Routine

Slab group	1
Slabs	1
Dist. factor	20 %
Position	Isocenter
Orientation	Transversal
Phase enc. dir.	A >> P
AutoAlign	---
Phase oversampling	0 %
Slice oversampling	20.0 %
Slices per slab	120
FoV read	256 mm
FoV phase	64.8 %
Slice thickness	1.00 mm
TR	400.90 ms
TE	1.77 ms
Averages	5
Concatenations	1
Filter	None
Coil elements	BO1-3;SP2-4

Contrast - Common

TR	400.90 ms
TE	1.77 ms
Magn. preparation	None
Flip angle	12 deg
Fat suppr.	None
Wrap-up Magn.	None

Contrast - Dynamic

Averages	5
Averaging mode	Short term
Reconstruction	Magnitude
Measurements	1
Multiple series	Off

Resolution - Common

FoV read	256 mm
FoV phase	64.8 %
Slice thickness	1.00 mm
Base resolution	256
Phase resolution	100 %
Slice resolution	100 %
Phase partial Fourier	Off
Slice partial Fourier	Off
Trajectory	Cartesian
Interpolation	Off

Resolution - iPAT

PAT mode	GRAPPA
Accel. factor PE	2
Ref. lines PE	24
Accel. factor 3D	1
Reference scan mode	Integrated

Resolution - Filter Image

Image Filter	Off
Distortion Corr.	Off
Prescan Normalize	Off
Normalize	Off
B1 filter	Off

Resolution - Filter Rawdata

Raw filter	Off
Elliptical filter	Off
POCS	Off

Geometry - Common

Slab group	1
Slabs	1
Dist. factor	20 %
Position	Isocenter
Orientation	Transversal
Phase enc. dir.	A >> P
Slice oversampling	20.0 %
Slices per slab	120
FoV read	256 mm
FoV phase	64.8 %
Slice thickness	1.00 mm
TR	400.90 ms
Multi-slice mode	Sequential
Series	Descending
Concatenations	1

Geometry - AutoAlign

Slab group	1
Position	Isocenter
Orientation	Transversal
Phase enc. dir.	A >> P
AutoAlign	---
Initial Position	Isocenter
L	0.0 mm
P	0.0 mm
H	0.0 mm
Initial Rotation	0.00 deg
Initial Orientation	Transversal

Geometry - Saturation

Fat suppr.	None
Wrap-up Magn.	None
Special sat.	None

Geometry - Navigator**System - Miscellaneous**

Positioning mode	FIX
Table position	H
Table position	0 mm
MSMA	S - C - T

System - Miscellaneous

Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Coil Combine Mode	Sum of Squares
Save uncombined	Off
Matrix Optimization	Off
Coil Focus	Flat
AutoAlign	---
Coil Select Mode	Default

System - Adjustments

B0 Shim mode	Cardiac
Adjust with body coil	On
Confirm freq. adjustment	Off
Assume Dominant Fat	Off
Assume Silicone	Off
Adjustment Tolerance	Auto

System - Adjust Volume

! Position	Isocenter
! Orientation	Transversal
! Rotation	0.00 deg
! A >> P	100 mm
! R >> L	100 mm
! F >> H	50 mm
Reset	Off

System - Tx/Rx

Frequency 1H	63.673772 MHz
Correction factor	1
Gain	Low
Img. Scale Cor.	1.000
Reset	Off
? Ref. amplitude 1H	0.000 V

Physio - Signal1

1st Signal/Mode	None
TR	400.90 ms
Concatenations	1
Segments	95

Physio - Cardiac

Tagging	None
Magn. preparation	None
Fat suppr.	None
Dark blood	Off
FoV read	256 mm
FoV phase	64.8 %
Phase resolution	100 %
Cine	Off
Trajectory	Cartesian
Dummy heartbeats	1

Physio - PACE

Resp. control	Off
Concatenations	1

Inline - Common

Subtract	Off
Measurements	1
StdDev	Off
Save original images	On

Inline - Cardiac

Inline Evaluation	Off
Magn. preparation	None
Contrasts	1
TE	1.77 ms
TR	400.90 ms
Save original images	On

Inline - MIP

MIP-Sag	Off
MIP-Cor	Off
MIP-Tra	Off
MIP-Time	Off
Save original images	On

Inline - Composing

Inline Composing	Off
Distortion Corr.	Off

Sequence - Part 1

Introduction	Off
Dimension	3D
Elliptical scanning	Off
Reordering	Linear
Asymmetric echo	Off
Contrasts	1
Flow comp.	No
Optimization	Min. TE TR
Multi-slice mode	Sequential
Echo spacing	4.2 ms
Sequence type	Gre
Bandwidth	592 Hz/Px

Sequence - Part 2

Define	Shots
Shots per slice	1
Segments	95
RF pulse type	Fast
Gradient mode	Fast
Excitation	Slab-sel.
Flip angle mode	Constant
RF spoiling	On
Phase Enc. Rewinder	On
Cine	Off

Sequence - Assistant

Mode	Off
Allowed delay	0 s

\\CMIC\Physics\T2MES\T2MES_low_res\FAMAP120HIGHBWTP_90

TA: 4:30 PM: REF Voxel size: 2.8×2.8×8.0 mmRel. SNR: 1.00 : fl

Properties

Prio recon	Off
Load images to viewer	On
Inline movie	Off
Auto store images	On
Load images to stamp segments	Off
Load images to graphic segments	Off
Auto open inline display	Off
Auto close inline display	Off
Start measurement without further preparation	On
Wait for user to start	Off
Start measurements	Single measurement

Routine

Slice group	1
Slices	1
Dist. factor	20 %
Position	Isocenter
Orientation	Transversal
Phase enc. dir.	A >> P
AutoAlign	---
Phase oversampling	0 %
FoV read	360 mm
FoV phase	75.0 %
Slice thickness	8.0 mm
TR	5000.0 ms
TE	4.2 ms
Averages	1
Concatenations	1
Filter	Distortion Corr.(2D)
Coil elements	BO1-3;SP2-4

Contrast - Common

TR	5000.0 ms
TE	4.2 ms
MTC	Off
Flip angle	90 deg
Fat suppr.	None
Water suppr.	None

Contrast - Dynamic

Averages	1
Averaging mode	Short term
Reconstruction	Magnitude
Measurements	1
Multiple series	Off

Resolution - Common

FoV read	360 mm
FoV phase	75.0 %
Slice thickness	8.0 mm
Base resolution	128
Phase resolution	75 %
Phase partial Fourier	6/8
Interpolation	Off

Resolution - Filter Image

Image Filter	Off
Distortion Corr.	On
Mode	2D

Resolution - Filter Image

Unfiltered images	Off
Prescan Normalize	Off
Normalize	Off
B1 filter	Off

Resolution - Filter Rawdata

Raw filter	Off
Elliptical filter	Off

Geometry - Common

Slice group	1
Slices	1
Dist. factor	20 %
Position	Isocenter
Orientation	Transversal
Phase enc. dir.	A >> P
FoV read	360 mm
FoV phase	75.0 %
Slice thickness	8.0 mm
TR	5000.0 ms
Multi-slice mode	Sequential
Series	Ascending
Concatenations	1

Geometry - AutoAlign

Slice group	1
Position	Isocenter
Orientation	Transversal
Phase enc. dir.	A >> P
AutoAlign	---
Initial Position	Isocenter
L	0.0 mm
P	0.0 mm
H	0.0 mm
Initial Rotation	0.00 deg
Initial Orientation	Transversal

Geometry - Saturation

Fat suppr.	None
Water suppr.	None
Special sat.	None

System - Miscellaneous

Positioning mode	REF
Table position	H
Table position	0 mm
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Coil Combine Mode	Sum of Squares
Save uncombined	Off
Matrix Optimization	Off
Coil Focus	Flat
AutoAlign	---
Coil Select Mode	Default

System - Adjustments

B0 Shim mode	Standard
Adjust with body coil	On

System - Adjustments

Confirm freq. adjustment	Off
Assume Dominant Fat	Off
Assume Silicone	Off
Adjustment Tolerance	Auto

System - Adjust Volume

! Position	Isocenter
! Orientation	Transversal
! Rotation	0.00 deg
! A >> P	100 mm
! R >> L	100 mm
! F >> H	25 mm
Reset	Off

System - Tx/Rx

Frequency 1H	63.673772 MHz
Correction factor	1
Gain	High
Img. Scale Cor.	1.000
Reset	Off
? Ref. amplitude 1H	0.000 V

Physio - Signal1

1st Signal/Mode	None
TR	5000.0 ms
Concatenations	1

Inline - Common

Subtract	Off
Measurements	1
StdDev	Off
Save original images	On

Inline - MIP

MIP-Sag	Off
MIP-Cor	Off
MIP-Tra	Off
MIP-Time	Off
Save original images	On

Inline - Composing

Inline Composing	Off
Distortion Corr.	On
Mode	2D
Unfiltered images	Off

Sequence - Part 1

Introduction	Off
Dimension	2D
Contrasts	1
Multi-slice mode	Sequential
Bandwidth	1560 Hz/Px

Sequence - Part 2

Gradient mode	Fast
RF spoiling	On

Sequence - Special

Online ICE	Off
Selection box	Second Choice
Spoil me!	On
Test Time	400 ms
dARRAY [1]	2.0 [UnitArr]

Sequence - Special

dARRAY [2]	12.0 [UnitArr]
dARRAY [3]	22.00 [UnitArr]

Sequence - Assistant

Mode	Off
------	-----

\\CMIC\Physics\T2MES\T2MES_low_res\FAMAP60HIGHBWTP_45

TA: 4:30 PM: FIX Voxel size: 2.8×2.8×8.0 mmRel. SNR: 1.00 : fl

Properties

Prio recon	Off
Load images to viewer	On
Inline movie	Off
Auto store images	On
Load images to stamp segments	Off
Load images to graphic segments	Off
Auto open inline display	Off
Auto close inline display	Off
Start measurement without further preparation	On
Wait for user to start	Off
Start measurements	Single measurement

Routine

Slice group	1
Slices	1
Dist. factor	20 %
Position	Isocenter
Orientation	Transversal
Phase enc. dir.	A >> P
AutoAlign	---
Phase oversampling	0 %
FoV read	360 mm
FoV phase	75.0 %
Slice thickness	8.0 mm
TR	5000.0 ms
TE	4.2 ms
Averages	1
Concatenations	1
Filter	Distortion Corr.(2D)
Coil elements	BO1-3;SP2-4

Contrast - Common

TR	5000.0 ms
TE	4.2 ms
MTC	Off
Flip angle	45 deg
Fat suppr.	None
Water suppr.	None

Contrast - Dynamic

Averages	1
Averaging mode	Short term
Reconstruction	Magnitude
Measurements	1
Multiple series	Off

Resolution - Common

FoV read	360 mm
FoV phase	75.0 %
Slice thickness	8.0 mm
Base resolution	128
Phase resolution	75 %
Phase partial Fourier	6/8
Interpolation	Off

Resolution - Filter Image

Image Filter	Off
Distortion Corr.	On
Mode	2D

Resolution - Filter Image

Unfiltered images	Off
Prescan Normalize	Off
Normalize	Off
B1 filter	Off

Resolution - Filter Rawdata

Raw filter	Off
Elliptical filter	Off

Geometry - Common

Slice group	1
Slices	1
Dist. factor	20 %
Position	Isocenter
Orientation	Transversal
Phase enc. dir.	A >> P
FoV read	360 mm
FoV phase	75.0 %
Slice thickness	8.0 mm
TR	5000.0 ms
Multi-slice mode	Sequential
Series	Ascending
Concatenations	1

Geometry - AutoAlign

Slice group	1
Position	Isocenter
Orientation	Transversal
Phase enc. dir.	A >> P
AutoAlign	---
Initial Position	Isocenter
L	0.0 mm
P	0.0 mm
H	0.0 mm
Initial Rotation	0.00 deg
Initial Orientation	Transversal

Geometry - Saturation

Fat suppr.	None
Water suppr.	None
Special sat.	None

System - Miscellaneous

Positioning mode	FIX
Table position	H
Table position	0 mm
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Coil Combine Mode	Sum of Squares
Save uncombined	Off
Matrix Optimization	Off
Coil Focus	Flat
AutoAlign	---
Coil Select Mode	Default

System - Adjustments

B0 Shim mode	Standard
Adjust with body coil	Off

System - Adjustments

Confirm freq. adjustment	Off
Assume Dominant Fat	Off
Assume Silicone	Off
Adjustment Tolerance	Auto

System - Adjust Volume

! Position	Isocenter
! Orientation	Transversal
! Rotation	0.00 deg
! A >> P	100 mm
! R >> L	100 mm
! F >> H	25 mm
Reset	Off

System - Tx/Rx

Frequency 1H	63.673772 MHz
Correction factor	1
Gain	High
Img. Scale Cor.	1.000
Reset	Off
? Ref. amplitude 1H	0.000 V

Physio - Signal1

1st Signal/Mode	None
TR	5000.0 ms
Concatenations	1

Inline - Common

Subtract	Off
Measurements	1
StdDev	Off
Save original images	On

Inline - MIP

MIP-Sag	Off
MIP-Cor	Off
MIP-Tra	Off
MIP-Time	Off
Save original images	On

Inline - Composing

Inline Composing	Off
Distortion Corr.	On
Mode	2D
Unfiltered images	Off

Sequence - Part 1

Introduction	Off
Dimension	2D
Contrasts	1
Multi-slice mode	Sequential
Bandwidth	1560 Hz/Px

Sequence - Part 2

Gradient mode	Fast
RF spoiling	On

Sequence - Special

Online ICE	Off
Selection box	Second Choice
Spoil me!	On
Test Time	400 ms
dARRAY [1]	2.0 [UnitArr]

Sequence - Special

dARRAY [2]	12.0 [UnitArr]
dARRAY [3]	22.00 [UnitArr]

Sequence - Assistant

Mode	Off
------	-----